Mycoplasma pneumoniae-IgA-ELISA medac

Deutsch/English/Français

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Mycoplasma pneumoniae-IgA-ELISA medac

Enzyme immunossay for quantitative determination of
Mycoplasma pneumoniae IgA antibodies in human serum

Cat.no.: 361

FOR IN VITRO DIAGNOSTIC USE ONLY

INTRODUCTION

Mycoplasmas are bacteria without cell walls (Class Mollicutes = soft-skinned, order Mycoplasmatales) having a particularly small genome. They are characterised by their small size (300 - 800 nm) and their capacity to assume widely differing shapes (pleomorphism). At the present time over 100 species of the Order Mycoplasmatales are known. Most of them occur exclusively in animals. Fourteen species have so far been isolated from human beings. They colonise the mucous membranes of the respiratory and urogenital tracts. Only a few of them are pathogenic. **Mycoplasma pneumoniae** is among the species pathogenic to human beings. Two variants of **M. pneumoniae** are known.

M. pneumoniae is an extracellular parasite of the mucosa of the respiratory tract. The pathogen is distinguished by its high level of host specificity. A specific adhesin is responsible for the adherence of the pathogen to the surface of the host cell. This adhesin also constitutes the virulence factor against which the crucial humoral response is directed. The transmission of **M. pneumoniae** takes place by droplet infection, but the incubation period of 10 to 20 days is relatively long.

M. pneumoniae is endemic, occurring with seasonal peaks in the spring and autumn. Epidemics at intervals of 3 - 4 years are not unusual. Infections with **M. pneumoniae** are among the most frequent causes of community-acquired tracheobronchitis and atypical pneumonia in children and young people. The highest incidence is found between the ages of five and fifteen years. Of all cases of atypical pneumonia, 5 - 10 % are attributable to **M. pneumoniae**. It also causes conditions such as pharyngitis, laryngitis, otitis media and myringitis. Following a **M. pneumoniae** infection which begins with an infection of the respiratory tract, various other clinical manifestations may occur in other organs and systems: these include pericarditis and myocarditis, reactive arthritis, meningitis, meningoencephalitis, and

polyneuritis together with a wide assortment of conditions affecting the skin. This at first seemingly harmless infection can therefore entail serious late consequences. The disease may run a particularly severe course in immunodeficient and immunosuppressed patients.

A previous attack of ***M. pneumoniae*** infection does not leave behind it any immunity to this pathogen. Frequent reinfections are therefore observed.

In children under five years of age the infection is predominantly asymptomatic. In such cases the disease runs a mild course. In older children and in young people and adults, however, the infection is accompanied by prostration, headache, fever and a dry, unproductive cough. But these clinical features are totally nonspecific and do not provide any clue to the cause of the infection. An effective diagnostic test is therefore necessary.

In the acute phase of the disease demonstration of the causative pathogen is the method of choice. Tests for the pathogen can be performed on throat swabs, specimens of sputum and bronchoalveolar lavage fluid (BAL). The classical technique for demonstrating ***M. pneumoniae*** by means of cell culture takes a long time (10 - 14 days). Nevertheless the isolation of the microorganism is successful in only 40 - 60 % of cases. The antigen-ELISA test requires material rich in cells, because the sensitivity of the test is relatively low. Effective qualitative DNA tests for this pathogen are not yet commercially available.

Serology constitutes the method of choice for patients with chronic infections, reinfections or extrapulmonary diseases due to ***M. pneumoniae***. Commercial test systems include the complement fixation test (CFT), the haemagglutination test (HAT) and enzyme immunoassay (ELISA). Unlike the CFT and the HAT, the ELISA techniques make it possible to differentiate between immunoglobulin classes IgG, IgA and IgM, so that acute infections can be distinguished from chronic infections and past infections.


The ***Mycoplasma pneumoniae*-IgA-ELISA medac** uses a recombinant antigen mixture for serological testing. This ensures that the results will have high specificity and sensitivity. Furthermore, the medac one point quantification (AU/ml) offers the best possible conditions for reproducible results and hence for the measurement of serial changes in titre.

In addition to the **Mycoplasma pneumoniae-IgA-ELISA medac**
Cat.no.: **361** for 96 determinations


we also distribute the following products:

Mycoplasma pneumoniae-IgG-ELISA medac
Cat.no.: **360** for 96 determinations,
Mycoplasma pneumoniae-IgM-ELISA medac
Cat.no.: **362** for 96 determinations.


TEST PRINCIPLE


The plate is coated with ***M. pneumoniae*** specific recombinant antigen.

The Mycoplasma-specific antibodies from the specimen are selectively bound to the antigen.

Peroxidase-conjugated anti-human IgA antibodies bind to the IgA antibodies (P = peroxidase).

Incubation with TMB-Substrate (*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

Advantages of the test

- ☞ The breakable microwell strips permit efficient use of the test.
- ☞ Suitable for automation on open ELISA devices.
- ☞ One-point quantification, no standard curve needed.

KIT CONTENTS

Cat. no.: 361

1.

MTP

Microplate: 12 x 8 wells, yellow-coded, (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with *M. pneumoniae*-specific antigen and FCS, ready to use.
2.

CONTROL	-
----------------	----------

Negative control: 2 vials with 0.75 ml each, human serum, ready to use, stained blue, contains NBCS, phenol, ProClin™ 300 and gentamycin sulfate.
3.

CONTROL	+
----------------	----------

Positive control: 2 vials with 0.75 ml each, human serum, ready to use, stained blue, contains NBCS, phenol, ProClin™ 300 and gentamycin sulfate.
4.

CAL

Calibrator: 2 vials with 0.75 ml each, human serum, ready to use, stained blue, contains NBCS, phenol, ProClin™ 300 and gentamycin sulfate.
5.

WB

Wash buffer: 1 bottle with 100 ml PBS/Tween (10 x), pH 7.2 -7.4, contains ProClin™ 300.
6.

BAC-DIL

Sample diluent: 1 bottle with 110 ml PBS/Tween/NBCS, pH 7.0 - 7.2, ready to use, stained blue, contains ProClin™ 300.
7.

CON

Conjugate: 3 vials with 5.0 ml each, goat anti-human IgA, HRP-conjugated, ready to use, stained yellow, contains BSA, phenol, ProClin™ 300 and gentamycin sulfate.
8.

TMB

TMB-substrate: 1 vial with 10 ml, ready to use.
9.

STOP

Stop solution: 2 vials with 14 ml each, 0.5 M sulfuric acid (H₂SO₄), ready to use.

1. STORAGE AND STABILITY

MATERIAL/REAGENT	STATE	STORAGE	STABILITY
Test kit	unopened	2 - 8 °C	until expiry date
Microplate	opened	2 - 8 °C in bag with desiccant	12 weeks
Controls/Calibrator	opened	2 - 8 °C	12 weeks
Wash buffer	diluted	2 - 8 °C	12 weeks
Sample diluent	opened	2 - 8 °C	12 weeks
Conjugate	opened	2 - 8 °C	12 weeks
TMB-substrate	opened	2 - 8 °C	12 weeks
Stop solution	opened	2 - 8 °C	until expiry date

Do not use the reagents after the expiry date.

2. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

- 2.1. Water for injection (H₂O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of wash buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37 °C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

3. PREPARATION OF THE REAGENTS

Before starting the test procedure all kit components must be equilibrated to room temperature.

Calculate the number of wells required.

3.1. Microplate

The aluminium bag has to be tightly resealed together with the desiccant after each removal of wells. Storage and stability of the wells are indicated under point 1.

3.2. Wash buffer

Mix one volume of wash buffer (10 x) with nine volumes of water for injection (e.g. 50 ml wash buffer (10 x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

Crystals in the wash buffer (10 x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.

Do not mix test specific reagents (microplate, controls, calibrator, conjugate) from different kit lots. In contrast to that, sample diluent, wash buffer, TMB-substrate and stop solution are generally exchangeable in all Chlamydia-and Mycoplasma-ELISA medac.

Reagents from other manufacturers must not be used in general.

Valid and reproducible results are only obtained if the test procedure is precisely followed.

4. SPECIMEN

4.1. The test is suitable for serum.

4.2. Pretreatment of sera, e.g. inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain any red blood cells.

4.3. Sera have to be diluted 1:100 with sample diluent.

5.A. TEST PROCEDURE

5.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

Microplate wells are ready to use and do not have to be pre-washed.

5.2. Pipette 50 µl sample diluent into well A1 as blank (see 6.A.). Add 50 µl each of the diluted sample as well as negative control, positive control and calibrator in duplicates to the wells.

If necessary the microplate wells can be kept in a humid chamber up to 30 min at RT before proceeding.

5.3. Incubate the microplate wells for 60 min (± 5 min) at 37 °C (± 1 °C) in a humid chamber or sealed with incubation cover foil.

5.4. After incubation wash the microplate wells three times with

200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

Do not allow the wells to dry out! Proceed immediately!

5.5. Add conjugate (coloured yellow) to each well.

50 µl of conjugate have to be pipetted into the wells if the test is done manually.

Please note:

When working with automated devices, 60 µl of conjugate have to be pipetted into each well due to a higher evaporation in the incubation chambers of the devices.

The suitability of the test for automated devices was confirmed during the evaluation of the test. Nevertheless we recommend to verify the compatibility of the test with the devices used in the lab.

5.6. Incubate again for 60 min (\pm 5 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil.

5.7. After incubation wash microplate wells again (see 5.4.).

5.8. Add 50 µl of TMB-substrate to each well and incubate for 30 min (\pm 2 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.

5.9. Stop the reaction by adding 100 µl of stop solution to each well. Positive samples turn yellow.

Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.

The reading must be done within 15 min after adding the stop solution!

5.B. TABLE FOR THE TEST PROCEDURE

	Blank (A1)	Negative Control	Positive Control	Calibrator	Sample
Sample diluent	50 µl	-	-	-	-
Neg. Control	-	50 µl	-	-	-
Pos. Control	-	-	50	-	-
Calibrator	-	-	-	50 µl	-
Sample	-	-	-	-	50 µl
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer					
Conjugate	50/60 µl *)	50/60 µl *)	50/60 µl *)	50/60 µl *)	50/60 µl *)
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer					
TMB-Substrate	50 µl	50 µl	50 µl	50 µl	50 µl
Incubate for 30 min at 37 °C in the dark					
Stop Solution	100 µl	100 µl	100 µl	100 µl	100 µl
Photometric reading at 450 nm (ref. 620 - 650 nm)					

*) manual/automatic procedure (see 5.5.)

6.A. CALCULATION OF RESULTS (VALIDITY)

- * The evaluation is performed using arbitrary units (AU).
- * Read OD values at 450 nm (reference wavelength 620 - 650 nm).
- * Subtract the OD value of the blank (well A1) from all other OD values.
- * Lot-specific data

The lot-specific data sheet provided with the kit contains the following information:

- Lot-specific calibration curve
- Curve parameters a and b
- Nominal OD value of the calibrator
- Lower OD limit of the calibrator
- Nominal concentration range (AU/ml) of the positive control

* Validity criteria

- The OD value of the **blank** has to be < 0.100.
- The mean OD value of the **negative control** has to be < 0.100.
- The unit value of the **positive control** has to be within the nominal range indicated in the lot-specific data sheet.
- The mean OD value of the calibrator has to be above the lower OD limit indicated in the lot-specific data sheet.

Repeat the run if the results do not meet the specification!

* Correction of the results

The measured OD values of the positive control and the samples have to be corrected as follows:

$$OD_{\text{corrected}} = \frac{\text{Nominal OD value of the calibrator}}{\text{Measured OD value of the calibrator}} \times OD_{\text{measured}}$$

* Quantification of the results

The corresponding concentrations of the corrected OD values in AU/ml can be read from the lot-specific calibration curve (see lot-specific data sheet).

Alternatively, the concentrations can be calculated using the following formula:

$$\text{Concentration [AU/ml]} = b \left(\frac{a}{OD_{\text{corrected}}} - 1 \right)$$

Most of the new ELISA readers allow to program the formula, thus enabling an automated data processing.

The measuring range spans from 9 to 125 AU/ml. Samples below this range have to be interpreted as < 9 AU/ml, those above as > 125 AU/ml. These values must not be extrapolated.

The cut-off is 10 AU/ml.

Grey zone = 9 - 11 AU/ml.

6.B. INTERPRETATION OF RESULTS/LIMITATIONS OF THE METHOD

- * Samples with Unit values below the lower limit of the grey zone are reported as **NEGATIVE**.
- * Samples with Unit values within the grey zone are reported as **EQUIVOCAL**.

These samples should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.
- * Samples with Unit values exceeding the upper limit of the grey zone are reported as **POSITIVE**.
- * The results should always be interpreted in connection with clinical data and additional diagnostic parameters.
- * High concentrations of hemoglobin in serum do not have an influence on the results. However, high concentrations of lipids may influence the test results.
- * Cross-reactivities with heterophilic antibodies cannot be excluded in individual cases.

6.C. SPECIFIC IgM-/IgA-/IgG-INTERPRETATION

Possible results			Interpretation
IgM	IgA	IgG	
+	-	-	1. Indication of early stage of infection. Retest IgM, IgA and IgG after 14 days. ^{1,2}
-	+	-	2. Indication of early stage of infection or solitary, persisting IgA. Retest IgA and IgG after 14 days.
+	+	-	3. Indication of acute infection. ^{1,2} Retest after 14 days.
+	-	+	4. Indication of acute infection. ^{1,2}
+	+	+	5. Indication of acute infection. ^{1,2}
-	+	+	6. Indication of current infection. ³
-	-	+	7. Indication of past infection. In case of clinical suspicion retest second serum sample after 14 days for IgA and IgG antibodies.
-	-	-	8. No serological indication of current or past infection. In case of clinical suspicion retest second serum sample after 14 days for IgM, IgA and IgG antibodies.

Comment:

- Borderline results may indicate commencing or abating infections. Retesting after 14 days is recommended.
- ¹ Current, acute infections are best detectable by parallel determination of IgM and IgA antibodies.
- ² Simultaneous detection of IgM and IgA antibodies is particularly frequent in children.
- ³ In adults, IgA antibodies are more reliable markers for current infections than IgM antibodies.

7. PERFORMANCE CHARACTERISTICS

We determined the following performance characteristics during the diagnostic evaluation.

7.A. SENSITIVITY AND SPECIFICITY

Patient Group	Specificity IgA
<i>M. pneumoniae</i> IgA negative sera (reference-ELISA/Aggl. Assay)	98 % (n=40)
Patient Group	Sensitivity IgA
Sera with IgA antibodies to <i>M. pneumoniae</i> (reference ELISA/Aggl. Assay) Patients with respiratory infections	59 % (n=35)

7.B. PRECISION

Sample	Intra-assay variation				Sample	Inter-assay variation			
	mean AU	SD	CV (%)	n		mean AU	SD	CV (%)	n
PC	12.3	0.7	6	22	PC	13	0.5	4	13
N° 1	20.6	0.6	3	22	N° 4	21.7	1.1	5	13
N° 2	21.7	0.7	3	22	N° 5	20.8	1.2	6	13
N° 3	116	3.2	3	22	N° 6	114.7	11.2	10	13

PC = positive control

GENERAL HANDLING ADVICES

- * To avoid cross contamination do not exchange the vials and their screw caps.
- * The reagents have to be sealed immediately after use to avoid evaporation and microbial contamination.
- * After use, the reagents have to be stored as indicated to guarantee the shelf life.
- * After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

HEALTH AND SAFETY INFORMATION

- * The local occupational safety and health regulations have to be regarded.
- * Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents), should be handled as potentially infectious and used with all necessary precautions.

DISPOSAL CONSIDERATIONS

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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LITERATUR/REFERENCES/LITTERATURE

Clyde, WA: Clinical overview of typical ***Mycoplasma pneumoniae*** infections. Clin Infect Dis 17 (Suppl 1), S32-6 (1993).

Dionisio D, Valassina M, Uberti M, Fabbri C, Parri F, Saffi EG: ***Mycoplasma pneumoniae*** non-pulmonary infection presenting with pharyngitis, polyarthritits and localized exanthem. Scand J Infect Dis 33, 782-783 (2001).

Döller G, Döller PC, Jacobs E, Schuy W: Zur Differentialdiagnostik von Infektionen des Respirationstrakts. Diagnose & Labor 43, 21-33 (1993).

Drasbek M, Nielsen PK, Persson K, Birkelund S, Christiansen G: Immune response to ***Mycoplasma pneumoniae*** P1 and P116 in patients with atypical pneumonia analyzed by ELISA. BMC Microbiol. 4, 7ff. (2004).

Ferwerda A, Moll HA, de Groot R: Respiratory tract infections by ***Mycoplasma pneumoniae*** in children: a review of diagnostic and therapeutic measures. Eur J Pediatr 160, 483-491 (2001).

Granstrom M, Holme T, Sjogren AM, Ortqvist A, Kalin M: The role of IgA determination by ELISA in the early serodiagnosis of ***Mycoplasma pneumoniae*** infection, in relation to IgG and mu-capture IgM methods. Med Microbiol. 40, 288-92 (1994).

Hammerschlag MR: ***Mycoplasma pneumoniae*** infections. Curr Opin Infect Dis 14, 181-186 (2001).

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Krause D: ***Mycoplasma pneumoniae*** cytoadherence: organization and assembly of the attachment organelle. Trends Microbiol 6, 15-18 (1998).

Seggev JS, Semak GV, Kurup VP: Isotype-specific antibody responses to acute ***Mycoplasma pneumoniae*** infection. Ann Allergy Asthma Immunol. 77, 66-73 (1996).

Sillis M: The limitations of IgM assays in the serological diagnosis of ***Mycoplasma pneumoniae*** infections. Med Microbiol. 33, 253-8 (1990).

Sotgiu S, Pugliatti M, Rosati G, Deiana GA, Sechi GP: Neurological disorders associated with ***Mycoplasma pneumoniae*** infection. Eur J Neurol 10, 165-168 (2003).