

CORMAY RF CALIBRATORS

CALIBRATORS FOR DETERMINATION OF RHEUMATOID FACTOR LEVELS



Kit name	Kit size	Cat. No
CORMAY RF CALIBRATORS	4 x 1 ml	4-277

INTRODUCTION

Rheumatoid factor (RF) is an autoantibody against human IgG commonly seen in sera at a high concentration in some conditions, particularly in patients with rheumatoid arthritis. The measurement of RF value is useful in evaluating the diagnosis, effects of therapy and prognosis of RA, systemic lupus erythematosus, chronic hepatopathy, etc.

CORMAY RF CALIBRATORS kit is prepared by diluting rheumatoid factor (RF) with buffer solution containing 1 w/v% bovine serum albumin at various concentrations. Concentrations are shown on each vial and in this kit insert.

CORMAY RF CALIBRATORS kit should be used for the preparation of the calibration curve for PRESTIGE 24i RF, CORMAY RF and CORMAY RF 500 kits.

PACKAGE

Cap colour	RF level	Size
BLUE	20 U/ml	1 x 1ml
GREEN	40 U/ml	1 x 1ml
RED	80 U/ml	1 x 1ml
ORANGE	120 U/ml	1 x 1ml

The calibrators are ready to use.

STABILITY AND STORAGE

Unopened calibrator is stable at 2-10⁰C until the expiration date printed on the vial. Protect from light.

WARNINGS AND NOTES

- Product for in vitro diagnostic use only.
- As the colours of caps vary according to RF levels in the vial, care should be taken not to interchange the caps.
- The calibrators should not be mixed with ones in the other vial.
- The calibrator should be used according to this kit insert.
- The calibrator is based on human serum that has been tested for the HIV antibody, HBsAg and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.

WASTE MANAGEMENT

Please refer to local legal requirements.

Date of issue: 10. 2004.

MANUFACTURER

P.Z. CORMAY
ul. Rapackiego 19, 20-150 Lublin, POLAND
P.O. Box 122, 20-954 Lublin 2
tel.: +48 (0) 81 749 44 00
fax: +48 (0) 81 749 44 34
<http://www.pzcormay.pl>

10/04/10/04