

CORMAY Lp(a) CONTROL N
CORMAY Lp(a) CONTROL P



CONTROL SERUM (LYOPHILISED)
FOR DETERMINATION
OF LIPOPROTEIN (a) CONCENTRATION

CORMAY Lp(a) CONTROL N **5 x 1 ml** **Cat. No 4-492**
CORMAY Lp(a) CONTROL P **5 x 1 ml** **Cat. No 4-493**

CORMAY Lp(a) CONTROL N and CORMAY Lp(a) CONTROL P are lyophilised human sera containing lipoprotein (a) at a concentration shown in the label attached on each bottle.

RECONSTITUTION

Carefully remove the rubber cap from vial and add 1 ml distilled water to the vial. Recap and leave at room temperature for 15 minutes, inverting gently to dissolve contents.

STABILITY AND STORAGE

1. Lyophilised serum remains stable when stored at 2-10°C until expiry date given on the product label.
2. The constituents in reconstituted control serum are stable for 2 weeks at 2-10°C.

NOTES

1. Product for in vitro diagnostic use only.
2. Once reconstituted control sera are ready for use.
3. When removing the rubber cap, avoid dispersing any lyophilized powder attached to the insides of the cap. Use appropriate laboratory gloves, clothing, etc. when handling and wash hands thoroughly after use.
4. This serum has been tested for the HIV antibody, HBsAg and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
5. The reagents contain sodium azide (< 0,1%) as a preservative. Avoid contact with skin and mucous membranes.
6. If the reagent is insufficient for measurement, addition of a new reagent should be avoided because that may cause measurement error.
7. This product should be used according to this kit insert.

WASTE MANAGEMENT

Please refer to local legal requirements.

Date of issue: 03. 2005.

MANUFACTURER

P.Z. CORMAY

ul. Rapackiego 19, 20-150 Lublin, POLAND
P.O. Box 122, 20-954 Lublin 2
tel.: +48 (0) 81 749 44 00
fax: +48 (0) 81 749 44 34
<http://www.pzcormay.pl>

03/05/03/05