

CORMAY MYOGLOBIN

DIAGNOSTIC KIT FOR DETERMINATION OF MYOGLOBIN CONCENTRATION



Kit name	Kit size	Cat. No
CORMAY MYOGLOBIN	1 x 84 ml	6-301

INTRODUCTION

Myoglobin (Mb) is a hemo-protein present in cardiac and skeletal muscle cells and is released into blood circulation when these cells are damaged. The determination of serum Mb level is useful in the diagnosis of myocardial infarction, muscular dystrophy, myositis and myopathy, and also for the assessment of treatment and disease prognosis.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between Mb in a sample and anti-Mb antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of Mb in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package	
1-Reagent	1 x 48 ml
2-Reagent	1 x 36 ml

Reagent preparation and stability

The reagents are ready to use.

The reagents are stable up to the kit expiry date printed on the package when stored at 2-10°C. Protect from light!

Concentrations in the test

suspension of latex particles sensitized with anti-Mb (rabbit) antibodies (pH 7.3) 0.12 w/v%
glycine buffer solution (pH 9.0)

Warnings and notes

- Product for in vitro diagnostic use only.
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- The reagents contain sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.

ADDITIONAL EQUIPMENT

- automated clinical chemistry analyser capable of accommodating two-reagent assays;
- general laboratory equipment;

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid).

It is recommended to perform the assay with freshly collected samples. If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided.

Samples which contain an excessive amount of Mb should be diluted with physiological saline and re-tested.

PROCEDURE

wavelength	572 nm
temperature	37°C

These reagents may be used in automatic analysers according to their service manual. These reagents may be used directly in HITACHI 911/912 analysers. Applications for analysers are available on request.

REFERENCE VALUES

serum, plasma	< 65 ng/ml
---------------	------------

It is recommended for each laboratory to establish its own reference ranges for local population. Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL II (Cat. No 4-290) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY MYOGLOBIN CALBRATORS kit (Cat. No 4-279) is recommended. Renewed calibration is recommended: after 1 month when using the reagent on the analyser, after lot change, as required. The standard solutions should be measured at least twice.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analyser TBA-30R. Results may vary if a different instrument is used.

- Analytical range:** 20 – 1000 ng/ml.
- Antigen excess:** up to 24000 ng/ml.
- Specificity / Interferences**
Haemoglobin up to 500 mg/dl, bilirubin up to 30 mg/dl, triglycerides up to 1500 mg/dl, RF up to 500 U/ml do not interfere with the test.

Precision

Repeatability (within run) n = 10	Mean [ng/ml]	SD [ng/ml]	CV [%]
level 1	48.6	2.3	4.63
level 2	86.2	3.3	3.87
level 3	323.5	3.8	1.17

Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using 43 samples gave following results:

$$y = 1.13x - 9.0 \text{ ng/ml};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Galvin J. P. et al.: Particle enhanced photometric immunoassay systems., Clin. Lab. Assays (Pap. Annu. Clin. Lab. Assays Conf.), 4th, 73 (1983).
2. Singer J. M. et al.: The latex fixation test. I. Application to the serologic diagnosis of rheumatoid arthritis, Amer. J. Med., 21, 888 (1956).

Date of issue: 10. 2004.

MANUFACTURER

P.Z. CORMAY

ul. Rapackiego 19, 20-150 Lublin, POLAND

P.O. Box 122, 20-954 Lublin 2

tel.: +48 (0) 81 749 44 00

fax: +48 (0) 81 749 44 34

<http://www.pzcormay.pl>

10/04/10/04