## Rubella-IgG-ELISA PKS medac

---

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s. S. 28

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**Rubella-IgG-ELISA PCS medac**

Enzyme immunoassay with Pipetting Control System (PCS) for the quantitative detection of IgG antibodies to rubella virus in serum and cerebrospinal fluid (CSF)

Cat. no.: 136-PKS

FOR IN VITRO DIAGNOSTIC USE ONLY

**INTRODUCTION**

The pathogen causing rubella (German measles), rubella virus, belongs to the family of *Togaviridae* and is assigned to the genus *Rubivirus*. The rubella virus has a single-stranded RNA genome and occurs only in one serotype.

Rubella appears mainly in childhood and takes normally a mild course. The disease is characterised by a small-spotted exanthema and is accompanied by lymphadenosis.


During the first trimenon of pregnancy, infections with the rubella virus can cause a foetal infection, which can lead to irreversible embryofoetal damages.

An immunization against a rubella virus infection can be achieved with suitable vaccines. With an immunization a protective immunity can be achieved in about 95 % of the vaccinees. Different test methods are available to measure the success of the vaccination as well as to determine the status of immunity.


The detection of rubella virus-specific IgG antibodies using the quantitative Rubella-IgG-ELISA PCS medac can be carried out fast and easily. By using test controls, which are calibrated against the international reference of the WHO, it is possible to quantify the specific anti-rubella virus IgG titer. Low titers in a hemagglutination inhibition assay (HAI) can be controlled by the quantitative detection of specific IgG antibodies.

The test is also suitable for the detection of pathogen-specific intrathecal antibody synthesis calculating a specific antibody index in serum-CSF pairs.


## TEST PRINCIPLE


The plate is coated with rubella virus antigen.

The rubella virus-specific antibodies from the specimen are selectively bound to the antigen.

Peroxidase-conjugated anti-human IgG antibodies bind to the rubella-specific antibodies (P = peroxidase).

Incubation with TMB-substrate (\*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

### Advantages of the test

- ☞ The Pipetting Control System allows to monitor visually each pipetting step through colour changes.
- ☞ The breakable microwell strips permit efficient use of the test.
- ☞ Suitable for automation on open ELISA devices.
- ☞ One-point quantification, no standard curve needed.
- ☞ No additional calibration curve for diagnostic of CSF needed.
- ☞ Suitable for parallel detection of immune status and pathogen-specific antibody index.

**KIT CONTENTS**

**Cat. no.: 136-PKS**

1. 

<b>MTP</b>
------------

Microplate: 12 x 8 wells, red-coded (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with rubella virus antigen, BSA and pH indicator, ready to use.
2. 

<b>CONTROL</b>	<b>-</b>
----------------	----------

Negative control: 2 vials with 0.75 ml each, human serum, ready to use, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
3. 

<b>CONTROL</b>	<b>+</b>
----------------	----------

Positive control: 2 vials with 0.75 ml each, human serum, ready to use, contains FCS, BSA, phenol, ProClin™ 300 and gentamycine sulfate.
4. 

<b>CAL</b>
------------

Calibrator: 2 vials with 0.75 ml each, human serum, ready to use, contains FCS, BSA, phenol, ProClin™ 300 and gentamycine sulfate.
5. 

<b>WB</b>
-----------

Wash buffer: 1 bottle with 100 ml, PBS/Tween (10 x), pH 7.2 - 7.4, contains ProClin™ 300.
6. 

<b>VIR-DIL</b>
----------------

Sample diluent: 1 bottle with 110 ml, PBS/Tween/BSA, pH 7.2 - 7.4, ready to use, contains ProClin™ 300.
7. 

<b>CON</b>
------------

Conjugate: 3 vials with 5.0 ml each, goat anti-human IgG, HRP-conjugated, ready to use, stained green, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
8. 

<b>TMB</b>
------------

TMB-substrate: 1 vial with 10 ml, ready to use.
9. 

<b>STOP</b>
-------------

Stop solution: 2 vials with 14 ml each, 0.5 M sulfuric acid (H<sub>2</sub>SO<sub>4</sub>), ready to use.

## **1. STORAGE AND STABILITY**

<b>MATERIAL/REAGENT</b>	<b>STATE</b>	<b>STORAGE</b>	<b>STABILITY</b>
Test kit	unopened	2...8°C	until expiry date
Microplate	opened	2...8°C in bag with desiccant	6 weeks
Controls/ Calibrator	opened	2...8°C	6 weeks
Wash buffer	diluted	2...8°C	6 weeks
Sample diluent	opened	2...8°C	6 weeks
Conjugate	opened	2...8°C	6 weeks
TMB-substrate	opened	2...8°C	6 weeks
Stop solution	opened	2...8°C	until expiry date

Do not use the reagents after the expiry date.

## **2. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED**

- 2.1. Water for injection (H<sub>2</sub>O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of wash buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37 °C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

## **3. PREPARATION OF THE REAGENTS**

Before starting the test procedure all kit components must be equilibrated to room temperature.

Calculate the number of wells required.

### **3.1. Microplate**

The aluminium bag has to be tightly resealed together with the desiccant after each removal of wells. Storage and stability of the wells are indicated in table 1.

**Note: The microplate wells have a light green colour. Eventually occurring greenish brown stains inside the wells are due to the production process and do not influence the test performance.**

### 3.2. Wash buffer

Mix one volume of wash buffer (10 x) with nine volumes of water for injection (e.g. 50 ml wash buffer (10 x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

**Crystals in the wash buffer (10 x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.**

**Do not mix test specific reagents (microplate, controls, conjugate, calibrator) from different kit lots.**

**In contrast to that sample diluent, wash buffer, TMB-substrate and stop solution are generally exchangeable in all virological ELISA medac.**

**Reagents from other manufacturers must not be used in general.**

**Valid and reproducible results are only obtained if the test procedure is precisely followed.**

## **4. SPECIMEN**

4.1. The test is suitable for serum and CSF samples (for the detection of CSF, see 8.).

4.2. Pretreatment of sera, e.g. inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain any red blood cells.

4.3. Sera have to be diluted 1 : 200 with sample diluent. We recommend to prepare an initial dilution of 1 : 50 (e.g. 10 µl serum + 490 µl sample diluent). For further 1 : 4 dilution just prepare the volume needed. Samples outside the measuring range can be diluted further.

4.4 **The diagnostic investigation of serum-CSF pairs is detailed described under 8.**

## **5.A. TEST PROCEDURE**

5.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

**Microplate wells are ready to use and do not have to be pre-washed.**

- 5.2. Leave well A1 empty as blank (see 6.A.). Add 50 µl each of the negative control, the positive control as well as the diluted samples in single determination, and 50 µl of the calibrator in duplicate to the wells.

**After pipetting the samples (pH neutral or basic fluid) the wells turn blue. A missing colour change in one well indicates that no sample or control has been added.**

**If necessary, the microplate wells can be kept in a humid chamber up to 60 min at 2 - 8 °C before proceeding.**

- 5.3. Incubate the microplate wells for 60 min ( $\pm$  5 min) at 37 °C ( $\pm$  1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.4. After incubation wash the microplate wells three times with 200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

**Do not allow the wells to dry out! Proceed immediately!**

- 5.5. Add conjugate (coloured green) to each well (except A1).

**50 µl of conjugate have to be pipetted into the wells if the test is done manually.**

**Please note:**

**When working with automated devices, 60 µl of conjugate have to be pipetted into each well due to a higher evaporation in the incubation chambers of the devices.**

**The suitability of the test for automated devices was confirmed during the evaluation of the test. Nevertheless we recommend to verify the compatibility of the test with the devices used in the lab.**

- 5.6. Incubate again for 60 min ( $\pm$  5 min) at 37 °C ( $\pm$  1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.7. After incubation wash microplate wells again (see 5.4.).
- 5.8. Add 50 µl of TMB-substrate to each well (also A1) and incubate for 30 min ( $\pm$  2 min) at 37 °C ( $\pm$  1 °C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.
- 5.9. Stop the reaction by adding 100 µl of stop solution to each well (also A1). Positive samples turn yellow.

Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.

The reading should be done within 15 min after adding the stop solution!

#### **5.B. TABLE FOR THE TEST PROCEDURE**

	Blank (A1)	Negative control	Positive control	Calibrator	Sample
Neg. control	-	50 µl	-	-	-
Pos. control	-	-	50 µl	-	-
Calibrator	-	-	-	50 µl	-
Sample	-	-	-	-	50 µl
Incubate for 60 min at 37°C, wash 3 x with 200 µl wash buffer					
Conjugate	-	50/60 µl *)	50/60 µl *)	50/60 µl *)	50/60 µl *)
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer					
TMB-substrate	50 µl	50 µl	50 µl	50 µl	50 µl
Incubate for 30 min at 37 °C in the dark					
Stop solution	100 µl	100 µl	100 µl	100 µl	100 µl
Photometric reading at 450 nm (ref. 620 - 650 nm)					

\*) manual/automatic procedure (see 5.6.)

#### **6.A. CALCULATION OF RESULTS (VALIDITY)**

- \* Read OD values at 450 nm (reference wavelength 620 - 650 nm).
- \* Subtract the OD value of the blank (well A1) from all other OD values.
- \* Lot-specific data  
 The lot-specific data sheet provided with the kit contains the following information:
  - Lot-specific calibration curve
  - Curve parameters a, b and c
  - Nominal OD value of the calibrator
  - Lower OD limit of the calibrator
  - Nominal concentration range (IU/ml) of the positive control
- \* Validity criteria
  - The OD value of the **negative control** has to be < **0.150**.

- The unit value of the **positive control** has to be within the nominal range indicated in the lot-specific data sheet.
- The mean OD value of the **calibrator** has to be above the lower OD limit indicated in the lot-specific data sheet.
- Additional validity criteria for the detection of serum-CSF pairs see 8.

**Repeat the run if the results do not meet the specification!**

\* Correction of the results

The measured OD values of the positive control and the samples have to be corrected as follows:

$$OD_{\text{corrected}} = \frac{\text{Nominal OD value of the calibrator}}{\text{Measured OD value of the calibrator}} \times OD_{\text{measured}}$$

\* Quantification of the results

The corresponding concentrations of the corrected OD values in IU/ml can be read from the lot-specific calibration curve (see lot-specific data sheet).

Alternatively, the concentrations can be calculated using the following formula:

$$\text{Concentration [IU/ml]} = b / \left( \frac{a}{OD_{\text{corrected}} - c} - 1 \right)$$

Most of the new ELISA readers allow to program the formula, thus enabling an automated data processing.

The measuring range spans from 5 to 200 IU/ml. Samples below this range have to be interpreted as < 5 IU/ml, those above as > 200 IU/ml. These values must not be extrapolated.

**The cut-off for immunity is 15 IU/ml.**

**Grey zone = cut-off ± 20 % (= 12 to 18 IU/ml).**

**6.B. INTERPRETATION OF RESULTS/LIMITATIONS OF THE METHOD**

- \* Samples with Unit values below the lower limit of the grey zone are reported as **NEGATIVE** with regard to immunity.

\* Samples with Unit values within the grey zone are reported as **EQUIVOCAL**.

Due to the fact that a protective immunity is not given with equivocal results these samples should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.

\* Samples with Unit values exceeding the upper limit of the grey zone are reported as **POSITIVE** with regard to immunity.

\* The results should always be interpreted in connection with clinical data and additional diagnostic parameters.

\* Very high concentrations of lipids in positive samples can influence the measured concentration of antibodies. High concentrations of hemoglobin and bilirubin do not influence the test results.

**7. PERFORMANCE CHARACTERISTICS**

We determined the following performance characteristics during the diagnostic evaluation.

**7.A. SENSITIVITY AND SPECIFICITY**

350 sera from patients, 100 of them detected as negative with the rubella hemagglutination inhibition assay, and 250 of them detected as positive with a reference ELISA, were measured with the Rubella-IgG-ELISA PCS medac. The results obtained are shown in the table below.

	Reference test	
	negative	positive
Rubella-IgG-ELISA PCS	negative	0
	99	0
	positive	250
	1*	250

\* Sample was detected as equivocal

Sensitivity = 100 %  
 Specificity = 99.00 %

Positive predictive value: 99.60 %  
 Negative predictive value: 100 %

## 7.B. PRECISION

Sample	Intra-assay variation				Sample	Inter-assay variation			
	mean OD	SD	CV (%)	n		mean IU	SD	CV (%)	n
NC	0.029	0.002	7	22	PC	13.6	0.8	6	13
PC	0.338	0.014	4	22	N° 2	17.9	1.5	8	13
Cal	0.952	0.046	5	22	N° 3	90.4	11.3	13	13
N° 1	1.764	0.070	4	22	N° 4	145.5	18.7	13	13

NC = Negative Control; PC = Positive Control; Cal = Calibrator

## 8. DIAGNOSTIC INVESTIGATIONS OF CEREBROSPINAL FLUID (CSF)

The detection of rubella-specific antibody synthesis in the central nervous system (CNS) during the diagnostic investigation of cerebrospinal fluid is an essential part of the differential diagnostic of acute infections, and also of chronic conditions involving the CNS.

Apart from fungi and bacteria, there are numerous viruses that can be responsible for infections with CNS involvement. Furthermore, in chronic diseases such as multiple sclerosis, high antibody indices for measles, rubella and/or varicella-zoster are found in about 95% of all cases.

The identification of rubella-specific intrathecal antibody synthesis is performed by estimation of the antibody index (AI) according to Reiber (Reiber 1987, 1999). In order to calculate the AI the following conditions must be fulfilled:

- estimation of the albumin quotient ( $Q_{alb}$ ) so as to assess the function of the blood-brain-barrier and to calculate the Limes value in patients with elevated IgG quotients ( $Q_{tot} > Q_{lim}$ )
- estimation of the total IgG quotient ( $Q_{tot}$ )

### 8.1. SPECIMEN

8.1.1. The test is suitable for serum-CSF paired samples.

8.1.2. Pre-treatment of sera and CSF samples, e.g., inactivation, is not necessary. However they should neither be contaminated with microorganisms nor contain any red blood cells.

#### 8.1.3. Serum

In addition to the 1 : 200 dilution for determination of serostatus, the sera are diluted 1 : 1500 with sample diluent. To calculate the AI, a dilution should be chosen so that its IU value will lie inside the assay range (see 6.A. and 8.2.). If the IU values for both dilutions are within the assay range of 5 - 200 IU, the IU value of the 1 : 1500 dilution should be chosen to calculate the AI. If the IU values for both dilutions are above 200 IU, the sample has to be further diluted.

#### 8.1.4. Cerebrospinal fluid

The CSF samples are diluted to a standard dilution 1 : 5 with sample diluent. If the measured antibody concentration is outside the measuring range, the sample has to be further diluted, or has to be retested in a lesser dilution (max. 1 : 2).

**Serum and cerebrospinal fluid have always to be assayed in parallel in the same test run (this also applies to repeat measurements)!**

### **8.2. TEST PROCEDURE**

The further test setup for rubella-IgG determination in serum-CSF pairs is performed as described under 5.A.

### **8.3. DETERMINATION OF TOTAL IgG AND ALBUMIN CONCENTRATIONS**

In addition to the rubella-specific IgG determinations, in each sample pair the total IgG concentration and the albumin concentration in serum and CSF have to be determined.

### **8.4. CALCULATION OF RESULTS/VALIDITY**

#### \* Validity

The validity criteria specified under 6.A. are applicable.

**Repeat the run if the results do not meet the specification.**

The following points also apply to CSF investigations:

- Sera with an antibody content < 5 IU in a dilution of 1 : 200 are considered as seronegative. In such cases the antibody index cannot be determined.
- In extremely rare cases, sero-negative patients may have demonstrable intrathecal antibodies because of rubella-associated encephalopathies.
- The assay range for CSF extends from 2 - 200 IU.
- CSF samples from patients with **positive serostatus** which at the dilution of 1 : 5 are below the assay range, should be retested in a lesser dilution (maximum 1 : 2).

\* Evaluation

- Calculation of IU values see 6.A.
- Calculation of the pathogen-specific IgG quotient ( $Q_{\text{spec}}$ ).

$$Q_{\text{spec}} = \frac{\text{IU CSF} \times \text{dilution CSF}}{\text{IU serum} \times \text{dilution serum}}$$

- Calculation of the antibody index

The pathogen-specific index is calculated from the fomulae:

1.  $AI = Q_{\text{spec}}/Q_{\text{tot}}$  (for  $Q_{\text{tot}} < Q_{\text{lim}}$ )
2.  $AI = Q_{\text{spec}}/Q_{\text{lim}}$  (for  $Q_{\text{tot}} > Q_{\text{lim}}$ )
3.  $Q_{\text{lim}} = 0,93 \times \sqrt{Q_{\text{alb}}^2 + 6 \times 10^{-6}} - 1,7 \times 10^{-3}$

## 8.5. INTERPRETATION OF RESULTS

- \* AI values from 0.6 - 1.3 are regarded as the normal range.
- \* AI values > 1.3 and ≤ 1.5 are regarded as borderline.
- \* **The abnormal range is defined as AI > 1.5.**
- \* AI values < 0.6 point to analytical errors and cannot be interpreted.
- \* The CSF diagnostic criteria for an acute, active disease of the CNS are raised cell count and a raised albumin quotient. These reflect obstruction to CSF flow due to some inflammatory condition.
- \* Elevated antibody indices are not reliable evidence of the acute phase of an infective CNS disease, because antibodies, even intrathecally, may persist for long periods, and because polyspecific CNS-intrinsic antibody synthesis may occur. In some circumstances it may be advisable to look for a significant change in the AI value by testing second paired serum-CSF samples. For this purpose further sample

collection will be necessary and should be performed after a time interval chosen in the light of the clinical circumstances.

#### **GENERAL HANDLING ADVICES**

- \* To avoid cross contamination do not exchange the vials and their screw caps.
- \* The reagents have to be sealed immediately after use to avoid evaporation and microbial contamination.
- \* After use, the reagents have to be stored as indicated to guarantee the shelf life.
- \* After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

#### **HEALTH AND SAFETY INFORMATION**

- \* The local occupational safety and health regulations have to be regarded.
- \* Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents), should be handled as potentially infectious and used with all necessary precautions.

#### **DISPOSAL CONSIDERATIONS**

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

#### **REFERENCES**

see page 28

**Date of issue: 06.04.2006**

## LITERATUR/REFERENCES

- Enders, G.: Röteln. Fortschr. Klin. Virol. München, 157 - 177, (1986).
- Hermann, K.L.: Available Rubella serological tests. Rev. Infect. Dis. 7, 108 - 112, (1985).
- Reiber, H. and Felgenhauer, K.: Protein transfer at the blood cerebrospinal fluid barrier and the quantitation of the humoral immune response within the central nervous system., Clin Chim Acta, 319 - 328, (1987).
- Reiber, H.: Liquordiagnostik, in: Klinische Neurologie, Berlitz, P. (Hrsg.), Springer Verlag, Heidelberg, 148 - 177, (1999).
- Pustowoit, B., Hofmann, J. und Trauer, H.: Rötelnvirus, in: Diagnostische Bibliothek, Bd. 1, Virusdiagnostik, Porstmann, T. (Hrsg.), Blackwell Wissenschafts-Verlag, 499 - 515, (1996).
- Selb, B.: Rötelnvirus, in: Medizinische Virusdiagnostik, Selb, B. (Hrsg.), Umschau Verlag, 198 - 202, (1992).
- Wolinsky, J.: Rubella, in: Virology, 2. Edition, Field, B. N., Kipe, D. M. et al. (Eds.), Raven Press, Ltd., New York, (1990).