CMV-IgG-ELISA PKS medac

Deutsch/English/Français

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s. S. 47

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CMV-IgG-ELISA PCS medac

Enzyme immunoassay with **Pipetting Control System** (PCS) for the quantitative detection of IgG antibodies to Cytomegalovirus (CMV) in serum, plasma and cerebrospinal fluid (CSF)

Cat. no.: 115/Q/PKS

FOR IN VITRO DIAGNOSTIC USE ONLY

INTRODUCTION

Cytomegalovirus (CMV) belongs to the family of human pathogenic Herpesviridae which is characterised by a double-stranded DNA genome. It is characteristic for these viruses to persist latently in the organism after primary infection. Reactivations of the virus can therefore occur under certain circumstances.

The course, CMV infections take in immunocompetent individuals, is normally without any symptoms. In individuals who suffer from immunosuppression (e.g. transplantation patients, HIV infected individuals, tumor patients, newborns) a variety of serious symptoms are observed.

CMV is one of the most frequent pathogens in prenatal infection which can cause a variety of serious disorders, also including late damages in children born without symptoms.

Depending on the geographic area and prevalence of infection, IgG antibodies to CMV are found in 50 - 100 % of adults.

The CMV-IgG-ELISA PCS medac can be applied for the screening of the immunity of patients with a clinically suspected CMV infection.

Furthermore, the determination of CMV antibodies in blood donors is of special importance, because the transmission of cytomegalovirus via blood transfusions to immunosuppressed patients, e.g. to transplant recipients, might have serious consequences, e.g. graft rejection.


The detection of specific IgG antibodies using the quantitative CMV-IgG-ELISA PCS medac can be carried out fast and easily. By using test controls, which are calibrated against the reference of the National Centre of Blood Products in France, the Bloodbank of Lille, it is possible to quantify the specific anti-CMV-IgG level.

In addition to the **CMV-IgG-ELISA PCS medac**
Cat. no.: **115/Q/PKS** for 96 determinations


we also distribute the following products:

CMV-IgG-ELA test PCS medac
Cat. no.: **115/PKS** for 96 determinations,
CMV-IgG-ELA test PCS medac
Cat. no.: **116/PKS** for 5 x 96 determinations,
CMV-IgA-ELA test PCS medac
Cat. no.: **112/PKS** for 96 determinations,
CMV-IgM-ELA test PCS medac
Cat. no.: **110/PKS** for 96 determinations.


TEST PRINCIPLE


The plate is coated with CMV antigen.

The CMV-specific antibodies from the specimen are selectively bound to the antigen.

Peroxidase-conjugated anti-human IgG antibodies bind to the IgG antibodies (P = peroxidase).

Incubation with TMB-substrate (*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

Advantages of the test

- ☞ The Pipetting Control System allows to monitor visually each pipetting step through colour changes.
- ☞ The breakable microwell strips permit efficient use of the test.
- ☞ Suitable for automation on open ELISA devices.
- ☞ One-point quantification, no standard curve needed.
- ☞ No additional calibration curve for diagnostic of CSF needed.

KIT CONTENTS

Cat. no.: 115/Q/PKS

1.

MTP

Microplate: 12 x 8 wells, orange-coded (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with CMV antigen, BSA and pH indicator, ready to use.
2.

CONTROL	-
----------------	----------

Negative control: 2 vials with 0.75 ml each, human serum, ready to use, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
3.

CONTROL	+
----------------	----------

Positive control: 2 vials with 0.75 ml each, human serum, ready to use, contains FCS, BSA, phenol, ProClin™ 300 and gentamycine sulfate.
4.

CAL

Calibrator: 2 vials with 0.75 ml each, human serum, ready to use, contains FCS, BSA, phenol, ProClin™ 300 and gentamycine sulfate.
5.

WB

Wash buffer: 1 bottle with 100 ml, PBS/Tween (10x), pH 7.2-7.4, contains ProClin™ 300.
6.

VIR-DIL

Sample diluent: 1 bottle with 110 ml, PBS/Tween/BSA, pH 7.2-7.4, ready to use, contains ProClin™ 300.
7.

CON

Conjugate: 3 vials with 5.0 ml each, goat anti-human IgG, HRP-conjugated, ready to use, stained green, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
8.

TMB

TMB-substrate: 1 vial with 10 ml, ready to use.
9.

STOP

Stop solution: 2 vials with 14 ml each, 0.5 M sulfuric acid (H₂SO₄), ready to use.

1. STORAGE AND STABILITY

MATERIAL/REAGENT	STATE	STORAGE	STABILITY
Test kit	unopened	2...8°C	until expiry date
Microplate	opened	2...8°C in bag with desiccant	6 weeks
Controls/ Calibrator	opened	2...8°C	6 weeks
Wash buffer	diluted	2...8°C	6 weeks
Sample diluent	opened	2...8°C	6 weeks
Conjugate	opened	2...8°C	6 weeks
TMB-substrate	opened	2...8°C	6 weeks
Stop solution	opened	2...8°C	until expiry date

Do not use the reagents after the expiry date.

2. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

- 2.1. Water for injection (H₂O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of Wash Buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37°C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

3. PREPARATION OF THE REAGENTS

Before starting the test procedure all kit components must be equilibrated to room temperature.

Calculate the number of wells required.

3.1. Microplate

The aluminium bag has to be tightly resealed together with the desiccant after each removal of wells. Storage and stability of the wells are indicated under point 1.

Note: The microplate wells have a light green colour. Eventually occurring greenish brown stains inside the wells are due to the production process and do not influence the test performance.

3.2. Wash buffer

Mix one volume of wash buffer (10 x) with nine volumes of water for injection (e.g. 50 ml wash buffer (10 x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

Crystals in the wash buffer (10 x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.

Do not mix test specific reagents (microplate, controls, conjugate, calibrator) from different kit lots. In contrast to that, sample diluent, wash buffer, TMB-substrate and stop solution are generally exchangeable in all virological ELISA medac.

Reagents from other manufactures must not be used in general.

Valid and reproducible results are only obtained if the test procedure is precisely followed .

4. SPECIMEN

- 4.1. The test is suitable for serum, plasma samples and CSF samples (for the detection of CSF see 8.).
- 4.2. Pretreatment of sera, e.g. inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain any red blood cells.
- 4.3. Serum or plasma samples have to be diluted 1 : 200 with sample diluent. We recommend to prepare an initial dilution of 1 : 50 (e.g. 10 µl serum + 490 µl sample diluent). For further 1 : 4 dilution just prepare the volume needed (e.g. 20 µl serum + 60 µl sample diluent). Samples outside the measuring range can be diluted further.
- 4.4. **The diagnostic investigation of serum-CSF is described in detail under 8.**

5.A. TEST PROCEDURE

- 5.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

Microplate wells are ready to use and do not have to be pre-washed.

- 5.2. Leave well A1 empty as blank (see 6.A.). Add 50 µl each of the diluted sample as well as negative control, positive control and Calibrator in duplicates to the wells.

After pipetting the samples (pH neutral or basic fluid) the wells turn blue. A missing colour change in one well indicates that no sample or control has been added.

If necessary, the microplate wells can be kept in a humid chamber up to 60 min at 2 - 8°C before proceeding.

- 5.3. Incubate the microplate wells for 60 min (\pm 5 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.4. After incubation wash the microplate wells three times with 200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

Do not allow the wells to dry out! Proceed immediately!

- 5.5. Add conjugate (coloured green) to each well (except A1).

50 µl of conjugate have to be pipetted into the wells if the test is done manually.

Please note:

When working with automated devices, 60 µl of conjugate have to be pipetted into each well due to a higher evaporation in the incubation chambers of the devices.

The suitability of the test for automated devices was confirmed during the evaluation of the test. Nevertheless we recommend to verify the compatibility of the test with the devices used in the lab.

- 5.6. Incubate again for 60 min (\pm 5 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.7. After incubation wash microplate wells again (see 5.4.).

5.8. Add 50 µl of TMB-substrate to each well (also A1) and incubate for 30 min (± 2 min) at 37 °C (± 1 °C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.

5.9. Stop the reaction by adding 100 µl of stop solution to each well (also A1). Positive samples turn yellow.

Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.

The reading should be done within 15 min after adding the Stop Solution!

5.B. TABLE FOR THE TEST PROCEDURE

	Blank (A1)	Negative control	Positive control	Calibrator	Sample
Neg. control	-	50 µl	-	-	-
Pos. control	-	-	50 µl	-	-
Calibrator	-	-	-	50 µl	-
Sample	-	-	-	-	50 µl
Incubate for 60 min at 37°C, wash 3x with 200 µl wash buffer					
Conjugate	-	50/60 µl *)	50/60 µl *)	50/60 µl *)	50/60 µl *)
Incubate for 60 min at 37°C, wash 3x with 200 µl wash buffer					
TMB-substrate	50 µl	50 µl	50 µl	50 µl	50 µl
Incubate for 30 min at 37°C in the dark					
Stop solution	100 µl	100 µl	100 µl	100 µl	100 µl
Photometric reading at 450 nm (ref. 620 - 650 nm)					

*) manual/automatic procedure (see 5.5.)

6.A. CALCULATION OF RESULTES (VALIDITY)

- * Read OD values at 450 nm (reference wavelength 620 - 650 nm).
- * Subtract the OD value of the blank (well A1) from all other OD values.

* Lot-specific data

The lot-specific data sheet provided with the kit contains the following information:

- Lot-specific calibration curve
- Curve parameters a, b and c
- Nominal OD value of the calibrator
- Lower OD limit of the calibrator
- Nominal concentration range (AU/ml) of the positive control

* Validity criteria

- The mean OD value of the **negative control** has to be **< 0.150**.
- The unit value of the **positive control** has to be within the nominal range indicated in the lot-specific data sheet.
- The mean OD value of the **calibrator** has to be above the lower OD limit indicated in the lot-specific data sheet.
- Additional validity criteria for the detection of serum-CSF pairs see 8.

Repeat the run if the results do not meet the specification!

* Correction of the results

The measured OD values of the positive control and the samples have to be corrected as follows:

$$OD_{\text{corrected}} = \frac{\text{Nominal OD value of the calibrator}}{\text{Measured OD value of the calibrator}} \times OD_{\text{measured}}$$

* Quantification of the results

The corresponding concentrations of the corrected OD values in AU/ml can be read from the lot-specific calibration curve (see lot-specific data sheet).

Alternatively, the concentrations can be calculated using the following formula:

$$\text{Concentration [AU/ml]} = b / \left(\frac{a}{OD_{\text{corrected}} - c} - 1 \right)$$

Most of the new ELISA readers allow to program the formula, thus enabling an automated data processing.

The measuring range spans from 0.45 to 15 AU/ml. Samples above have to be interpreted as > 15 AU/ml. These values must not be

extrapolated. The samples have to be retested in higher dilutions.

The cut-off is 0.55 AU/ml

Grey zone = 0.45 - 0.65 AU/ml

Attention! Important Note!

Due to the mathematical algorithm of the quantification negative or not defined AU values can be obtained in the following cases:

- If the corrected OD values are $< c$, negative AU values are indicated. These samples have to be interpreted as negative.
- If the corrected OD value is $= c$ (not allowed division by 0) the samples have to be interpreted as negative.
- Highly positive samples with corrected OD values $\geq a + c$ are calculated as negative or not defined AU values (not allowed division by 0). These samples have to be retested in higher dilutions or have to be interpreted as > 15 AU/ml.

6.B. INTERPRETATION OF RESULTS/LIMITATIONS OF THE METHOD

- * Samples with unit values below the lower limit of the grey zone are reported as **NEGATIVE** with regard to immunity.
- * Samples with unit values within the grey zone are reported as **EQUIVOCAL**.

These samples should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.
- * Samples with unit values exceeding the upper limit of the grey zone are reported as **POSITIVE**.
- * The results should always be interpreted in connection with clinical data and additional diagnostic parameters.
- * Low false positive reactions, caused by antibodies against other Herpesviridae or by heterophilic antibodies, cannot be excluded in single cases.
- * Very high concentrations of lipids in positive samples can influence the measured concentration of antibodies. High concentrations of hemoglobin or bilirubin do not influence the test results.

* A lower antibody concentration has to be expected in citrate plasma in comparison to serum due to additional dilution.

7. PERFORMANCE CHARACTERISTICS

We determined the following performance characteristics during the diagnostic evaluation.

7.A. SENSITIVITY AND SPECIFICITY

519 samples were measured in comparison to an established reference assay. 423 samples were obtained from blood donors, 43 samples from children and 53 samples were derived from 17 pregnant women with symptoms of a primary CMV infection (12 follow ups).

The results obtained are shown in the table below:

		Reference assay		
		negative	equivocal	positive
CMV-IgG- ELISA PCS medac	negative	167	0	0
	equivocal	4	0	0
	positive	1	9	338

Specificity = 97.1 %
 Sensitivity = 100 %
 Negative predictive value : 100 %
 Positive predictive value : 97.1 %
 Concordance : 97.3 %

7.B. PRECISION

Sample	Intra-assay variation				Sample	Inter-assay variation			
	mean AU	SD	CV (%)	n		mean AU	SD	CV (%)	n
PK	1.42	0.13	9.1	21	PC	1.20	0.06	5.0	11
N° 1	0.86	0.07	8.1	21	N° 4	2.67	0.09	3.4	11
N° 2	5.86	0.54	9.2	21	N° 5	4.24	0.38	9.0	11
N° 3	11.81	0.94	8.0	21	N° 6	10.45	0.73	7.0	11

PC = Positive Control; N° 4 = CSF 1:4

8. DIAGNOSTIC INVESTIGATIONS OF CEREBROSPINAL FLUID (CSF)

The detection of CMV-specific antibody synthesis in the central nervous system (CNS) during the diagnostic investigation of cerebrospinal fluid is an essential part of the differential diagnostic of infections involving the CNS.

Apart from fungi and bacteria, there are numerous viruses that can be responsible for infections with CNS involvement.

The identification of CMV-specific intrathecal antibody synthesis is performed by estimation of the antibody index (AI) according to Reiber (Reiber 1987, 1999). In order to calculate the AI the following conditions must be fulfilled:

- estimation of the albumin quotient (Q_{alb}) so as to assess the function of the blood-brain-barrier and to calculate the Limes value in patients with elevated IgG quotients ($Q_{tot} > Q_{lim}$)
- estimation of the total IgG quotient (Q_{tot})

8.1. SPECIMEN

- 8.1.1. The test is suitable for serum-CSF paired samples.
- 8.1.2. Pre-treatment of sera and CSF samples, e.g., inactivation, is not necessary. However they should neither be contaminated with microorganisms nor contain any red blood cells.
- 8.1.3. Serum
In addition to the 1 : 200 dilution for determination of serostatus, the sera are diluted 1 : 800 with sample diluent. To calculate the AI, a dilution should be chosen so that its AU value will be positive and lie inside the measuring range (see 6.A. and 8.2.). If the AU values for both dilutions are within this range of 0.65 - 15 AU, the AU value of the 1 : 800 dilution should be chosen to calculate the AI. If the AU values for both dilutions are above 15 AU, the sample has to be further diluted.
- 8.1.4. Cerebrospinal fluid
The CSF samples are diluted to a standard dilution 1 : 4 with sample diluent. If the measured antibody concentration is outside the measuring range (0.45 - 15 AU), the sample has to be further diluted, or has to be retested in a lesser dilution (max. 1 : 2).

Serum and cerebrospinal fluid have always to be assayed in parallel in the same test run (this also applies to repeat measurements)!

8.2. TEST PROCEDURE

The further test setup for CMV-IgG determination in serum-CSF pairs is performed as described under 5.A.

8.3. DETERMINATION OF TOTAL IgG AND ALBUMIN CONCENTRATIONS

In addition to the CMV-specific IgG determinations, in each sample pair the total IgG concentration and the albumin concentration in serum and CSF have to be determined.

8.4. CALCULATION OF RESULTS/VALIDITY

* Validity

The validity criteria specified under 6.A. are applicable.

Repeat the run if the results do not meet the specification.

The following points also apply to CSF investigations:

- Only positive serum samples (antibody content > 0.65 AU) can be used for the calculation of the AI.
- Sera with an antibody content < 0.45 AU in a dilution of 1 : 200 are considered as sero-negative. In such cases the antibody index cannot be determined.
- The AI cannot be calculated in sero-negative and grey zone samples.
- In extremely rare cases, sero-negative patients may have demonstrable intrathecal antibodies because of CMV-associated encephalopathies.
- The assay range for CSF extends from 0.45 - 15 AU.
- CSF samples from patients with **positive serostatus** which at the dilution of 1 : 4 are below the assay range, should be re-tested in a lesser dilution (maximum 1 : 2).

* Evaluation

- Calculation of AU values see 6.A.
- Calculation of the pathogen-specific IgG quotient (Q_{spec}).

$$Q_{\text{spec}} = \frac{\text{AU CSF} \times \text{dilution CSF}}{\text{AU serum} \times \text{dilution serum}}$$

- Calculation of the antibody index

The pathogen-specific index is calculated from the fomulae:

1. $AI = Q_{\text{spec}}/Q_{\text{tot}}$ (for $Q_{\text{tot}} < Q_{\text{lim}}$)
2. $AI = Q_{\text{spec}}/Q_{\text{lim}}$ (for $Q_{\text{tot}} > Q_{\text{lim}}$)
3. $Q_{\text{lim}} = 0.93 \times \sqrt{Q_{\text{alb}}^2 + 6 \times 10^{-6}} - 1.7 \times 10^{-3}$

8.5. INTERPRETATION OF RESULTS

- * AI values from 0.6 - 1.3 are regarded as the normal range.
- * AI values > 1.3 and ≤ 1.5 are regarded as borderline.
- * **The abnormal range is defined as AI > 1.5 .**
- * AI values < 0.6 point to analytical errors and cannot be interpreted.
- * The CSF diagnostic criteria for an acute, active disease of the CNS are raised cell count and a raised albumin quotient. These reflect obstruction to CSF flow due to some inflammatory condition.
- * Elevated antibody indices are not reliable evidence of the acute phase of an infective CNS disease, because antibodies, even intrathecally, may persist for long periods, and because polyspecific CNS-intrinsic antibody synthesis may occur. In some circumstances it may be advisable to look for a significant change in the AI value by testing second paired serum-CSF samples. For this purpose further sample collection will be necessary and should be performed after a time interval chosen in the light of the clinical circumstances.

GENERAL HANDLING ADVICES

- * To avoid cross contamination do not exchange the vials and their screw caps.
- * The reagents have to be sealed immediately after use to avoid evaporation and microbial contamination.
- * After use, the reagents have to be stored as indicated to guarantee the shelf life.

- * After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

HEALTH AND SAFETY INFORMATION

- * The local occupational safety and health regulations have to be regarded.
- * Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents) should be handled as potentially infectious and used with all necessary precautions.

DISPOSAL CONSIDERATIONS

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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Flik, J., A comparison of 2 quantitative ELISAs used to detect IgG antibodies against HCMV in transplant recipients, Poster presented at the 5th CMV-Congress in Stockholm, (1995).

Hengster, P. and Dierich, M. P., Cytomegalovirus serology: Comparison of different ELISAs with complement fixation and indirect immunofluorescence for detection of antibodies. Lab. med. 12, 363-367, (1988).

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Pass, R. F. et al., Outcome of Symptomatic Congenital Cytomegalovirus Infection: Results of Longterm Longitudinal Follow up, Pediatrics 66, No. 5, 758-762, (1980).

Peterson, P. K. et al., Cytomegalovirus Disease in Renal Allograft Recipients: A Prospective Study of Clinical Features, Risk Factors and Impact on Renal Transplantation, Medicine 59, No. 4, 283-300, (1980).

Plotkin, S. A., Immunology of Cytomegalovirus, Comprehensive Immunology, Immunology of Human Infection, Plenum Publishing Corporation, 233 Spring Street, New York 89, 108, (1982).

Reynolds, D. W., Stagno, S. and Alford, C. A., in Laboratory Diagnosis of Cytomegalo Infections, Editors: Lennette, E. H. and Schmidt, N. J., American Public Health Association, 5th ed., Chapter 13, 399-439, (1979).

Sachers, M., Emmerich, P., Mohr, H. and Schmitz, H., Simple Detection of Antibodies to different Viruses using Rheumatoid Factor and Enzyme Labelled Antigen (ELA), J. Virol. Methods 10, 99-110, (1985).

Schmitz, H., von Deimling, U. and Flehmig, B., Detection of IgM Antibodies to Cytomegalovirus (CMV) Using an Enzyme Labelled Antigen (ELA), J. Gen. Virol. 50, 59-68, (1980).