Chlamydien-IgG-rELISA medac

Deutsch/English/Français

Chlamydien-IgG-rELISA medac

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Chlamydia-IgG-rELISA medac

Recombinant enzyme immunoassay for the quantitative detection of specific IgG antibodies to chlamydial LPS in human serum

Cat. no.: 480/TMB

FOR IN VITRO DIAGNOSTIC USE ONLY

INTRODUCTION

Chlamydiae are gram-negative bacterial pathogens. They have an obligate intracellular life cycle in mucosal surfaces, endothelial cells, smooth muscle cells, and according to recent findings in certain tissue structures of the central nervous system. Chlamydiae depend on energy-rich phosphates of their host cells and are therefore energy parasites.

The genus chlamydia comprises four species: *C. pneumoniae*, *C. trachomatis*, *C. psittaci*, and *C. pecorum*. *C. pneumoniae* and *C. trachomatis* are obligate pathogens of humans. *C. psittaci* is pathogenic for humans and a variety of animal species. To date, *C. pecorum* has been isolated from animals only.

Chlamydia trachomatis is one of the most frequent sexually transmitted pathogens worldwide. It causes infections of the urogenital tract and the eye. In most cases *C. trachomatis* infections are asymptomatic. This results in a variety of chronic diseases, which are sustained by the ascended and persisting agents. The clinical pictures include:

In women endometritis, adnexitis, periappendicitis, perihepatitis, and reactive arthritis. As a consequence of repeated adnexitis the tubes occlude which results in sterility.

In men the pathogens may ascend to the epididymis (→ epididymitis) and into the prostate gland (→ prostatitis) after incomplete or unsuccessful treatment of urethritis. A reduction of fertility has been discussed in these cases. Furthermore, post-urethritic reactive arthritis is known in men.

In newborns *C. trachomatis* can be transferred during delivery from the infected birth channel to the infant and can cause neonatal conjunctivitis and/or pneumonia.

Chlamydia pneumoniae infections occur worldwide. In addition to flu-like illness, the clinical picture includes sinusitis, pharyngitis, bronchitis, chronic obstructive pulmonary disease, atypical pneumonia, and reactive arthritis. A causal involvement in infection-conditioned asthma, sarcoidosis, lung cancer, atherosclerosis, acute myocardial infarction, stroke, multiple sclerosis, and late-onset of Alzheimer's disease has been discussed.

Chlamydia psittaci infects birds and mammals, which may transmit the pathogen to humans (ornithosis). The resultant frequently severe pneumonias may lead to life-threatening conditions if no targeted antibiotic intervention is started rapidly.

The usual laboratory methods for the detection of chlamydia infections include antigen and antibody determinations. The antigen detection (IFA, ELISA, PCR, LCR) is of value for the diagnosis of peripherally localized infections. In ascended courses of disease the detection of the pathogen is limited and has to be complemented by serology. Serological methods include complement fixation test, IFA, MIF, genus- and species-specific ELISA systems.

Chlamydiae contain as common immunodominant antigen the lipopolysaccharide (LPS), to which the first immune reaction is directed. The corresponding antibodies are already detectable within a few days after infection and, thus, allow early diagnosis. The use of paired sera enables the discrimination of current infections, reinfections, and reactivations (defined titer increases) from chronic persistent ones (constant antibody titers). Because of the limited persistence of the LPS antibodies after successful eradication of the pathogens, the diagnosis of current infections is obviously not influenced by past infections.

The Chlamydia-IgG-, IgA-, and IgM-rELISA medac are based upon a molecularly defined, genetechnologically produced antigen. It is an exclusively chlamydia-specific fragment from the LPS which has not been found in any other bacterial LPS.

The comparison of sera from blood donors with sera from patients with urogenital or respiratory infections has shown, that antibodies to chlamydia-LPS are detected more frequently in patients than in blood donors.

As the reaction is genus-specific, the test result will not differentiate between the chlamydia species. However, as the clinical picture already will hint to the suspected chlamydia species and additionally all the chlamydiae show an equal sensitivity towards specific antibiotics, the genus-specific reactivity does not hamper the consequences.

In addition to the

Chlamydia-IgG-rELISA medac

Cat. no.: **480/TMB** for 96 determinations,

we also distribute the following products:

Chlamydia-IgA-rELISA medac

Cat. no.: **490/TMB** for 96 determinations,

Chlamydia-IgM-rELISA medac

Cat. no.: **485/TMB** for 96 determinations

Chlamydia pneumoniae IgG sELISA medac

Cat. no.: **430/TMB** for 96 determinations,

Chlamydia pneumoniae IgA sELISA medac

Cat. no.: **431/TMB** for 96 determinations,

Chlamydia pneumoniae IgM sELISA medac

Cat. no.: **432/TMB** for 96 determinations,

Chlamydia trachomatis-IgG-pELISA medac

Cat. no.: **497/TMB** for 96 determinations,


Chlamydia trachomatis-IgA-pELISA medac

Cat. no.: **498/TMB** for 96 determinations,


cHSP60-IgG-ELISA medac

Cat. no.: **435** for 96 determinations.


TEST PRINCIPLE


The plate is coated with chlamydia-specific LPS fragment.

The chlamydia-specific LPS antibodies from the specimen bind to the antigen.

Peroxidase conjugated anti-human IgG antibodies bind to the IgG antibodies (P = peroxidase).

Incubation with TMB-substrate (*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

Advantages of the test

- ☞ Chemically defined, chlamydia-specific antigen.
- ☞ Non-infectious antigen material.
- ☞ The breakable microwell strips permit efficient use of the test.

KIT CONTENTS

Cat. no.: 480/TMB

1.

MTP

Microplate: 12 x 8 wells (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with chlamydia-specific, recombinant antigen and NBCS, ready to use.
2.

CONTROL	-
----------------	----------

Negative control: 2 vials with 0.75 ml each, human serum, ready to use, stained blue, contains NBCS, phenol, ProClin™ 300 and gentamycin sulfate.
3.

CONTROL	+
----------------	----------

Positive control: 2 vials with 0.75 ml each, human serum, ready to use, stained blue, contains BSA, phenol, ProClin™ 300 and gentamycin sulfate.
4.

WB

Wash buffer: 1 bottle with 100 ml PBS/Tween (10 x), pH 7.2 - 7.4, contains ProClin™ 300.
5.

BAC-DIL

Sample diluent: 1 bottle with 110 ml PBS/Tween/NBCS, pH 7.0 - 7.2, ready to use, stained blue, contains ProClin™ 300.
6.

CON

Conjugate: 4 vials with 5.0 ml each, goat anti-human IgG, HRP-conjugated, ready to use, stained green, contains BSA, phenol, ProClin™ 300 and gentamycin sulfate.
7.

TMB

TMB-substrate: 1 vial with 10 ml, ready to use.
8.

STOP

Stop solution: 2 vials with 14 ml each, 0.5 M sulfuric acid (H₂SO₄), ready to use.

1. STORAGE AND STABILITY

MATERIAL/REAGENT	STATE	STORAGE	STABILITY
Test kit	unopened	2 - 8 °C	until expiry date
Microplate	opened	2 - 8 °C in bag with desiccant	6 weeks
Controls	opened	2 - 8 °C	6 weeks
Wash buffer	diluted	2 - 8 °C	6 weeks
Sample diluent	opened	2 - 8 °C	6 weeks
Conjugate	opened	2 - 8 °C	6 weeks
TMB-substrate	opened	2 - 8 °C	6 weeks
Stop solution	opened	2 - 8 °C	until expiry date

Do not use the reagents after the expiry date.

2. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

- 2.1. Water for injection (H₂O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of wash buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37 °C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

3. PREPARATION OF THE REAGENTS

Before starting the test procedure all kit components must be equilibrated to room temperature (RT).

Calculate the number of wells required.

3.1. Microplate

After each removal of wells the aluminium bag has to be tightly resealed together with the desiccant. Storage and stability of the wells are indicated under point 1.

3.2. Wash buffer

Mix one volume of wash buffer (10 x) with nine volumes of water for injection (e.g. 50 ml wash buffer (10 x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

Crystals in the wash buffer (10 x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.

Do not mix test specific reagents (microplate, controls, conjugate) from different kit lots. In contrast to that, sample diluent, wash buffer, TMB-substrate and stop solution are generally exchangeable in all Chlamydia-and Mycoplasma-ELISA medac.

Reagents from other manufacturers must not be used in general.

Valid and reproducible results are only obtained if the test procedure is precisely followed.

4. SPECIMEN

- 4.1. The test is suitable for serum samples, but not for plasma.
- 4.2. Pretreatment of sera, e.g., inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain red blood cells.
- 4.3. Sera have to be diluted 1:100 with sample diluent.

5.A. TEST PROCEDURE

- 5.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

Microplate wells are ready for use and do not have to be prewashed.

- 5.2. Pipette 50 µl of sample diluent into the well A1 as blank (see 6.A.), in duplicate each 50 µl of the negative control, and each 50 µl of the positive control and the diluted patients' samples for single determinations.

If necessary, the microplate wells can be kept in a humid chamber up to 30 min at 2 - 8 °C before proceeding.

- 5.3. Incubate the microplate wells for 60 min (± 5 min) at 37 °C (± 1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.4. After incubation wash the microplate wells three times with each 200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

Do not allow the wells to dry out! Proceed immediately!

5.5. Add conjugate (coloured green) to each well.

50 µl of conjugate have to be pipetted into the wells if the test is done manually.

Please note:

When working with automated instruments, 60 µl of conjugate have to be pipetted into each well due to a higher evaporation in the incubation chambers of the automates.

The principal suitability of the test for automation could be shown during validation. Nevertheless we recommend to verify the compatibility with the automate employed.

5.6. Incubate the microplate wells again for 60 min (\pm 5 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil.

5.7. After incubation wash microplate wells again (see 5.4.).

5.8 Add 50 µl of TMB-substrate to each well and incubate the microplate wells for 30 min (\pm 2 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.

5.9. Stop the reaction by adding 100 µl of stop solution to each well. Positive samples turn yellow.

Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.

The reading should be done within 15 min after adding the stop solution!

5.B. TABLE FOR THE TEST PROCEDURE

	Blank (A1)	Negative control	Positive control	Sample
Sample diluent	50 µl	-	-	-
Negative control	-	50 µl	-	-
Positive control	-	-	50 µl	-
Sample	-	-	-	50 µl
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer.				
Conjugate	50/60 µl *)	50/60 µl *)	50/60 µl *)	50/60 µl *)
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer.				
TMB-substrate	50 µl	50 µl	50 µl	50 µl
Incubate for 30 min at 37 °C in the dark.				
Stop solution	100 µl	100 µl	100 µl	100 µl
Photometric reading at 450 nm (ref. 620 - 650 nm)				

*) manual/automatic procedure (see 5.5.)

6.A. CALCULATION OF RESULTS (VALIDITY)

- * Read OD values at 450 nm (reference wavelength 620 - 650 nm).
- * Subtract the OD value of the blank (well A1) from all other OD values.
- * The OD value of the blank has to be **< 0.100**.
- * The mean OD value of the **negative control** has to be **< 0.200**.
- * The OD value of the **positive control** has to be **> 0.800**.
- * **Cut-off = mean OD value of the negative control + 0.320**
- * **Grey zone = Cut-off ± 10%**

6.B. EVALUATION OF RESULTS:

6.B.1. QUALITATIVE

Result	Valuation
OD < Grey zone	negative
OD of Cut-off +/- 10%	equivocal
OD > Grey zone	positive

6.B.2. QUANTITATIVE

The endtiter (= last positive or borderline dilution) of each sample can be calculated.

Firstly, the cut-off index has to be calculated.

$$\text{Cut-off Index} = \frac{\text{OD Sample}}{\text{Cut-off}}$$

Secondly, the calculation of the endtiter is possible in two ways:

1. **1:Endtiter = 111 x (OD/Cut-off)**

The result has to be rounded down to the next lower whole titer step, i.e., to 1:100, 1:200, 1:400, etc.

This calculation is valid for OD values ≤ 2.0 . Samples with OD values > 2.0 should be retested employing a higher starting dilution (e.g., 1:400).

If a higher starting dilution than 1:100 has been used the above mentioned formula has to be multiplied with the corresponding dilution factor.

Example: Employed sample dilution 1:400 means an fourfold higher starting dilution.

OD = 1.78
Cut-off = 0.42
1:Endtiter = 111 x (1.78/0.42) x 4 = 1882 \Rightarrow 1:1600

2. The endtiter which corresponds to the cut-off index can also be determined from the table below.

Cut-off Index	IgG	
	Result	Titer
< 0.9	Negative	< 1:100
0.9 - ≤ 1.1	Equivocal	
> 1.1 - ≤ 1.8	Positive	1:100
> 1.8 - ≤ 3.6	Positive	1:200
> 3.6 - ≤ 7.2	Positive	1:400

This calculation is only valid for OD values ≤ 2.0; otherwise see above mentioned example.

6.C. INTERPRETATION OF THE RESULTS

6.C.1. CRITERIA FOR SIGNIFICANT TITER INCREASES

For the diagnosis of current, fresh infections, reinfections, reactivations (significant titer increases) and dissociation from persistent infections (constant antibody titers) we principally recommend the use of paired sera. They should be obtained 10 - 14 days apart.

The following criteria have been described as a possibility for the assessment of significant titer increases (Persson et al., Verkooyen et al.):

Threefold or greater increase in
 specific IgG **or** IgA antibody titers
 or
 twofold or greater increase in
 specific IgG **and** IgA antibody titers
 or
 twofold or greater increase in
 specific IgM antibody titers

6.C.2. SPECIFIC IgG and IgA INTERPRETATION

Possible Results		Interpretation
IgG	IgA	
+	+	1. IgG and IgA positive; indication of active infection. Further judgement dependent on symptoms and titer course.
+	-	2. Only IgG positive; indication of a past infection. In case of clinical suspicion determine a threefold IgG titer increase between two serum samples (interval 10 - 14 days) and retest for IgA.
+	+/-	3. IgG positive, IgA grey zone; possibility of a fresh or a reactivated or subsiding infection; retest IgG and IgA after 10 - 14 days.
+/-	-	4. Only IgG in grey zone; past infection cannot be excluded. In case of clinical suspicion retest IgG and IgA after 10 - 14 days.
-	+	5. Only IgA positive; possibility of an early stage of active infection or solitary persisting IgA (see p.29). Retest IgA and IgG after 10 - 14 days.
-	+/-	6. Only IgA in the grey zone; possibility of a very early stage of active infection; retest IgG and IgA after 10 - 14 days.
+/-	+	7. IgG in the grey zone, IgA positive; possibility of an early stage of active infection. Retest IgG and IgA after 10 - 14 days.
-	-	8. IgG and IgA negative; no serological indication of infection. In case of justified clinical suspicion perform direct antigen detection and retest IgG and IgA after 10 - 14 days.

This table should be used as guideline only.

- * The IgG results should always be interpreted in connection with IgA and IgM, the clinical picture and additional diagnostic parameters.
- * In individual cases solitary IgA antibodies may persist. This immunological phenomenon occurs in various bacterial and viral infections. A clinical relevance cannot be assessed.
- * Samples with OD values within the grey zone should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.
- * High concentrations of hemoglobin and of bilirubin in serum do not have an influence on the results.
- * Scarcely, high lipid concentrations in serum may influence the test results.
- * Cross-reactivities with antibodies to parvovirus-B19 may not be excluded in individual cases.

Comment:

In cases of fresh acute chlamydial infections the serological antibody results may be negative despite clinics and positive antigen detection. If a serological confirmation of a positive antigen result or if a follow-up is desired we recommend to test after 10 - 14 days for seroconversion.

7. PERFORMANCE CHARACTERISTICS

The following performance characteristics were determined during the diagnostic evaluation.

7.A. PREVALENCE

Sera from various patients cohorts (culture-positive, -negative STD patients, prostitutes, infertile females, infertile males, patients with rheumatic diseases, patients with chronic obstructive respiratory diseases [COPD]), and controls (pregnant women, 2 blood donor cohorts) were investigated for the **determination of LPS antibody prevalence.**

Cohorts	Prevalence		
	IgG	IgA	IgM
Culture-positive STD patients (medac-internal investigations)	67% (88/132)	56% (74/132)	14% (18/132)
Culture-negative STD patients (medac-internal investigations)	33% (41/125)	13% (16/125)	2% (3/125)
Prostitutes (Schmitz et al.)	86% (255/295)	47% (141/295)	25% (63/295)
Infertile females (Schmitz et al.)	75% (62/83)	45% (37/83)	10% (8/83)
Infertile males (Schmitz et al.)	71% (144/203)	35% (71/203)	9% (18/203)
Patients with rheumatic diseases (medac-internal investigations)	83% (158/191)	63% (120/191)	7% (13/191)
COPD patients (Verkooyen et al. 1997)	53% (144/271)	32% (88/271)	3% (9/271)
Pregnant women (medac-internal investigations)	32% (62/192)	17% (32/192)	8% (16/192)
Blood donor group 1 (Verkooyen et al. 1998)	29% (325/1104)	n.d*	n.d*
Blood donor group 2 (medac-internal investigations)	37% (154/416)	13% (54/416)	3% (6/240)

*: n.d = not determined

7.B. PRECISION

Sample	Intra-assay variation				Sample	Interassay variation (n = 11)		
	Ø OD	SD	CV (%)	n		Ø OD	SD	CV (%)
NC	0.123	0.010	8	24	NC	0.109	0.015	14
BC	0.431	0.021	5	24	BC	0.431	0.027	6
PC	1.419	0.050	4	24	PC	1.417	0.071	5
N° 1	0.398	0.012	3	23	N° 3	0.567	0.034	6
N° 2	2.224	0.094	4	24	N° 4	1.960	0.086	4

NC = negative control; BC = weak positive control (not included in the kit);
PC = positive control

GENERAL HANDLING ADVICES

- * Do not exchange the vials and their screw caps in order to avoid cross-contamination.
- * The reagents have to be sealed immediately after use in order to avoid evaporation and microbial contamination.
- * After use, the reagents have to be stored as indicated to guarantee the shelf life.
- * After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

HEALTH AND SAFETY INFORMATION

- * The local occupational safety and health regulations have to be regarded.
- * Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents), should be handled as potentially infectious and used with all necessary precautions.

DISPOSAL CONSIDERATIONS

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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Brade, L., Holst, O., Kosma, P., Zhang, Y.-X., Paulsen, H., Krausse, R., Brade, H.: Characterization of murine monoclonal and murine, rabbit, and human polyclonal antibodies against chlamydial lipopolysaccharide. *Infect. Immun.* 58, 205-213 (1990).

Brade, L., Brunnemann, H., Ernst, M., Fu, Y., Kosma, P., Näher, H., Persson, K., Brade, H.: Occurrence of antibodies against chlamydial lipopolysaccharide in human sera as measured by ELISA using an artificial glycoconjugate antigen. *FEMS Immunol. Med. Microbiol.* 8, 27-42 (1994).

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Falck, G., Gnarpe, J., Gnarpe, H.: Persistent *Chlamydia pneumoniae* infections in a Swedish family. *Scand. J. Infect. Dis.* 28, 271-273 (1996)

Gérard, H.C., Schumacher R.H., El-Gabalawi, H., Goldbach-Mansky, R., Hudson, A.P.: ***Chlamydia pneumoniae*** present in the human synovium are viable and metabolically active. *Microbiol. Pathol.* 29, 17-24 (2000).

Grayston, J.T., Wang, S.P., Kuo, C.C., Campbell L.A.: Current knowledge of ***Chlamydia pneumoniae***, strain **TWAR**, an important cause of pneumonia and other respiratory diseases. *Eur. J. Clin. Microbiol. Infect. Dis.* 8, 191-202 (1989)

Gupta, S., Leatham, E.W., Carrington, D., Mendall, M.A., Kaski, J.C., Camm, A.J.: Elevated ***Chlamydia pneumoniae*** antibodies, cardiovascular events, and azithromycin in male survivors of myocardial infarction. *Circulation* 96, 404-407 (1997).

Hahn, D.L., Peeling, R.W., Dillon, E., McDonald, R., Saikku, P.: Serologic markers for ***C. pneumoniae*** in asthma. *Ann. Allergy Asthma Immunol.* 84, 227-233 (2000).

- Köhler, M., Jendro, C.: Bedeutung der Persistenz von ***Chlamydia trachomatis*** im Gelenk für die Pathogenese der Chlamydien-induzierten Arthritis unter Berücksichtigung der therapeutischen Implikationen. Akt. Rheumatol. 22, 170-175 (1997).
- Land, J.A., Evers, J.L., Goossens, J.V.: How to use Chlamydia antibody testing in subfertility patients. Hum. Reprod. 13, 1094-1098 (1998).
- Maass, M., Bartels, C., Engel, P., Mamat, U., Sievers, H.H.: Endovascular presence of viable ***Chlamydia pneumoniae*** is a common phenomenon in coronary artery disease. JACC 31, 827-32 (1997).
- Morré, S.A., De Groot, C.J., Killestein, J., Meijer, C.J., Polman, C.H., Van der Falk, P., Van den Brule, A.J.: Is ***Chlamydia pneumoniae*** present in the central nervous system of multiple sclerosis patients? Ann. Neurol. 48, 399 (2000).
- Paavonen, J.: ***Chlamydia trachomatis***: A major cause of mucopurulent cervicitis and pelvic inflammatory disease in women. In: Sexually Transmitted Diseases. Advances in Diagnosis and Treatment. Curr. Probl. Dermatol. Elsner, P., Eichmann, A. (eds.). Basel, Karger, 24, 110-122 (1996).
- Patton, D., Askienazy-Elbhar, M., Henry-Suchet, J., Campbell, L.A., Cappucio, A., Tannous, W., Wang, S.P., Kuo, C.C.: Detection of ***Chlamydia trachomatis*** in fallopian tube tissue in women with post-infectious tubal infertility. Am. J. Obstet. Gynecol. 171, 95-101 (1994).
- Persson, K., Haidl, S.: Evaluation of a commercial test for antibodies to the chlamydial lipopolysaccharide (Medac™) for serodiagnosis of acute infections by ***Chlamydia pneumoniae (TWAR)*** and ***Chlamydia psittaci***. APMIS 108, 131-138 (2000).
- Petersen, E.E., Clad, A.: Genitale Chlamydieninfektionen. Deutsches Ärzteblatt 92, A-277-282 (1995).
- Saikku, P.: Diagnosis of acute and chronic ***Chlamydia pneumoniae*** infections. In: Orfila, J. et al. (eds). Chlamydial Infections. Proc. Eighth Int. Symp. on Human Chlamydial Infect. Società Editrice Esculapio, Bologna, Italy, p 163-170 (1994).
- Schachter, J.: The intracellular life of Chlamydia. Curr. Top. Microbiol. Immun. 138, 109-39 (1988).

Schmitz, F.J., Küppers, B., Reinartz, R., Vossel, R., Schuppe, D., Bielfeld, D., Heinz, H.P.: Prevalence of ***Chlamydia trachomatis*** antigen and Chlamydia antibodies in prostitutes, infertile female and infertile male patients - Comparison of different methods. In: Stary, A. (ed.) Proc. Europ. Soc. Chlamydia Res. Società Editrice Esculapio, Bologna, Italy, p. 339 (1996).

Sriram, S., Stratton, C.W., Yao S.: ***Chlamydia pneumoniae*** infections in the central nervous system in multiple sclerosis. Ann. Neurol. 46, 6-14 (1999).

Verkooyen, R.P., van Lent, N.A., Mousavi Joulandan, S.A., Snijder, R.J., van den Bosch, J.M., van Helden, H.P., Verbrugh, H.A.: Diagnosis of ***Chlamydia pneumoniae*** infection in patients with chronic obstructive pulmonary disease by micro-immunofluorescence and ELISA. J. Med. Microbiol. 46, 959-964 (1997).

Verkooyen, R.P., Willemsse, D., Hiep-van Casteren S.C.A.M., Mousavi Joulandan, S.A., Snijder, R.J., van den Bosch, J.M.M., van Helden, H.P.T., Peeters, M.F., Verbrugh, H.A.: Evaluation of PCR, culture, and serology for diagnosis of ***Chlamydia pneumoniae*** respiratory infections. J. Clin. Microbiol., 36, 2301-2307 (1998).

Ward, M.E.: The immunobiology and immunopathology of chlamydial infections. APMIS 103, 769-796 (1995).

Weström, L.V.: Chlamydia and its effect on reproduction. J. Brit. Fertil. Soc. 1, 23-30 (1996).

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