Chlamydia pneumoniae-IgA-sELISA medac

Deutsch/English/Français

Chlamydia pneumoniae-IgA-sELISA medac

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Chlamydia pneumoniae-IgA-sELISA medac

Enzyme immunoassay for the detection of IgA antibodies
to ***Chlamydia pneumoniae***

Cat.no.: 431/TMB

FOR IN VITRO DIAGNOSTIC USE ONLY

INTRODUCTION

Chlamydiae are gram-negative bacterial pathogens. They have an obligate intracellular life cycle in mucosal surfaces, endothelial cells, smooth muscle cells, and according to recent findings in certain tissue structures of the central nervous system. Chlamydiae depend on energy-rich phosphates of their host cells and are therefore energy parasites.

The genus chlamydia comprises four species: ***C. pneumoniae***, ***C. trachomatis***, ***C. psittaci***, and ***C. pecorum***. ***C. pneumoniae*** and ***C. trachomatis*** are obligate pathogens of humans. ***C. psittaci*** is pathogenic for humans and a variety of animal species. To date, ***C. pecorum*** has been isolated from animals only.

Infections with ***C. pneumoniae*** occur worldwide. The spectrum of diseases, in addition to flu-like illness, includes sinusitis, pharyngitis, bronchitis, chronic obstructive pulmonary diseases, pneumonia, and reactive arthritis. A causal involvement of ***C. pneumoniae*** in infectious asthma, sarcoidosis, lung cancer, atherosclerosis, acute myocardial infarction, brain stroke, multiple sclerosis, and late onset of Alzheimer's disease remains an area of actual investigation.

According to Grayston and Saikku (1989), who firstly had described this chlamydia species, almost everybody is infected and reinfected with ***C. pneumoniae*** throughout his life. The mainly weak and/or diffuse symptomatic course of ***C. pneumoniae*** infections makes their detection more difficult; undetected infections may lead to chronic courses of disease with serious sequelae.

The prevalence of ***C. pneumoniae*** infections has been estimated at less than 10% in children not yet attending school, followed by a strong increase to 50% beyond the age of 20 years and up to 80-100% beyond the age of 70 years.

The diagnosis of ***C. pneumoniae*** infections is based upon isolation of the pathogen from cell culture, direct antigen detection, nucleic acid amplification tests, and serology.

Isolation of the pathogen from cell culture normally needs various subsequent passages; it is time consuming, restricted to special laboratories, and only in few cases successful. Direct immunofluorescence assays (IFA) and antigen enzyme immunoassays (EIA) have not been widely used; they have been reported to be of low sensitivity and also specificity. A commercially PCR or LCR, to date, are not available.

The shortcomings of culture and non-culture methods for the diagnosis of ***C. pneumoniae*** infections have made serology the method of choice.

For the detection of species-specific antibodies the microimmunofluorescence (MIF) has been considered as the "gold standard". MIF is laborious, subjective, and requires much experience. This test system is not standardized; the use of different antigens and different cut off criteria for past and recent or current infections creates significant lab-to-lab variations of results.

The ***Chlamydia pneumoniae*-sELISA medac** employs a highly purified and specific antigen. The detection of IgM, IgA, and IgG antibodies allows the appraisal of the infection's status and the follow up after treatment.


The ***Chlamydia pneumoniae*-sELISA medac** fulfills the needs for standardization, objectivity, reproducibility, and automation in the routine laboratory.

In addition to the **Chlamydia pneumoniae-IgA-sELISA medac**
Cat. no.: **431/TMB** for 96 determinations


we also distribute the following products:

Chlamydia pneumoniae-IgG-sELISA medac
Cat. no.: **430/TMB** for 96 determinations,
Chlamydia pneumoniae-IgM-sELISA medac
Cat. no.: **432/TMB** for 96 determinations,
Chlamydia trachomatis-IgG-pELISA medac
Cat. no.: **497/TMB** for 96 determinations,
Chlamydia trachomatis-IgA-pELISA medac
Cat. no.: **498/TMB** for 96 determinations,
Chlamydia-IgG-rELISA medac
Cat. no.: **480/TMB** for 96 determinations,
Chlamydia-IgA-rELISA medac
Cat. no.: **490/TMB** for 96 determinations,
Chlamydia-IgM-rELISA medac
Cat. no.: **485/TMB** for 96 determinations,
cHSP60-IgG-ELISA medac
Cat. no.: **435** for 96 determinations.


TEST PRINCIPLE


The plate is coated with a highly purified *C. pneumoniae*-specific antigen preparation.

The *C. pneumoniae*-specific antibodies from the specimen bind to the antigen.

Peroxidase-conjugated anti-human IgA antibodies bind to the IgA antibodies (P = peroxidase).

Incubation with TMB-substrate (*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

Advantages of the test

- ☞ High sensitivity and specificity.
- ☞ The breakable microwell strips permit efficient use of the test.

KIT CONTENTS

Cat. no.: 431/TMB

1.

MTP

Microplate: 12 x 8 wells, pink-coded (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with **C. pneumoniae**-specific antigen and FCS, ready for use.
2.

CONTROL	-
----------------	----------

Negative control: 2 vials with 0.75 ml each, human serum, ready for use, stained blue, contains NBCS, phenol, ProClin™ 300 and gentamycin sulfate.
3.

CONTROL	+
----------------	----------

Positive control: 2 vials with 0.75 ml each, human serum, ready for use, stained blue, contains BSA, phenol, ProClin™ 300 and gentamycin sulfate.
4.

WB

Wash buffer: 1 bottle with 100 ml PBS/Tween (10 x), pH 7.2-7.4, contains ProClin™ 300.
5.

BAC-DIL

Sample diluent: 1 bottle with 110 ml PBS/Tween/NBCS, pH 7.0-7.2, ready for use, stained blue, contains ProClin™ 300.
6.

CON

Conjugate: 4 vials with 5.0 ml each, goat anti-human IgA, HRP-conjugated, ready for use, stained yellow, contains BSA, phenol, ProClin™ 300 and gentamycin sulfate.
7.

TMB

TMB-substrate: 1 vial with 10 ml, ready for use.
8.

STOP

Stop solution: 2 vials with 14 ml each, 0.5 M sulfuric acid (H₂SO₄), ready for use.

1. STORAGE AND STABILITY

MATERIAL/REAGENT	STATE	STORAGE	STABILITY
Test kit	unopened	2 - 8 °C	Until expiry date
Microplate	opened	2 - 8 °C in bag with desiccant	12 weeks
Controls	opened	2 - 8 °C	12 weeks
Wash buffer	diluted	2 - 8 °C	12 weeks
Sample diluent	opened	2 - 8 °C	12 weeks
Conjugate	opened	2 - 8 °C	12 weeks
TMB-substrate	opened	2 - 8 °C	12 weeks
Stop solution	opened	2 - 8 °C	Until expiry date

Do not use the reagents after the expiry date.

2. REAGENTS AND MATERIALS REQUIRED BUT NOT SUPPLIED

- 2.1. Water for injection (H₂O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of wash buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37 °C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

3. PREPARATION OF THE REAGENTS

Before starting the test procedure all kit components must be equilibrated to room temperature (RT).

Calculate the number of wells required.

3.1. Microplate

After each removal of wells the aluminium bag has to be tightly resealed together with the desiccant. Storage and stability of the wells are indicated under point 1.

3.2. Wash buffer

Mix one volume of wash buffer (10 x) with nine volumes of water for injection (e.g. 50 ml wash buffer (10 x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

Crystals in the wash buffer (10 x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.

Do not mix test specific reagents (microplate, controls, conjugate) from different kit lots. In contrast to that, sample diluent, wash buffer, TMB-substrate and stop solution are generally exchangeable in all Chlamydia-and Mycoplasma-ELISA medac.

Reagents from other manufacturers must not be used in general.

Valid and reproducible results are only obtained if the test procedure is precisely followed.

4. SPECIMEN

- 4.1. The test is suitable for serum samples, but not for plasma.
- 4.2. Pretreatment of sera, e.g. inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain red blood cells.
- 4.3. Sera have to be diluted 1:50 with sample diluent.

5.A. TEST PROCEDURE

- 5.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

Microplate wells are ready for use and do not have to be prewashed.

- 5.2. Pipette 50 µl of sample diluent into the well A1 as blank (see 6.A.), in duplicate each 50 µl of the negative control, and each 50 µl of the positive control and the diluted patients' samples for single determinations.

If necessary, the microplate wells can be kept in a humid chamber up to 30 min at RT before proceeding.

- 5.3. Incubate the microplate wells for 60 min (\pm 5 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.4. After incubation wash the microplate wells three times with each 200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

Do not allow the wells to dry out! Proceed immediately!

5.5. Add conjugate (coloured yellow) to each well.

50 µl of conjugate have to be pipetted into the wells if the test is done manually.

Please note:

When working with automated instruments, 60 µl of conjugate have to be pipetted into each well because of a higher evaporation in the incubation chambers of the automates.

The suitability of the test for automated devices was confirmed during the evaluation of the test. Nevertheless we recommend to verify the compatibility of the test with the devices used in the lab.

5.6. Incubate the microplate wells again for 60 min (\pm 5 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil.

5.7. After incubation wash microplate wells again (see 5.4.).

5.8. Add 50 µl of TMB-substrate to each well and incubate the microplate wells for 30 min (\pm 2 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.

5.9. Stop the reaction by adding 100 µl of stop solution to each well. Positive samples turn yellow.

Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.

The reading should be done within 15 min after adding the stop solution!

5.B. TABLE FOR THE TEST PROCEDURE

	Blank (A1)	Negative control	Positive control	Sample
Sample diluent	50 µl	-	-	-
Negative control	-	50 µl	-	-
Positive control	-	-	50 µl	-
Sample	-	-	-	50 µl
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer.				
Conjugate	50/60 µl *)	50/60 µl *)	50/60 µl *)	50/60 µl *)
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer.				
TMB-substrate	50 µl	50 µl	50 µl	50 µl
Incubate for 30 min at 37 °C in the dark.				
Stop solution	100 µl	100 µl	100 µl	100 µl
Photometric reading at 450 nm (ref. 620 - 650 nm)				

*) manual/automatic procedure (see 5.5.)

6.A. CALCULATION OF RESULTS (VALIDITY)

- * Read OD values at 450 nm (reference wavelength 620 - 650 nm).
- * Subtract the OD value of the blank (well A1) from all other OD values.
- * The OD value of the blank has to be **< 0.100**.
- * The mean OD value of the **negative control** has to be **< 0.100**.
- * The OD value of the **positive control** has to be **> 0.800**.
- * **Cut-off = mean OD value of the negative control + 0.390**
- * **Grey zone = Cut-off ± 10%**

6.B. EVALUATION OF RESULTS

6.B.1. QUALITATIVE

Result	Valuation
OD < Grey zone	negative
OD of Cut-off +/- 10%	equivocal
OD > Grey zone	positive

6.B.2. QUANTITATIVE

Cut-off-Index: $\frac{\text{OD Sample}}{\text{OD Cut-off}}$	Valuation
< 0.9	negative
0.9 - 1.1	equivocal
> 1.1	positive

- * Samples with OD values within the grey zone should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.
- * The results should always be interpreted in connection with IgG and IgM and with the clinical data and additional diagnostic parameters.
- * High concentrations of hemoglobin, of bilirubin and of lipids in serum do not have an influence on the results.
- * Cross-reactivities with antinuclear antibodies, heterophilic antibodies and antibodies to *C. psittaci* as well as to *C. trachomatis* cannot be excluded in individual cases.

6.C. SPECIFIC IgG/IgA INTERPRETATION

Possible Results		Interpretation
IgG	IgA	
+	+	1. IgG and IgA positive; indication of infection. Further judgement dependent on symptoms and titer course.
+	-	2. Only IgG positive; indication of a past infection. In case of clinical suspicion determine a fourfold IgG titer increase between two serum samples and retest for IgA.
+	+/-	3. IgG positive, IgA grey zone; possibility of a beginning or a subsiding infection; retest IgA after 10 - 14 days.
+/-	-	4. Only IgG in grey zone; past infection cannot be excluded. In case of clinical suspicion retest IgG and IgA after 10 - 14 days.
-	+	5. Only IgA positive; possibility of an early stage of infection or solitary, persisting IgA; retest IgA and IgG after 10 - 14 days.
-	+/-	6. Only IgA in the grey zone; possibility of a early stage of infection; retest IgA and IgG after 10 - 14 days.
+/-	+	7. IgG in the grey zone, IgA positive; possibility of an early stage of infection; retest IgA and IgG after 10 - 14 days.
-	-	8. IgG and IgA negative; no indication of infection. In case of justified clinical suspicion perform direct antigen detection and retest IgA and IgG after 10 - 14 days.

Comment:

In cases of fresh acute chlamydial infections the serological antibody results may be negative despite clinics and positive antigen detection. If a serological confirmation of a positive antigen result or if a follow-up is desired we recommend to test after 10 - 14 days for seroconversion.

7. PERFORMANCE CHARACTERISTICS

We determined the following performance characteristics during the diagnostic evaluation.

7.A. SPECIFICITY AND SENSITIVITY

Sera from patients without symptoms of respiratory infection were investigated for the **determination of specificity**.

By MIF no antibodies to *C. pneumoniae* were detectable.

Sera from patients suspected for respiratory infection were investigated for the **determination of sensitivity**.

By MIF antibodies to *C. pneumoniae* were detectable in all sera.

Patient Group	Specificity	
	IgA	IgG
Patients without respiratory infections; no antibodies to <i>C. pneumoniae</i> by MIF.	93% (n=86)	95% (n=42)

Patient Group	Sensitivity	
	IgA	IgG
Patients with respiratory infections; all sera with antibodies to <i>C. pneumoniae</i> by MIF.	95% (n=74)	99% (n=117)

7.B. PRECISION

Sample	Intra-Assay Variation				Sample	Inter-Assay Variation (n = 11)		
	mean OD	SD	CV (%)	n		mean OD	SD	CV (%)
NC	0.026	0.011	42	21	NC	0.027	0.007	28
BC	0.538	0.022	4	21	BC	0.533	0.045	8
PC	1.251	0.037	3	21	PC	1.243	0.110	9
N° 1	0.169	0.012	7	21	N° 3	0.663	0.048	7
N° 2	1.419	0.057	4	21	N° 4	0.975	0.092	9

NC = negative control; BC = weak positive control (not included in the kit);
PC = positive control

GENERAL HANDLING ADVICES

- * To avoid cross-contamination do not exchange the vials and their screw caps.
- * The reagents have to be sealed immediately after use to avoid evaporation and microbial contamination.
- * After use, the reagents have to be stored as indicated to guarantee the shelf life.
- * After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

HEALTH AND SAFETY INFORMATION

- * The local occupational safety and health regulations have to be regarded.
- * Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents), should be handled as potentially infectious and used with all necessary precautions.

DISPOSAL CONSIDERATIONS

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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Danesh, J., Collins., Peto, R.: Chronic infections and coronary heart disease: is there a link? *Lancet* 350, 430-436 (1997).

Elkind, M.S., Lin, I.F., Grayston, J.T., Sacco, R.L.: ***Chlamydia pneumoniae*** and the risk of first ischemic stroke. The Northern Manhattan stroke study. *Stroke* 31, 1521-1525 (2000).

Gérard, H.C., Schumacher, R.H., El-Gabalawy, H., Goldbach-Mankys, R., Hudson, A.P.: ***Chlamydia pneumoniae*** present in the human synovium are viable and metabolically active. *Microbiol. Pathol.* 29, 17-24 (2000).

Grayston, J.T.: Epidemiology of ***Chlamydia pneumoniae*** (TWAR). In: *Chlamydia Research*. Angelika Stary (ed.). Proceedings of the Third Meeting of the European Society for Chlamydia Research, Vienna, Austria, 11.-14. September, 211-214 (1996).

Grayston, J.T., Aldous, M.B., Easton, A., Wang, S.P., Kuo C.C., Campbell, L.A. Altman, J.: Evidence that ***Chlamydia pneumoniae*** causes pneumonia and bronchitis. *J. Infect. Dis.* 168, 1231-1235 (1993).

Gupta, S., Leatham, E.W., Carrington, D., Mendall, M.A., Kaski, J.C., Camm, A.J.: Elevated ***Chlamydia pneumoniae*** antibodies, cardiovascular events, and azithromycin in male survivors of myocardial infarction. *Circulation* 99, 401-407 (1997).

Gurfinkel, E., Bozovich, G., Caroca, A., Beck, E, Mautner, B.: Randomised trial of roxithromycin in nin-Q-wave coronary syndromes: ROXIS pilot study. *Lancet* 350, 404-407 (1997).

Hahn, D.L., Peeling, R.W., Dillon, E., McDonald, R., Saikku, P.: Serologic Markers for ***C. pneumoniae*** in asthma. *Ann. Allergy Asthma Immunol.* 84, 227-233 (2000).

Heinzl, S.: Chlamydien, Arteriosklerose, Asthma und MS. *Med. Monatsschr. Pharm.* 23, 37 (2000).

Hunter, S.F., Hafler, D.A.: Ubiquitous pathogens: links between infection and autoimmunity in MS? *Neurology* 55, 164-165 (2000).

Kuo, C.-C., Jackson, L.A., Campbell, L.A., Grayston, J.T.: ***Chlamydia pneumoniae***(TWAR). Clin. Microbiol. Rev. 8, 451-461 (1995).

Layh-Schmitt, G., Bendl, C., Hildt, U., Dong-Si, T., Juttler, E., Schnitzler, P., Grond-Ginsbach, C., Grau, A.J.: Evidence for infection with ***Chlamydia pneumoniae*** in a subgroup of patients with multiple sclerosis. Ann. Neurol. 47, 652-655 (2000).

Maass, M., Bartels, C., Engel, P.M., Mamat, U., Sievers, H.H.: Endovascular presence of viable ***Chlamydia pneumoniae*** is a common phenomenon in coronary artery disease. J. Am. Coll. Cardiol. 31, 827-832 (1998).

Morré, S.A., De Groot, C.J., Killestein, J., Meijer, C.J., Polman, C.H., Van der Valk, P., Van den Brule, A.J.: Is ***Chlamydia pneumoniae*** present in the central nervous system of multiple sclerosis patients? Ann. Neurol. 48, 399 (2000).

Muhlestein, J.B.: The link between ***Chlamydia pneumoniae*** and atherosclerosis. Infect. Med. 14, 380-382, 392, 426 (1997).

Saikku, P.: Chronic ***Chlamydia pneumoniae*** infections. In: Chlamydia Research. Angelika Stary (ed.). Proceedings of the Third Meeting of the European Society for Chlamydia Research, Vienna, Austria, 11.-14. September. 215-218 (1996).

Saikku, P., Leinonen, M., Mattila, K., Ekman, M.R., Nieminen, M.S., Mäkela, P.H., Huttunen, J.K., Valtonen, V.: Serological evidence of an association of a novel Chlamydia, TWAR, with chronic coronary heart disease and acute myocardial infarction. Lancet 2, 983-986 (1988).

Saikku, P., Leinonen, M., Tenkanen, L., Linnanmäki, E., Ekman, M.R., Manninen, V., Manttari, M., Frick, M.H., Huttunen, J.K.: Chronic ***Chlamydia pneumoniae*** infection as a risk factor for coronary heart disease in the Helsinki Heart Study. Ann. Int. Med. 116, 272-278 (1992).

Samra, Z., Soffer, Y.: IgA antichlamydial antibodies as a diagnostic tool for monitoring of active chlamydial infection. Eur. J. Epidemiol. 8, 882-884 (1992).

Schumacher, H.R.: Chlamydial Arthritis. In: Chlamydia Research. Pekka Saikku (ed.). Proceedings of the 4th Meeting of the European Society for Chlamydia Research, Helsinki, Finland, 20.-23. August. 229 (2000).

Sriram, S., Stratton, C.W., Yao, S.: ***Chlamydia pneumoniae*** infection in the central nervous system in multiple sclerosis. Ann. Neurol. 46, 6-14 (1999).

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