Chlamydia trachomatis-IgG-pELISA medac

Deutsch/English/Français

Chlamydia trachomatis-IgG-pELISA medac

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Chlamydia trachomatis-IgG-pELISA medac

Enzyme immunoassay for the detection of IgG antibodies
to ***Chlamydia trachomatis***

Cat. no.: 497/TMB

FOR IN VITRO DIAGNOSTIC USE ONLY

INTRODUCTION

Chlamydiae are gram-negative bacterial pathogens. They have an obligate intracellular life cycle in mucosal surfaces and according to recent findings in endothelial cells and smooth muscle cells. Chlamydiae depend on energy-rich phosphates of their host cells and are therefore called energy parasites.

The genus Chlamydia comprises four species: ***C. pneumoniae***, ***C. psittaci***, ***C. pecorum*** and ***C. trachomatis***.

C. pneumoniae and ***C. trachomatis*** are obligate pathogens of humans. ***C. psittaci*** is pathogenic for humans and a variety of animal species. Up to now ***C. pecorum*** has been isolated from animals only.

C. trachomatis is one of the most frequent agents causing sexually transmitted infections of the urogenital tract and the eye. During delivery ***C. trachomatis*** can be transferred to the infant and can cause neonatal conjunctivitis and pneumonia.

C. trachomatis is subdivided into 18 serovars causing different clinical manifestations. Serovars A-C are responsible for the endemic trachoma, an ocular disease, which can lead to blindness. It affects hundreds of million people in developing countries. Lymphogranuloma venereum, a tropical sexually transmitted disease, is the result of infections with serotypes L₁-L₃. In industrial countries serotypes D-K are the cause of the most frequent sexually transmitted oculo-genital infections.

Most ***C. trachomatis*** infections in women are asymptomatic. This results in the many chronic diseases sustained by the ascended, persisting agents. Initial site of infections is usually the cervix with the clinical manifestation of cervicitis. Urethritis is less common in women. Between a primary infection and the appearance of chronic sequelae often months or even years will pass. Asymptomatic ***C. trachomatis*** infections in women may lead to endometritis, adnexitis,

periappendicitis, perihepatitis, peritonitis and reactive arthritis. This arthritis form is also known in men as post-urethritic reactive arthritis.

Repeated adnexitis may lead to secondary infertility. In men ascending infections following asymptomatic urethritis may result in epididymitis and prostatitis. Reduction of male fertility due to **C. trachomatis** infections is also discussed.

Peripheral **C. trachomatis** pathogens can be detected with high sensitivity and specificity by laboratory methods like ELISA, IFT, PCR, LCR. However, unrecognized or insufficiently treated **C. trachomatis** infections often lead to chronic manifestations. In these cases the antigen determination is limited and have to be completed by serology.

Chlamydiae contain common immunodominant antigens like the genus-specific lipopolysaccharide (LPS). To perform species specific serology **C. trachomatis**-specific immunogenic structures are needed. These structures are located in the variable domains of the major outer membrane protein (MOMP).

Until now the microimmunofluorescence test (MIF) is considered to be the so-called "gold standard" of species-specific serology. The MIF antigens consist of purified fixed total elementary bodies. Therefore a crossreactivity between the species cannot be excluded.

The **C. trachomatis-pELISA medac** uses a synthetic peptide from the immuno-dominant region of the MOMP. With this highly specific antigen a discrimination of **C. trachomatis**-specific antibodies from the whole anti-chlamydiae antibody response is possible with high specificity.

In infections with mainly chronic course serology is the method of choice.


Here, the detection of IgA antibodies in combination with IgG antibodies plays the decisive role. Serology also is indicated as additional method for the diagnosis of suspected peripheral **C. trachomatis** infections and for post-treatment follow up.

In addition to the **Chlamydia trachomatis-IgG-pELISA medac**
Cat.-no.: **497/TMB** for 96 determinations


we also distribute the following products:

Chlamydia trachomatis-IgA-pELISA medac
Cat.-no.: **498/TMB** for 96 determinations,
Chlamydia-IgG-rELISA medac
Cat.-no.: **480/TMB** for 96 determinations,
Chlamydia-IgA-rELISA medac
Cat.-no.: **490/TMB** for 96 determinations,
Chlamydia-IgM-rELISA medac
Cat.-no.: **485/TMB** for 96 determinations,
Chlamydia pneumoniae-IgG-rELISA medac
Cat.-no.: **430/TMB** for 96 determinations,
Chlamydia pneumoniae-IgA-rELISA medac
Cat.-no.: **431/TMB** for 96 determinations,
Chlamydia pneumoniae-IgM-rELISA medac
Cat.-no.: **482/TMB** for 96 determinations,
cHSP60-IgG-ELISA medac
Cat. no.: **435** for 96 determinations.


TEST PRINCIPLE


The plate is coated with synthetic *C. trachomatis*-specific peptid.

The *C. trachomatis*-specific antibodies from the specimen are bound to the antigen.

Peroxidase-conjugated anti-human IgG antibodies bind to the IgG antibodies (P = peroxidase).

Incubation with TMB-substrate (*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

Advantages of the test

- ☞ No false positive results caused by crossreactivities with other Chlamydiae.
- ☞ Non infectious antigen material.
- ☞ Chemically defined structure of the antigen.
- ☞ The breakable microwell strips permit efficient use of the test.

KIT CONTENTS

Cat. no.: 497/TMB

1.

MTP

Microplate: 12 x 8 wells, purple-coded (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with **C. trachomatis**-specific, synthetic peptid and NBCS, ready to use.
2.

CONTROL	-
----------------	----------

Negative control: 2 vials with 0.75 ml each, human serum, ready to use, stained blue, contains NBCS, phenol, ProClin™ 300 and gentamycine sulfate.
3.

CONTROL	+
----------------	----------

Positive control: 2 vials with 0.75 ml each, human serum, ready to use, stained blue, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
4.

WB

Wash buffer: 1 bottle with 100 ml, PBS/Tween (10 x), pH 7.2 - 7.4, contains ProClin™ 300.
5.

BAC-DIL

Sample diluent: 1 bottle with 110 ml, PBS/Tween/NBCS, pH 7.0 -7.2, ready to use, stained blue, contains ProClin™ 300.
6.

CON

Conjugate: 4 vials with 5.0 ml each, goat anti-human IgG, HRP-conjugated, ready to use, stained green, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
7.

TMB

TMB-substrate: 1 vial with 10 ml, ready to use.
8.

STOP

Stop solution: 2 vials with 14 ml each, 0.5 M sulfuric acid (H₂SO₄), ready to use.

1. STORAGE AND STABILITY

MATERIAL/REAGENT	STATE	STORAGE	STABILITY
Test kit	unopened	2 - 8 °C	until expiry date
Microplate	opened	2 - 8 °C in bag with desiccant	12 weeks
Controls	opened	2 - 8 °C	12 weeks
Wash buffer	diluted	2 - 8 °C	12 weeks
Sample diluent	opened	2 - 8 °C	12 weeks
Conjugate	opened	2 - 8 °C	12 weeks
TMB-substrate	opened	2 - 8 °C	12 weeks
Stop solution	opened	2 - 8 °C	until expiry date

Do not use the reagents after the expiry date.

2. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

- 2.1. Water for injection (H₂O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of Wash Buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37 °C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

3. PREPARATION OF THE REAGENTS

Before starting the test procedure all kit components must be equilibrated to room temperature.

Calculate the number of wells required.

3.1. Microplate

The aluminium bag has to be tightly resealed together with the desiccant after each removal of wells. Storage and stability of the wells are indicated under point 1.

3.2. Wash buffer

Mix one volume of wash buffer (10 x) with nine volumes of water for injection (e.g. 50 ml wash buffer (10 x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

Crystals in the wash buffer (10 x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.

Do not mix test specific reagents (microplate, controls, conjugate) from different kit lots. In contrast to that, sample diluent, wash buffer, TMB-substrate and stop solution are generally exchangeable in all Chlamydia-and Mycoplasma-ELISA medac.

Reagents from other manufacturers must not be used in general.

Valid and reproducible results are only obtained if the test procedure is precisely followed.

4. SPECIMEN

- 4.1. The test is suitable for serum samples.
- 4.2. Pretreatment of sera, e.g. inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain any red blood cells.
- 4.3. Sera have to be diluted 1:50 with sample diluent. They can be diluted further for titer determination.

5.A. TEST PROCEDURE

- 5.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

Microplate wells are ready to use and do not have to be prewashed.

- 5.2. Pipette 50 µl of the sample diluent into the well A1 as blank (see 6.A.), and 50 µl of the negative control (in duplicate), positive control, and the diluted patient samples.
- 5.3. Incubate the microplate wells for 60 min (± 5 min) at 37 °C (± 1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.4. After incubation wash the microplate wells three times with 200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

Do not allow the wells to dry out! Proceed immediately!

5.5. Add conjugate (coloured green) to each well.

50 µl of conjugate have to be pipetted into the wells if the test is done manually.

Please note:

When working with automated devices, 60 µl of conjugate have to be pipetted into each well due to a higher evaporation in the incubation chambers of the devices.

The suitability of the test for automated devices was confirmed during the evaluation of the test. Nevertheless we recommend to verify the compatibility of the test with the devices used in the lab.

5.6. Incubate again for 60 min (\pm 5 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil.

5.7. After incubation wash microplate wells again (see 5.4.).

5.8. Add 50 µl of TMB-substrate to each well and incubate for 30 min (\pm 2 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.

5.9. Stop the reaction by adding 100 µl of stop solution to each well. Positive samples turn yellow.

Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.

The reading should be done within 15 min after adding the stop solution!

5.B. TABLE FOR THE TEST PROCEDURE

	Blank (A1)	Negative control	Positive control	Sample
Sample diluent	50 µl	-	-	-
Negative control	-	50 µl	-	-
Positive control	-	-	50 µl	-
Sample	-	-	-	50 µl
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer.				
Conjugate	50/60 µl *)	50/60 µl *)	50/60 µl *)	50/60 µl *)
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer.				
TMB-substrate	50 µl	50 µl	50 µl	50 µl
Incubate for 30 min at 37 °C in the dark.				
Stop solution	100 µl	100 µl	100 µl	100 µl
Photometric reading at 450 nm (ref. 620 - 650 nm)				

*) manual/automatic procedure (see 5.5.)

6.A. CALCULATION OF RESULTS (VALIDITY)

- * Read OD values at 450 nm (reference wavelength 620 - 650 nm).
- * Subtract the OD value of the blank (well A1) from all other OD values.
- * The OD value of the blank has to be **< 0.100**.
- * The mean OD value of the **negative control** has to be **< 0.100**.
- * The OD value of the **positive control** has to be **> 0.800**.
- * **Cut-off = mean OD value of the negative control + 0.340**
- * **Grey zone = Cut-off ± 10%**

6.B. EVALUATION OF RESULTS:

6.B.1. QUALITATIVE

Result	Valuation
OD < Grey zone	negative
OD of Cut-off +/- 10%	equivocal
OD > Grey zone	positive

6.B.2. QUANTITATIVE

Cut-off-Index: $\frac{\text{OD Sample}}{\text{OD Cut-off}}$	Valuation
< 0.9	negative
0.9 - 1.1	equivocal
> 1.1	positive

- * Samples with OD values within the grey zone should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.
- * The results should always be interpreted in connection with IgA and with the clinical data and additional diagnostic parameters.
- * High concentrations of hemoglobin, of bilirubin and of lipid in serum do not have an influence on the results.

6.C. SPECIFIC IgG/IgA INTERPRETATION FOR SERUM

Possible Results		Interpretation
IgG	IgA	
+	+	1. IgG and IgA positive; indication of infection. Further judgement dependent on symptoms and titer course.
+	-	2. Only IgG positive; indication of a past infection. In case of clinical suspicion determine a fourfold IgG titer increase between two serum samples and retest for IgA.
+	+/-	3. IgG positive, IgA grey zone; possibility of a beginning or a subsiding infection; retest IgA after 10 - 14 days.
+/-	-	4. Only IgG in grey zone; past infection cannot be excluded. In case of clinical suspicion retest IgG and IgA after 10 - 14 days.
-	+	5. Only IgA positive; possibility of an early stage of infection or solitary, persisting IgA; retest IgA and IgG after 10 - 14 days.
-	+/-	6. Only IgA in the grey zone; possibility of a early stage of infection; retest IgA and IgG after 10 - 14 days.
+/-	+	7. IgG in the grey zone, IgA positive; possibility of an early stage of infection; retest IgA and IgG after 10 - 14 days.
-	-	8. IgG and IgA negative; no indication of infection. In case of justified clinical suspicion perform direct antigen detection and retest IgA and IgG after 10 - 14 days.

Comment:

In cases of fresh acute chlamydial infections the serological antibody results may be negative despite clinics and positive antigen detection. If a serological confirmation of a positive antigen result or if a follow-up is desired we recommend to test after 10 - 14 days for seroconversion.

7. PERFORMANCE CHARACTERISTICS

We determined the following performance characteristics during the diagnostic evaluation.

7.A. SPECIFICITY AND SENSITIVITY

Sera from children aged < 10 years and patients with serologically proven (MIF) *C. pneumoniae* infection were measured **to determine the specificity**.

Sera from three different patient groups were tested for the **determination of sensitivity/prevalence**.

Patient group	Specificity		
	IgA	IgG	
Children < 10 years	98% (n=100)	99% (n=100)	
Patients with <i>C. pneumoniae</i> infection (MIF positive)	100% (n=47)	100% (n=47)	
	Sensitivity/Prevalence		
STD patients culture negative for <i>C. trachomatis</i>	<i>C. trachomatis</i> - pELISA	6.2% (n=112)	25.0% (n=96)
	MIF	3.6% (n=112)	27.1% (n=96)
STD patients culture positive for <i>C. trachomatis</i>	<i>C. trachomatis</i> - pELISA	24.6% (n=114)	65.1% (n=106)
	MIF	20.2% (n=114)	67.0% (n=106)
Blood donors	<i>C. trachomatis</i> - pELISA	6.4% (n=299)	15.1% (n=299)
Arthritis patients ¹	<i>C. trachomatis</i> - pELISA	41.0% (n=83)	48.2% (n=83)

¹The sera were preselected according to a positive IgG result in the genus-specific Chlamydia rELISA.

7.B. PRECISION

Sample	Intra-assay variation				Sample	Inter-assay variation (n = 13)		
	mean OD	SD	CV (%)	n		mean OD	SD	CV (%)
NC	0.037	0.005	13.2	21	NC	0.038	0.011	30.4
BC	0.653	0.034	5.2	21	BC	0.371	0.045	12.0
PC	1.789	0.061	3.4	21	PC	1.502	0.094	6.2
N° 1	0.050	0.009	17.2	21	N° 3	1.724	0.102	5.9
N° 2	2.111	0.064	3.0	24	N° 4	0.525	0.070	13.4

NC = negative control; BC = weak positive control (not included in the kit);
PC = positive control

GENERAL HANDLING ADVICES

- * To avoid cross contamination do not exchange the vials and their screw caps.
- * The reagents have to be sealed immediately after use to avoid evaporation and microbial contamination.
- * After use, the reagents have to be stored as indicated to guarantee the shelf life.
- * After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

HEALTH AND SAFETY INFORMATION

- * The local occupational safety and health regulations have to be regarded.
- * Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents), should be handled as potentially infectious and used with all necessary precautions.

DISPOSAL CONSIDERATIONS

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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LITERATUR/REFERENCES/LITTÉRATURE

Christiansen, G., Pedersen, L., Clausen, J.D., Birkelund, S.: Cell and molecular biology of Chlamydia. In: Chlamydia Research. Angelika Stary (ed.). Proceedings of the Third Meeting of the European Society for Chlamydia Research, Vienna, Austria, 11.-14. September. pp. 3-6 (1996).

DeSchreyer, A., Meheus, A.: Epidemiology of sexually transmitted diseases: The global picture. Bull. World Health Organization 68, 639-654 (1990).

Hoyme, U.B., Spitzbart, H.: Past and current prevalence of **Chlamydia trachomatis** in women in Germany. In: Chlamydia Research. Angelika Stary (ed.). Proceedings of the Third Meeting of the European Society for Chlamydia Research, Vienna, Austria 11-14 September. p. 391 (1996).

Keck, C., Gerber-Schafer, C., Clad, A., Wilhelm, C., Breckwoldt, M.: Seminal tract infections: impact on male fertility and treatment options. Hum. Reprod. Update 4(6), 891-903 (1998).

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Nickel, J.C.: Chronic Prostatitis: An infectious disease? Infect. Urol. 13(2), 31-38 (2000).

Ochsendorf, F.R., Ozdemir, K., Rabenau, H., Fenner, T., Oremek, R., Milbradt, R., Doerr, H.W.: **Chlamydia trachomatis** and male infertility: chlamydia-IgA antibodies in seminal plasma are **C. trachomatis** specific and associated with an inflammatory response. Eur. Acad. Dermatol. Venereol. 12(2), 143-152 (1999).

Paavonen, J.: **Chlamydia trachomatis**: A major cause of mucopurulent cervicitis and pelvic inflammatory disease in women. In: Sexually Transmitted Diseases. Advances in Diagnosis and Treatment. Curr. Probl. Dermatol. Elsner. P., Eichmann, A. (eds.). Basel, Karger, 24, 110-122 (1996).

Patton, D., Askienazy-Elbhar, M., Henry-Suchet, J., Campbell, L.A., Cappucio, A., Tannous, W., Wang, S.-P., Kuo, S.-C.: Detection of **Chlamydia trachomatis** in fallopian tube tissue in women with post-infectious tubal infertility. Am. J. Obstet. Gynecol. 171, 95-101 (1994).

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Petersen, E.E., Clad, A., Mendel, R., Prillwitz, J., Hintz, K.: Prevalence of chlamydial infections in Germany: Screening of asymptomatic women and men by testing first void urine by ligase chain reaction (LCR). In: Chlamydia Research. Angelika Stary (ed.). Proceedings of the Third Meeting of the European Society for Chlamydia Research, Vienna, Austria, 11-14 September. p. 415 (1996).

Saikku, P.: Diagnosis of acute and chronic **Chlamydia pneumoniae** infections. In: Chlamydial Infections. Orfila, J., Byrne, G.I., Chernesky, M.A., Grayston, J.T., Jones, R.B., Ridgway, G.L., Saikku, P., Schachter, J., Stamm, W.E., Stephens, R.S. (eds.). Proceedings of the Eighth International Symposium on Human Chlamydial Infections, Chateau de Montvillargenne, 60270 Gouvieux-Chantilly, France, 19-24 June. p. 163-170 (1994).

Schachter, J.: The intracellular life of Chlamydia. Curr. Top. Microbiol. Immun. 138, 109-39 (1988).

Schachter, J.: Evolution of diagnostic tests for **Chlamydia trachomatis** infections. In: Chlamydia Research. Angelika Stary (ed.). Proceedings of the Third Meeting of the European Society for Chlamydia Research, Vienna, Austria, 11-14 September. pp. 243-146 (1996).

Schuppe, H.C., Bispink, G., Peet, D.J., Propping, D., Böttcher, M., Dettlaff, S.: The significance of antibodies against **C. trachomatis** in seminal plasma. IVth European Chlamydia Congress, 2000 Helsinki, im Druck

Ward, M.E.: The immunobiology and immunopathology of chlamydial infections. APMIS 103, 769-796 (1995).

Weström, L.V.: Chlamydia and its effect on reproduction. J. Brit. Fertil. Soc. 1, 23-30 (1996).

Wolff, H.: The biologic significance of white blood cells in semen. Fertil. Steril. 63(6), 1143-1157 (1995).

Wolff, H., Neubert, U., Zebhauser, M., Bezold, G., Korting, H.C., Meurer, M.: **Chlamydia trachomatis** induces an inflammatory response in the male genital tract and is associated with altered semen quality. Fertil. Steril. 55(5), 1017-1019 (1991).

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