**Mycoplasma pneumoniae-IgM-ELISA medac**

Deutsch/English/Français

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## **Mycoplasma pneumoniae-IgM-ELISA medac**

Enzyme immunoassay for qualitative determination of  
***Mycoplasma pneumoniae*** IgM antibodies in human serum

Cat.no.: 362

FOR IN VITRO DIAGNOSTIC USE ONLY

### **INTRODUCTION**

Mycoplasmas are bacteria without cell walls (Class Mollicutes = soft-skinned, order Mycoplasmatales) having a particularly small genome. They are characterised by their small size (300-800 nm) and their capacity to assume widely differing shapes (pleomorphism). At the present time over 100 species of the Order Mycoplasmatales are known. Most of them occur exclusively in animals. Fourteen species have so far been isolated from human beings. They colonise the mucous membranes of the respiratory and urogenital tracts. Only a few of them are pathogenic. ***Mycoplasma pneumoniae*** is among the species pathogenic to human beings. Two variants of ***M. pneumoniae*** are known.

***M. pneumoniae*** is an extracellular parasite of the mucosa of the respiratory tract. The pathogen is distinguished by its high level of host specificity. A specific adhesin is responsible for the adherence of the pathogen to the surface of the host cell. This adhesin also constitutes the virulence factor against which the crucial humoral response is directed. The transmission of *M. pneumoniae* takes place by droplet infection, but the incubation period of 10 to 20 days is relatively long.

***M. pneumoniae*** is endemic, occurring with seasonal peaks in the spring and autumn. Epidemics at intervals of 3-4 years are not unusual. Infections with ***M. pneumoniae*** are among the most frequent causes of community-acquired tracheobronchitis and atypical pneumonia in children and young people. The highest incidence is found between the ages of five and fifteen years. Of all cases of atypical pneumonia, 5-10% are attributable to ***M. pneumoniae***. It also causes conditions such as pharyngitis, laryngitis, otitis media and myringitis. Following a ***M. pneumoniae*** infection which begins with an infection of the respiratory tract, various other clinical manifestations may occur in other organs and systems: these include pericarditis and myocarditis, reactive arthritis, meningitis, meningoencephalitis, and polyneuritis together with a wide assortment of conditions affecting the skin. This at first seemingly harmless infection can therefore entail serious la-

te consequences. The disease may run a particularly severe course in immunodeficient and immunosuppressed patients.

A previous attack of *M. pneumoniae* infection does not leave behind it any immunity to this pathogen. Frequent reinfections are therefore observed.

In children under five years of age the infection is predominantly asymptomatic. In such cases the disease runs a mild course. In older children and in young people and adults, however, the infection is accompanied by prostration, headache, fever and a dry, unproductive cough. But these clinical features are totally nonspecific and do not provide any clue to the cause of the infection. An effective diagnostic test is therefore necessary.

In the acute phase of the disease demonstration of the causative pathogen is the method of choice. Tests for the pathogen can be performed on throat swabs, specimens of sputum and bronchoalveolar lavage fluid (BAL). The classical technique for demonstrating *M. pneumoniae* by means of cell culture takes a long time (10-14 days). Nevertheless the isolation of the microorganism is successful in only 40-60% of cases. The antigen-ELISA test requires material rich in cells, because the sensitivity of the test is relatively low. Effective qualitative DNA tests for this pathogen are not yet commercially available.

Serology constitutes the method of choice for patients with chronic infections, reinfections or extrapulmonary diseases due to *M. pneumoniae*. Commercial test systems include the complement fixation test (CFT), the haemagglutination test (HAT) and enzyme immunoassay (ELISA). Unlike the CFT and the HAT, the ELISA techniques make it possible to differentiate between immunoglobulin classes IgG, IgA and IgM, so that acute infections can be distinguished from chronic infections and past infections.


The **Mycoplasma pneumoniae-IgM-ELISA medac** employs a native, purified antigen for the serological detection.

In addition to the **Mycoplasma pneumoniae-IgM-ELISA medac**  
Cat.no.: **362** for 96 determinations


we also distribute the following products:

**Mycoplasma pneumoniae-IgG-ELISA medac**  
Cat.no.: **360** for 96 determinations,  
**Mycoplasma pneumoniae-IgA-ELISA medac**  
Cat.no.: **361** for 96 determinations.


## TEST PRINCIPLE


The plate is coated with a highly purified *M. pneumoniae*-specific antigen preparation.

The *M. pneumoniae*-specific antibodies from the specimen bind to the antigen.

Peroxidase-conjugated anti-human IgM antibodies bind to the IgM antibodies (P = peroxidase).

Incubation with TMB-substrate (\*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

### Advantages of the test

- ☞ High sensitivity and specificity.
- ☞ The breakable microwell strips permit efficient use of the test.

## **KIT CONTENTS**

**Cat. no.: 362**

1. 

<b>MTP</b>
------------

Microplate: 12 x 8 wells, (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with *M. pneumoniae*-specific antigen and FCS, ready to use.
2. 

<b>CONTROL</b>	<b>-</b>
----------------	----------

Negative control: 2 vials with 0.75 ml each, human serum, ready to use, stained blue, contains NBCS, phenol, ProClin™ 300 and gentamycin sulfate.
3. 

<b>CONTROL</b>	<b>+</b>
----------------	----------

Positive control: 2 vials with 0.75 ml each, human serum, ready to use, stained blue, contains BSA, phenol, ProClin™ 300 and gentamycin sulfate.
4. 

<b>WB</b>
-----------

Wash buffer: 1 bottle with 100 ml PBS/Tween (10x), pH 7.2 - 7.4, contains ProClin™ 300.
5. 

<b>BAC-DIL</b>
----------------

Sample diluent: 1 bottle with 110 ml PBS/Tween/NBCS, pH 7.0 - 7.2, ready to use, stained blue, contains ProClin™ 300.
6. 

<b>CON</b>
------------

Conjugate: 3 vials with 5.0 ml each, goat anti-human IgM, HRP-conjugated, ready to use, stained red, contains BSA, phenol, ProClin™ 300 and gentamycin sulfate.
7. 

<b>TMB</b>
------------

TMB-substrate: 1 vial with 10 ml, ready to use.
8. 

<b>STOP</b>
-------------

Stop solution: 2 vials with 14 ml each, 0.5 M sulfuric acid (H<sub>2</sub>SO<sub>4</sub>), ready to use.
9. 

<b>RF-ABS</b>
---------------

IgG/Rf-absorbent: 1 vial with 4 ml, goat anti-human IgG antibody, ready to use, contains < 0.1% sodium azide.

## **1. STORAGE AND STABILITY**

<b>MATERIAL/REAGENT</b>	<b>STATE</b>	<b>STORAGE</b>	<b>STABILITY</b>
Test kit	unopened	2 - 8 °C	until expiry date
Microplate	opened	2 - 8 °C in bag with desiccant	12 weeks
Controls	opened	2 - 8 °C	12 weeks
Wash buffer	diluted	2 - 8 °C	12 weeks
Sample diluent	opened	2 - 8 °C	12 weeks
Conjugate	opened	2 - 8 °C	12 weeks
TMB-substrate	opened	2 - 8 °C	12 weeks
Stop solution	opened	2 - 8 °C	until expiry date
IgG/Rf-absorbent	opened	2 - 8 °C	12 weeks

Do not use the reagents after the expiry date.

## **2. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED**

- 2.1. Water for injection (H<sub>2</sub>O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of wash buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37 °C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

## **3. PREPARATION OF THE REAGENTS**

Before starting the test procedure all kit components must be equilibrated to room temperature (RT).

Calculate the number of wells required.

### 3.1. Microplate

After each removal of wells the aluminium bag has to be tightly resealed together with the desiccant. Storage and stability of the wells are indicated under point 1.

### 3.2. Wash buffer

Mix one volume of wash buffer (10x) with nine volumes of water for injection (e.g. 50 ml wash buffer (10x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

**Crystals in the wash buffer (10x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.**

**Do not mix test specific reagents (microplate, controls, conjugate) from different kit lots. In contrast to that, sample diluent, wash buffer, TMB-substrate, IgG/Rf-absorbent and stop solution are generally exchangeable in all Chlamydia-and Mycoplasma-ELISA medac.**

**Reagents from other manufacturers must not be used in general.**

**Valid and reproducible results are only obtained if the test procedure is precisely followed.**

## **4. SPECIMEN**

- 4.1. The test is suitable for serum samples.
- 4.2. In order to avoid interferences with high IgG titers and rheumatoid factors (RF) an absorption of IgG/RF shall be performed.
- 4.3. An additional pretreatment of sera, e.g. inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain red blood cells.

## **5. TEST PROCEDURE**

### **5.A. IgG/RF-ABSORPTION**

**Attention:**

- \* **The controls are ready to use (no absorption necessary).**
  - \* **The following volumes indicated apply for single determinations only.**
- 5.A.1. Serum: 10 µl serum are diluted with 490 µl sample diluent (dilution 1:50).
  - 5.A.2. Absorption: 30 µl IgG/Rf-Absorbent and 30 µl diluted serum are mixed (dilution 1:100) and incubated for 15 min at room temperature.

Alternative: The absorption can be performed overnight at 2 - 8 °C.

5.A.3. Now the test dilution is 1:100.

## **5.B. TEST RUN**

5.B.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

**Microplate wells are ready to use and do not have to be prewashed.**

5.B.2. Pipette 50 µl of the sample diluent into the well A1 as blank (see 6.A.), and 50 µl of the negative control (in duplicate), positive control, and the diluted patient samples.

**If necessary the microplate wells can be kept for max. 30 min in a humid chamber at RT before proceeding.**

5.B.3. Incubate the microplate wells for 60 min (± 5 min) at 37 °C (± 1 °C) in a humid chamber or sealed with incubation cover foil.

5.B.4. After incubation wash the microplate wells three times with 200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

**Do not allow the wells to dry out! Proceed immediately!**

5.B.5. Add conjugate (coloured red) to each well.

**50 µl of conjugate have to be pipetted into the wells if the test is done manually.**

### **Please note:**

**When working with automated devices, 60 µl of conjugate have to be pipetted into each well due to a higher evaporation in the incubation chambers of the devices.**

**The suitability of the test for automated devices was confirmed during the evaluation of the test. Nevertheless we recommend to verify the compatibility of the test with the devices used in the lab.**

5.B.6. Incubate again for 60 min (± 5 min) at 37 °C (± 1 °C) in a humid chamber or sealed with incubation cover foil.


- 5.B.7. After incubation wash microplate wells again (see 5.B.4.).
- 5.B.8. Add 50  $\mu$ l of TMB-substrate to each well and incubate for 30 min ( $\pm$  2 min) at 37  $^{\circ}$ C ( $\pm$  1  $^{\circ}$ C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.
- 5.B.9. Stop the reaction by adding 100  $\mu$ l of stop solution to each well. Positive samples turn yellow.

**Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.**

**The reading must be done within 15 min after adding the stop solution!**

**5.C. TABLE FOR THE IgG/RF-ABSORPTION**

Indication per determination

#### **5.D. TABLE FOR THE TEST PROCEDURE**

	Blank (A1)	Negative control	Positive control	Sample
Sample diluent	50 µl	-	-	-
Negative control	-	50 µl	-	-
Positive control	-	-	50 µl	-
Absorbed sample	-	-	-	50 µl
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer.				
Conjugate	50/60 µl *)	50/60 µl *)	50/60 µl *)	50/60 µl *)
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer.				
TMB-substrate	50 µl	50 µl	50 µl	50 µl
Incubate for 30 min at 37 °C in the dark.				
Stop solution	100 µl	100 µl	100 µl	100 µl
Photometric reading at 450 nm (ref. 620 - 650 nm)				

\*) manual/automatic procedure (see 5.B.5.)

#### **6.A. CALCULATION OF RESULTS (VALIDITY)**

- \* Read OD values at 450 nm (reference wavelength 620 - 650 nm).
- \* Subtract the OD value of the blank (well A1) from all other OD values.
- \* The OD value of the blank has to be **< 0.100**.
- \* The mean OD value of the **negative control** has to be **< 0.100**.
- \* The OD value of the **positive control** has to be **> 0.800**.
- \* **Cut-off = mean OD value of the negative control + 0.380**
- \* **Grey zone = Cut-off ± 10%**.

**6.B. EVALUATION OF RESULTS:**

**6.B.1. QUALITATIVE**

Result	Valuation
OD < Grey zone	negative
OD of Cut-off +/- 10%	equivocal
OD > Grey zone	positive

**6.B.2. SEMIQUANTITATIVE**

Cut-off-Index: $\frac{\text{OD Sample}}{\text{OD Cut-off}}$	Valuation
< 0.9	negative
0.9 - 1.1	equivocal
> 1.1	positive

- \* Samples with OD values within the grey zone should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.
- \* The results should always be interpreted in connection with IgG and IgA, and with the clinical data and additional diagnostic parameters.
- \* High concentrations of hemoglobin in serum do not have an influence on the results. However, high concentrations of lipids might increase the OD values of negative samples.
- \* Cross-reactivities with heterophilic antibodies cannot be excluded in individual cases.
- \* In recent EBV-infections false positive IgM results may occur in certain circumstances.

**6.C. SPECIFIC IgM-/IgA-/IgG-INTERPRETATION**

Possible results			Interpretation
IgM	IgA	IgG	
+	-	-	1. Indication of early stage of infection or solitary, persisting IgM. Retest IgM, IgA and IgG after 14 days. <sup>1,2</sup>
-	+	-	2. Indication of early stage of infection or solitary, persisting IgA. Retest IgA and IgG after 14 days.
+	+	-	3. Indication of acute infection. <sup>1,2</sup> Retest IgG after 14 days.
+	-	+	4. Indication of acute infection.
+	+	+	5. Indication of acute infection.
-	+	+	6. Indication of primary current or reinfection. <sup>3</sup>
-	-	+	7. Indication of past infection. In case of clinical suspicion retest after 14 days for IgA and IgG antibodies.
-	-	-	8. No serological indication of current or past infection. In case of clinical suspicion retest after 14 days for IgM, IgA and IgG antibodies.

**Comment:**

- Borderline results may indicate commencing or abating infections. Retesting after 14 days is recommended.
- <sup>1</sup> Current, acute infections are best detectable by parallel determination of IgM and IgA antibodies.
- <sup>2</sup> Simultaneous detection of IgM and IgA antibodies is particularly frequent in children.
- <sup>3</sup> In adults, IgA antibodies are more reliable markers for current infections than IgM antibodies.

## 7. PERFORMANCE CHARACTERISTICS

We determined the following performance characteristics during the diagnostic evaluation.

### 7.A. SPECIFICITY AND SENSITIVITY

<b>Patient Group</b>	<b>Specificity IgM</b>
<i>M. pneumoniae</i> IgM negative sera (reference ELISA/ Aggl. Assay)	<b>100% (n=40)</b>
<b>Patient Group</b>	<b>Sensitivity IgM</b>
Sera <b>with</b> IgM antibodies to <i>M. pneumoniae</i> (reference ELISA/Aggl. Assay) Patients with respiratory infections	<b>71% (n=44)</b>

### 7.B. PRECISION

Sample	Intra-assay variation				Sample	Inter-assay variation (n = 13)		
	mean OD	SD	CV (%)	n		mean OD	SD	CV (%)
<b>NC</b>	0.028	0.008	29	22	<b>NC</b>	0.017	0.004	24
<b>BC</b>	0.838	0.019	2	22	<b>BC</b>	0.792	0.050	6
<b>PC</b>	1.141	0.020	2	22	<b>PC</b>	0.990	0.035	4
<b>N° 1</b>	0.033	0.004	12	22	<b>N° 1</b>	0.021	0.004	19
<b>N° 2</b>	0.651	0.015	2	22	<b>N° 2</b>	0.542	0.017	3
<b>N° 3</b>	1.611	0.033	2	22	<b>N° 3</b>	1.480	0.062	4
					<b>N° 4</b>	1.409	0.093	7

NC = negative control; BC = weak positive control (not included in the kit);  
PC = positive control

### **GENERAL HANDLING ADVICES**

- \* To avoid cross contamination do not exchange the vials and their screw caps.
- \* The reagents have to be sealed immediately after use to avoid evaporation and microbial contamination.
- \* After use, the reagents have to be stored as indicated to guarantee the shelf life.
- \* After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

### **HEALTH AND SAFETY INFORMATION**

- \* The local occupational safety and health regulations have to be regarded.
- \* Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents), should be handled as potentially infectious and used with all necessary precautions.

### **DISPOSAL CONSIDERATIONS**

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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## LITERATUR/REFERENCES

Clyde, WA: Clinical overview of typical *Mycoplasma pneumoniae* infections. Clin Infect Dis 17 (Suppl 1), S32-6 (1993).

Dionisio D, Valassina M, Uberti M, Fabbri C, Parri F, Saffi EG: *Mycoplasma pneumoniae* non-pulmonary infection presenting with pharyngitis, polyarthritits and localized exanthem. Scand J Infect Dis 33, 782-783 (2001).

Döller G, Döller PC, Jacobs E, Schuy W: Zur Differentialdiagnostik von Infektionen des Respirationstrakts. Diagnose & Labor 43, 21-33 (1993).

Drasbek M, Nielsen PK, Persson K, Birkelund S, Christiansen G: Immune response to *Mycoplasma pneumoniae* P1 and P116 in patients with atypical pneumonia analyzed by ELISA. BMC Microbiol. 4, 7ff. (2004).

Ferwerda A, Moll HA, de Groot R: Respiratory tract infections by *Mycoplasma pneumoniae* in children: a review of diagnostic and therapeutic measures. Eur J Pediatr 160, 483-491 (2001).

Granstrom M, Holme T, Sjogren AM, Ortqvist A, Kalin M: The role of IgA determination by ELISA in the early serodiagnosis of *Mycoplasma pneumoniae* infection, in relation to IgG and mu-capture IgM methods. Med Microbiol. 40, 288-92 (1994).

Hammerschlag MR: *Mycoplasma pneumoniae* infections. Curr Opin Infect Dis 14, 181-186 (2001).

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Krause D: *Mycoplasma pneumoniae* cytoadherence: organization and assembly of the attachment organelle. Trends Microbiol 6, 15-18 (1998).

Seggev JS, Semak GV, Kurup VP: Isotype-specific antibody responses to acute *Mycoplasma pneumoniae* infection. Ann Allergy Asthma Immunol. 77, 66-73 (1996).

Sillis M: The limitations of IgM assays in the serological diagnosis of *Mycoplasma pneumoniae* infections. Med Microbiol. 33, 253-8 (1990).

Sotgiu S, Pugliatti M, Rosati G, Deiana GA, Sechi GP: Neurological disorders associated with *Mycoplasma pneumoniae* infection. Eur J Neurol 10, 165-168 (2003).