recomBlot HCV IgG 2.0

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recomBlot HCV IgG 2.0

Immunoblot test with recombinant antigens separated by gel electrophoresis, for the determination of IgG antibodies against the **Hepatitis C virus (HCV)** in human serum or plasma.

1. General aspects

recomBlot HCV IgG 2.0 is a qualitative *in vitro* test for the determination and identification of IgG antibodies against Hepatitis C virus in human serum or plasma. In contrast to ELISAs or Spottests the test permits a statement about the specific antibodies assigned to different stages of infection by the electrophoretical separation of the antigens. The **recomBlot HCV IgG 2.0** is used to confirm results for samples that have tested positive in screening assays.

2. Hepatitis C virus

Before its identification, the hepatitis C virus was classified as the so-called nonA/nonB hepatitis viruses that cause most of the cases of hepatitis transmitted by blood transfusions and by other blood products. Cases of nonA/nonB hepatitis are characterised by an incubation period of 2 to 26 weeks (average: 8 weeks) and, in most cases, a relatively mild course in the acute phase. As many as 80% of patients develop a chronic hepatitis, which in some cases results in a liver cirrhosis or primary hepatic cell carcinoma. The number of chronically infected persons world-wide is estimated at 170 million. The prevalence level in Germany is 0.6%.

The causative infective agent of this hepatitis type was identified in 1989 as a single-strand RNA virus, which was then designated as the hepatitis C virus (HCV). This pathogen belongs to the family of flaviviruses, which also includes the yellow fever and TBE (tick-borne encephalitis) viruses. Its genome has a length of 9.5 Kb and codes for one capsid protein (core), two membrane proteins (E1 and E2) and four non-structural proteins (NS2 to NS5).

2.1 Epidemiology and transmission

Until suitable testing methods were available, the hepatitis C virus caused 80 to 90% of nonA/nonB hepatitis infections. This virus infects only humans and is usually transmitted via blood transmissions or other blood products. One in every three to five thousand blood donations is positive. The risk of developing a clinical infection after receiving a dose of positive blood is about 75%.

80% of all new transmissions are now observed in drug addicts due to syringe and needle sharing. Other potential transmission routes include sexual intercourse and household contact with infected persons under poor hygienic conditions. Hospital staff are at risk from needle injuries. The virus can be transmitted from mother to child during pregnancy. The transmission route is unknown in about 30% of cases.

2.2 Pathogenesis and clinical aspects

The virus enters the bloodstream directly, usually in contaminated blood or blood products. It is transported to the liver by infected macrophages, where it infects and destroys hepatocytes. The result is a liver inflammation (hepatitis) with cellular necrosis.

About 75% of infections take an inapparent course, severe courses being rare exceptions. Chronic persistent or chronic reactive hepatitis develops in as many as 80% of all infected persons. In 10 to 20% of the chronic cases, cirrhosis of the liver eventually develops, and about 4% of these cases develop a primary hepatocellular carcinoma.

3. Diagnosis

The first step in laboratory diagnostics of a hepatitis C infection is the detection of specific antibodies against HCV proteins using ELISA (enzyme-linked immunosorbent assay) screening tests. Then ELISA reactivity is confirmed by the more highly specific Western Blot test. The higher level of specificity of the latter test is achieved by separating the HCV antigens. Thereby their respective reactivity levels can be evaluated separately.

A positive result for HCV antibodies recommends the viral RNA detection using the NAT (nucleic acid amplification technique), for example RT/PCR. NAT should confirm or exclude virus persistence, which may result in therapeutic consequences. In the case of negative PCR coupled with confirmed antibody

detection, the RT/PCR should be repeated at intervals of a few months. Only then a relatively reliable conclusion can be drawn regarding the chronicity of a HCV infection.

Antibodies to HCV-specific antigens are usually detectable three to four weeks after a HCV infection. An exception is the "diagnostic window" during the incubation and very early acute phase of an infection. In that time the very low viral load delays the antibody response. Immunosuppressed patients and neonates with perinatal infections also show a "distorted" antibody picture (neonates due to antibodies from the mother). A confirmed diagnosis of these HCV infections is only possible by repeated positive detection of the HCV genome using the NAT.

4. Test principle

The purified recombinant antigens are separated by molecular weight via SDS polyacrylamide gel electrophoresis (SDS-PAGE) and subsequently transferred to nitrocellulose membrane electrophoretically (western blotting). Free binding sites on the membrane are saturated with a solution of proteins. Following this matrix is washed and cut into strips.

To facilitate detection of HCV-specific antibodies, the strips are incubated with the diluted serum sample, whereby the antibodies bind to the antigens on the strips. Unbound antibodies are then flushed away and the strips are incubated in a second step with anti-human IgG coupled with horse radish peroxidase. Specifically bound antibodies are detected by means of a colour reaction catalysed by the peroxidase. If reaction against one of the HCV-specific proteins has taken place, a dark band appears at the corresponding locus on the strip.

As reaction control, a band with anti-human immunoglobulin is also applied at the upper end of the strips (below the number) that must show a reaction to every serum.

The present immunoblot uses five recombinant produced HCV antigens:

NS-3:	Non-structural protein with activity of viral helicase and protease
Helicase:	Part of NS-3 with activity of viral helicase
NS-5-12:	Part of NS5 (RNA polymerase) of HCV
NS-4:	Unknown function, N-terminal part of NS-4
Core:	Capsid protein, forming the major part of the nucleocapsid

In order to avoid cross reactions and false positive results, the recombinant antigens are produced without any foreign protein.

5. Package contents

recomBlot HCV IgG 2.0

The reagents in a pack are sufficient for 20 determinations.

Each reagent set contains:

WASHBUF 10 X	100 ml	Wash and dilution buffer (ten times the concentration) contains phosphate buffer, NaCl, KCl, detergent and preservative
SUBS TMB	45 ml	Substrate solution Tetramethylbenzidin (TMB, ready-to-use)
MILKPOW	5 g	Skim milk powder
INCUTRAY	2 pieces	Incubation trays with 10 wells each
CONTRSTR	1	Predeveloped control strip (kit specific)
INSTRU	1	Instructions for use
EVALFORM	1	Evaluation form
TESTSTR	2 pieces	Test tubes with 10 consecutively numbered western blot test strips coated with HCV antigens
CONTROL ± IgG	80 µl	Weak positive serum control (white screw cap) Human origin, anti HIV-1/2 and HBs-Ag negative, contains 0,1% MIT and 0,1% Oxypyrion
CONTROL - IgG	80 µl	Negative serum control (blue screw cap) Human origin, anti HIV-1/2, anti-HCV and HBs-Ag negative, contains 0,1% MIT and 0,1% Oxypyrion
CONJ IgG	60 µl	Anti-human IgG conjugate (green screw cap, thousand times the concentration) from rabbit, contains NaN ₃ (<0,1%)

6. Additional reagents and accessory equipment required

Deionised water, vacuum extraction system with disinfection trap, micro pipettes, plastic forceps, shaker, metric cylinder with graduations, scales.

7. Information on test and reagents

7.1 Precautions

- ☞ Control sera are from blood donors verified for the absence of antibodies to HIV 1/2 and Hepatitis Bs-antigen. Since an infection cannot be excluded with absolute certainty despite this precaution, the product must be treated with the same care as the patient samples.
- ☞ Suitable single-use gloves must be worn during the entire test procedure.
- ☞ The conjugates contain sodium azide. Avoid contact with skin or mucosa.
- ☞ All reagents and materials contaminated with potentially infectious samples must be treated with suitable disinfectants or autoclaved at 121°C for at least 1 hour.
- ☞ Use incubation trays once only.
- ☞ Handle strips carefully with a plastic forceps.

7.2 Handling information

Store reagents before and after use at 2°C - 8°C, **do not freeze**. Temper all components before starting the test for at least 30 minutes to 18°C - 25°C (room temperature). Both the test and incubation procedures are carried out at room temperature.

Mix the control sera and concentrated conjugates well before use. Mix the patient sera well before use as well.

Do not open the tube containing the test strips until just before use to avoid condensation of water. The strips that are not used must be left in the tube and stored further at 2°C - 8°C (reclose tube tightly, test strips must not become moist before testing!).

Weak positive and negative controls must always be run parallel to the sera being tested, regardless of how many are being tested. This is the only way to ensure correct test interpretation.

A quality guarantee can only be given up to the expiry date on the packages.

Only reagents with lot numbers corresponding to the respective lot number on the label of the kit package may be used.

The test must be performed by well-trained and authorised qualified personnel.

In case of substantial modifications of the product or the instructions for use, the application of the test might differ from the purpose intended by MIKROGEN.

Automation is possible, please refer to MIKROGEN for details.

7.3 Preparation of solutions

7.3.1 Preparation of ready-to-use wash buffer

This buffer is used for both serum and conjugate dilution as well as for the washing steps.

Prior to dilution, the volume of wash buffer required for the corresponding number of tests must be determined.

The skim milk powder is first dissolved in wash buffer concentrate. This mixture is then filled up with deionised water to the final volume (dilution: 1 + 9). See Table 1 for the amounts required. Volumes for numbers of test strips not listed in the table have to be determined by calculation.

Ready-to-use wash buffer can be stored at 2°C - 8°C for four weeks.

Table 1: Wash buffer per test strip used

Test strips used	Skim milk powder	Wash buffer concentrate	Deionised water	Ready-to-use wash buffer
1	0,1 g	2 ml	18 ml	20 ml
2	0,2 g	4 ml	36 ml	40 ml
3	0,3 g	6 ml	54 ml	60 ml
5	0,5 g	10 ml	90 ml	100 ml
10	1 g	20 ml	180 ml	200 ml
20	2 g	40 ml	360 ml	400 ml

7.3.2 Preparation of conjugate solution

The conjugate solution for the IgG test is to be prepared immediately before use. It is not possible to store the ready-to-use conjugate solution.

One part of IgG conjugate concentrate is diluted with 1000 parts of ready-to-use wash buffer (1 + 1000).

The amounts used are listed in Table 2. Volumes for numbers of test strips not listed in the table have to be determined by calculation.

Table 2: Volumes for anti-human IgG conjugate dilution

Test strips used	IgG conjugate concentrate	Ready-to-use wash buffer
1	2 µl	2 ml
2	4 µl	4 ml
3	6 µl	6 ml
5	10 µl	10 ml
10	20 µl	20 ml
15	30 µl	30 ml
20	40 µl	40 ml

* The volumes have been calculated without dead volume. Depending on handling (manually or by machine processing) please prepare conjugate solution for additional 1 to 3 strips.

7.3.3 Substrate solution

The substrate is ready for use! Bring to room temperature (18°C - 25°C) before starting the colour reaction.

Avoid contamination of the unused substrate solution by nonsterile pipette tips etc. at all costs since this may affect test sensitivity.

7.4 Storage and stability

Store reagents at 2°C - 8°C before and after use.

Ready-to-use wash buffer can be stored at 2°C - 8°C for four weeks.

The conjugate solution is to be prepared always freshly.

8. Sample material

The sample material can be serum or plasma that is separated from the coagulum as soon as possible after sampling. Heat-inactivated samples may result in raised background reaction levels. A microbial contamination of the sample is to be avoided at all costs. Insoluble substances are to be removed from the sample prior to incubation by means of centrifugation.

Use of lipaemic, haemolytic or icteric samples did not result in incorrectly evaluated data (see 12). Therefore these sera can be analysed by the *recomBlot* HCV IgG 2.0. Please consider the possibility of dark background which may be generated by these sera on the whole nitrocellulose strip.

Important!

If the tests are not carried out immediately, the samples can be stored for up to 2 weeks at 2°C-8°C. Longer storage of the samples is possible at -20°C or colder. Repeated freezing and thawing of the samples is not recommended because this may affect the quality of the results.

9. Test procedure

9.1 General

The reproducibility of the results depends to a great extent on consistent washing of the strips. The washing frequencies described under 9.3. and 9.5 should therefore always be maintained.

9.2 Incubation of samples

1. One well in the incubation tray (see 6) is required per sample to be tested. **2 ml** of ready-to-use wash buffer is pipetted into each well. Then one test strip is carefully dipped into each of the wells filled with ready-to-use wash buffer using a plastic forceps. The strip number must face upwards.

Important!

The strip must be completely wet and immersed in the wash buffer.

Record the tube and strip numbers used on the evaluation sheet.

2. Adding the samples

IgG test procedure: 20 µl of an undiluted sample (human serum or plasma) or of the corresponding weak positive control resp. negative control, are pipetted into the proper wells (**dilution 1 + 100**).

Please be sure to add the sample at one end of the immersed strips into the wash buffer and mix as soon as possible by shaking the tray carefully.

Record the sample numbers on the evaluation sheet.

Cover the incubation tray with the plastic lid and incubate for **1 hour** at room temperature while shaking gently.

Important!

Make sure the incubation solutions are not carried over into other wells; be particularly careful to avoid splashing when opening and closing the lid (risk of cross-contamination).

9.3 Washing

1. Following incubation, the plastic lids are carefully removed from the incubation trays.
2. The serum dilution is carefully aspirated from the individual wells.

Important!

When the solutions have been aspirated from a well, the pipette tip has to be changed or rinsed well with deionised water after each aspiration procedure to prevent cross-contamination.

In the case of machine processing the references of the equipment manufacturer have to be considered.

3. Then place **2 ml** of the ready-to-use wash buffer in each well and wash on the shaker while shaking gently for **5 minutes**. The wash buffer is aspirated after the washing procedure.
4. Carry out the washing step under 3. a total of **three times**.

9.4 Incubation with peroxidase conjugate

After washing the strips, place **2 ml** of the appropriately prepared **conjugate solution** (see Table 2). into each well and incubate for **45 minutes** while shaking gently at room temperature, whereby the incubation tray is covered with the plastic lid.

9.5 Washing

The conjugate solutions are aspirated from the wells and the strips are washed again (see 9.3).

9.6 Substrate reaction

Add **1.5 ml substrate solution** into each well and incubate for **5 - 15 minutes**, shaking gently and under observation, at room temperature.

Important!

Observe the colour reaction process and leave all strips in the substrate solution until the strips incubated with the weak positive serum control shows the core band at least. Further band can also appear on the wet strip (s. 11.2)!

The colour reaction may be excessive with highly positive sera. Recommendation: Stop the colour reaction earlier in such strips.

9.7 Stopping the reaction

1. After the solution is aspirated, wash the strips briefly **three times with deionised water**.
2. Use a plastic forceps to remove the strips carefully from the water and place them between 2 layers of absorbent paper to dry for about 2 hours. The strips may be adhesively attached to the enclosed evaluation sheet and the results may be recorded.
3. The strips should be stored protected from exposure to light.

10. Summary of the test procedure

1.	bring all reagents to room temperature
2.	deposit the strips in 2 ml ready-to-use wash buffer and take care that they are completely immersed in the ready-to-use wash buffer
3.	add 20 µl of patient sample resp. control serum
4.	incubate for 1 hour at room temperature while shaking gently
5.	wash on the shaker three times with 2 ml of ready-to-use wash buffer for 5 minutes each time
6.	add 2 ml appropriately prepared conjugate solution
7.	incubate at room temperature for 45 minutes while shaking gently
8.	wash on the shaker three times with 2 ml of ready-to-use wash buffer for 5 minutes each time
9.	add 1.5 ml of the substrate solution, incubate while shaking at room temperature for 5 – 15 minutes. Bands see 9.6
10.	wash the strips at least three times with deionised water
11.	dry the strips between 2 layers of absorbent paper for 2 hours and read off the result

11. Evaluation

11.1 Evaluation of band intensity

1. A weak positive control and negative control must be included in each test run, irrespective of the number of samples. The reactivity of the weak positive control serves as cutoff and as a discrimination from the background reactivity.
2. On the attached evaluation sheet, record the date, batch and tube number, and the detected antibody class.
3. Enter the sample identification number on the protocol sheet.
4. Now attach the corresponding test strips with a glue stick into the corresponding fields of the evaluation sheet. To do this, place the test strips with the reaction control band on the marking lines. Then use clear adhesive tape to attach the test strips left from the marking line. Complete sticking of the test strip with glue or adhesive tape may lead to aberrations of the colouring.
5. Place the predeveloped control test strip with the labelled antigen bands on the fixed test strips in such a way that the reaction control bands are aligned. The predeveloped test strip is kit specific and comes from the same nitrocellulose sheet as the test strips in the test kit. The molecular weights and names of the antigen bands are marked.
6. Band intensities are assessed by comparison with the weak positive control (refer to 11.2 resp. Table 4). Please record the values in the evaluation sheet.


11.2 Control results

Weak positive controls must always be run parallel to the sera being tested. The test conditions (sensitivity) are to be examined and reactivity with the individual antigen bands are to be confirmed. The localisation of the antigen bands must agree with the attached control strip (kit specific).

The control band (reaction control) below the strip number is used to check the correct performance of the test. This band must show a distinct colouring.

The test strips incubated with the weak positive serum control show the following band patterns.

- Weak positive control**

The reactivity of the weak positive control serves as sensitivity check and as cutoff. The dry strip must show the reaction control band, Helicase and core. Further bands may also appear. In the *wet state* and during the *staining process*, at least the core band should be weakly visible. All of the other bands may appear (see 9.6).

- Negative control:**

The negative control does not show any band patterns. The reaction control must only be seen (circa 2 - 6 mm below the identification number of each strip). The reaction control serves as a control for the proper incubation of serum, conjugate and substrate.

Table 3: Band intensity

Bands	Intensity
Negative	-
Weak intensity (visibly, but weaker than the core antigen of the weak positive control)	+/-
Intensity well visible (same intensity as with the core-antigen of the weak positive control)	+
Strong intensity (stronger intensity as with the weak positive control)	2+

11.3 Test results and test interpretation

As a basis for reliable and simplified test evaluation, band specific point values for the determination of the test results were established. The point values of the HCV antigens in the *recomBlot HCV IgG 2.0* based on the clinical evaluations and mathematical analysis.

The test result is obtained by the addition of the point values of the separate bands, that are determined as ±, + or ++ (Table 4 and Table 5). The resulting total of points is entered into the column marked with the sign for sums. The positive, questionable or negative evaluation of the sample can then be directly determined and entered into the evaluation column.

Table 4: Point values of the HCV antigens of the *recomBlot HCV IgG 2.0*

	+/-	+ and ++
NS3	0	3
Helicase	0	3
NS5-12	0	2
NS4	0	4
Core	5	8

Table 5: Interpretation of the test results

Sum of points	Evaluation
0 - 5	negative
6 - 9	borderline
> 9	positive

11.4 Directions for Interpretation

It should be noted, that deviations from the test procedure may lead to misinterpretation. The given incubation times must be strictly observed.

Strong positive samples may show an overreaction with the antigens on the strips, resulting in a strong colour reaction. Additional bands may occur as well. It is therefore recommended to repeat the test in cases of strong positive samples with a dilution of 1:500. Dilution in ready-to-use wash buffer.

Samples with borderline results have to be retested with a freshly obtained serum sample after three to four weeks. If samples with isolated reactions against NS-antigens do not show any change, it might be caused by cross reactions. Furthermore it may indicate a residual titer of a past infection or an unspecific reaction with the recombinant antigens.

For confirmation of a positive test result, NAT (e.g. PCR) is recommended for the detection of the HCV genome.

There is no correlation between positive detection of HCV antibodies and infectiousness.

Negative test results do not exclude a virus infection. It is possible that HCV antibodies may not be detected in haemodialysed patients despite positive test results for HCV RNA (Hinrichsen et al., 2002; Bukh et al., 1993).

Dark test strips: A few sera of patients may generate a dark general or patterned colouring on the whole nitrocellulose strip (e.g. sera of patients with milk protein allergies). Different factors according to the patient serum are responsible for this phenomenon. The evaluation of these strips is usually possible only with some restrictions. E.g. "inverse" bands (white bands on dark background) are assessed as negative. The corresponding serum should strictly be retested by use of other serological methods.

12. Clinical results

A total of 1,137 sera were tested with the *recomBlot* HCV IgG 2.0 for clinical evaluation. The collective comprised different panels of sera: 399 blood donor sera, 299 HCV-positive sera, pre-tested samples from routine diagnostic procedures, 18 seroconversions (with a total of 156 sera), 50 potentially interfering samples (icteric, haemolytic, lipaemic and RF-positive samples), 33 sera with antibodies against the recombinant antigen expression system and 200 clinical sera from patients with EBV infections (n = 20), HEV infections (n = 22), pregnancy (n = 50), clinically defined sera (n = 30) and routine diagnostics sera (n = 78).

	recomBlot HCV IgG 2.0			
	total	negative	indeterminate	positive
Blood donors	399	380	19	0
Clinical sera EBV infections, pregnancy sera, HEV patients, clinically determined sera, sera of routine diagnostics	200	191	6	3 ¹⁾
Potentially interfering sera (icteric, haemolytic, lipaemic and RF-positive samples)	50	47	3	0
Sera with antibodies against the recombinant antigen expression system	33	33	0	0
HCV-positive, pretested samples	299	1 ²⁾	1	297
Seroconversions (18 courses)	156	39	20	97
Total	1137	691	49	397

1) Positive results were confirmed by comparative assays.

2) The serum is weak positive in the ELISA, positive in PCR and negative in several confirmatory tests. According to the external referee, the sera from routine diagnostics also include haemodialysis patients with positive PCR results and negative HCV antibody titer (11.4).

Testing of the 18 seroconversions, with an average of 7-8 samplings, showed that *recomBlot* HCV IgG 2.0 results in indeterminate or positive interpretations, whereas ELISA are still negative or only weakly positive. The *recomBlot* HCV IgG 2.0 shows the sensitivity level of comparable confirmatory tests at least.

Sera with a confirmed HCV status (i.e. with corresponding comparative test results) were used in calculation of sensitivity and specificity.

recomBlot HCV IgG 2.0	HCV status well defined		HCV status not defined	total
	negative	positive		
negative	605	1 ²⁾	85	691
positive	0	300	97	397
borderline	22	1	26	49
total	627	302	208	1137

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We will be pleased to send you further literature on HCV diagnostics at your request.

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MIKROGEN GmbH