recomLine Campylobacter IgG

recomLine Campylobacter IgA

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recomLine Campylobacter IgG

recomLine Campylobacter IgA

Line immunoassay based on recombinant antigens for the determination of IgG and IgA antibodies against *Campylobacter jejuni* and *Campylobacter coli* in human serum and plasma.

1. General aspects, intended use

recomLine Testkit Campylobacter is a qualitative in-vitro-Test for the detection of IgG- and IgA-antibodies against protein antigens of *Campylobacter jejuni* and *Campylobacter coli*.

In cases of previous or persistent Campylobacter infection with primary diagnosis based on stool sample culturing, recomLine Campylobacter provides for identification of specific Campylobacter antibodies by means of separate lineups of different antigens produced by recombinant engineering for the purpose of clarifying post-infection complications (see Sequelae).

2. *Campylobacter jejuni* and *Campylobacter coli*

The genus Campylobacter comprises gram-negative, spiral-shaped, microaerophilic, mesophilic to thermophilic bacteria with bipolar flagella. In 1963, Sebald and Vernon named a bacterium Escherich had described as early as 1889 *Campylobacter jejuni*. Isolation from stool samples was achieved in 1972 by Dekeyser et al. Jones et al. made out a food association in 1981/1982. In taxonomic terms, Campylobacter is classified with the epsilon subdivision of the Proteobacteria.

The pathogen reservoir is mainly the intestinal tract of warm-blooded wild, domestic and pet animals (birds and mammals).

Intestinal Campylobacter infections are the second most frequent enteric bacterial infections reported in Germany after enteric salmonellosis (Robert Koch Institute 2004 *RKI Infektionsepidemiologisches Jahrbuch meldepflichtiger Krankheiten 2004: Salmonella 56.947, Campylobacter 55.745*), whereby *Campylobacter jejuni* is much more frequent, accounting for over 90% of cases as compared to *Campylobacter coli* at approx. 9%. According to the Robert Koch Institute, Campylobacter infections have now become more frequent than Salmonella infections in 2005, an indication that incidence of such infections is increasing.

Contaminated foods (mainly poultry) and drinking water (tropical countries) constitute the main sources of infection. Unreported cases not reflected in the statistics probably outnumber reported cases many times over. Annual incidence rates vary seasonally (higher rates in the summer months) as well as regionally, with diagnoses following travel to warm countries throughout the year due to considerably higher incidence in Third World countries.

Acute disease

The course of infection with Campylobacter is almost exclusively oral, the pathogen is highly adapted to the intestinal tract. Systemic infections up to and including meningitis are observed less frequently (especially in immunocompromised patients, infants and the elderly) (Skirrow & Blaser, 2000). The incubation period is brief (1-7 days) and presumably depends on the infective dose. Besides nearly asymptomatic (clinically inapparent) courses, infected persons suffer from painful gastrointestinal symptoms with sometimes bloody diarrhoea, fever, meningism and myalgias. The acute clinical picture persists for only a few days in most patients. The course of the infection is self-limiting in most cases, although severe losses of fluids and electrolytes must be replaced as required.

Sequelae:

Reactive arthritis (ReA) is among the known sequelae to enteral, and frequently urogenital infections as well caused by certain bacterial and viral pathogens, with onset a few weeks after the primary infection. In addition to urogenital mycoplasmas and chlamydiae, enteral salmonellae, shigellae and yersiniae, *Campylobacter spp.* are prominent among secondary pathogen diagnoses (Locht & Krogfeldt, 2002; Cox et al., 2003; Hannu et al., 2004). The background pathophysiological mechanism is presumably molecular mimicry of cross-reacting antibodies to antigens of the synovial membrane. Some cases of arthritis of unclear genesis can be correctly diagnosed as post-infection Reactive Arthritis, and not rheumatoid arthritis, on the basis of the appropriate laboratory parameters.

Similarly to Reactive Arthritis, **Guillain-Barré-Syndrom (GBS)** may also develop a few weeks after infections, usually enteral or respiratory, caused by certain bacterial and viral pathogens. The presumed level of incidence is between 1:1,000 and 1:10,000. GBS is an acute, immunomediated polyradicular

neuropathy caused by an abnormal humoral immune response to the peripheral myelin sheath and/or the neural axon. A number of different viruses (e.g. CMV, EBV, VZV, measles virus, mumps virus), *Mycoplasma pneumoniae* and bacteria (*Borrelia burgdorferi*, *Haemophilus influenzae*) can cause a GBS, but various case-control studies have implicated that mainly *Campylobacter spp.*, and in particular *Campylobacter jejuni* (Nachamkin et al., 2000; Prendergast et al., 2004; Gilbert et al., 2004; Leonard et al., 2004) as prior agent with a frequency level of 30%, based mainly on serological finding. Direct pathogen detection methods produced negative results in most cases. In pathophysiological terms, pathogen-specific antibodies cross-react with neuronal antigens (e.g. Moran et al., 2000; Yuki et al., 2004), thereby inducing inflammation mediators to invade macrophages, resulting in subsequent local micro-oedemas with transitory or permanent failure of the affected neurones. Neurological failures in GBS therefore include purely motor ("AMAN" = acute motor axonal neuropathy, "AIDP" = acute inflammatory demyelinating polyneuropathy) as well as sensory and mixed types ("ANSAM" = acute motor-sensory axonal neuropathy). Mild and reversible courses are observed as well as severe courses with permanent paralysis. About 5% of GBS cases are terminal. The course of disease appears to be directly proportional to the specific antibody titre. Clinical observations and initial findings from surveillance studies on "flaccid paralysis" (differential diagnosis: poliomyelitis) suggest that, in pathophysiological terms, discrete, transitory paralysis symptoms that take the form of signs of exhaustion in the extremities correspond much more frequently to an inapparent GBS than has been assumed to date.

3. Diagnosis

In suspected *Campylobacter* infections, the bacteriological stool culture, preferably determined on selective mediums, is the standard diagnostic method (MIQ, Kist et al., 2000). However, the result or species-specific diagnosis may require two to three days due to the special demands of the pathogens. As is the case with other enteral bacterial infections, the results of serological antibody detection methods are not very definitive in acutely ill patients and are at most useful in aetiological clarification after the pathogens are no longer detectable. As a rule, however, *Campylobacter spp.* are still excreted sporadically during convalescence, in some cases over a period of several weeks, whereby negative and positive stool samples may alternate. Chronic shedders over periods of many months are observed, who may form a partial collective among patients with infection sequelae.

Campylobacter infections are notifiable according to the (German) Infection Protection Law (IfSG).

Due to the relatively high frequency of occurrence of these pathogens among the causes of post-infection complications, serological identification of antibodies to *Campylobacter spp.* is a necessary standard procedure in patients suffering from infection sequelae. The necessary tests can be based on complement-binding reactions (CBR), ELISA methods or immunoblots.

Campylobacter infections are still underestimated. Testing rates are generally low and usually involve CBR tests, whereby the use of the total protein does not ensure a very high level of specificity in view of the cross-reactivity involved. CBR can detect IgG and IgM antibodies. IgA antibodies cannot be detected, nor can IgG and IgM antibodies be differentiated.

The *recomLine Campylobacter* is based on recombinant *Campylobacter* antigens which contains the MOMP (major outer membrane protein), the PEB4 (cytoplasmic protein membrane protein), the PEB2 (periplasmic protein), the PEB1 (periplasmic protein, cell binding factor) and OMP18 (outer membrane protein, peptidoglycan-associated lipoprotein).

Recombinant Antigen	Abbreviation	Size of rec. Antigens [kDa]
major outer membrane protein	MOMP	46
cytoplasmic protein	PEB4	31
periplasmic protein	PEB2	27
periplasmic protein	PEB1	28
outer membrane protein	OMP18	18

The antigens used in this test facilitate the detection and identification of alternatively IgG resp. IgA antibodies specific to *Campylobacter jejuni* and *Campylobacter coli*. This allows a view on the reactivities against the most important antigens in only one test run.

4. Test principle

The recombinant, highly purified proteins are applied to a nitrocellulose membrane. This matrix is then cut into strips.

To facilitate detection of Campylobacter specific antibodies, the strips are incubated with the diluted serum or plasma sample, whereby the antibodies bind to the antigens on the strips. Unbound antibodies are then flushed away and the strips are incubated in a second step with anti-human IgG or IgA coupled with horse radish peroxidase. Specifically bound antibodies are detected by means of a colour reaction catalysed by the peroxidase. If an antigen-antibody reaction has taken place, a dark band appears at the corresponding locus on the strip.

Four control bands are arranged side-by-side at the upper end of the test strip:

1. The reaction control band under the strip number, which must show a reaction to every serum.
2. Two conjugate control bands: IgG/IgA. These bands serve as controls for the antibody class detected in each case. For instance, if the test strip is used to detect IgG antibodies, the conjugate control band IgG shows a clearly defined band. This applies analogously to the conjugate band IgA.
3. "Cutoff control": For control of the coloration process and evaluation of the test strip. The intensity of this band provides a basis for evaluation of the antibody reactivity as positive, borderline or negative.

5. Package contents

The reagents in a pack are sufficient for 25 determinations.

Each reagent set contains:

WASHBUF 10 X	100 ml	Wash buffer (ten times the concentration) Contains phosphate buffer, NaCl, KCl, detergent, preservatives: methylisothiazolone and oxyprione
SUBS TMB	45 ml	Substrate Solution Tetramethylbenzidin (TMB, ready-to-use)
MILKPOW	5 g	Skim milk powder
INSTRU	1	Instructions for use
EVALFORM	1	Evaluation form

5.1 recomLine Campylobacter IgG

Additionally to the components listed under point 5 each reagent set contains:

TESTSTR	1 piece	Test tubes with 25 consecutively numbered test strips, coated with recombinant Campylobacter -antigens
CONJ IgG	6 ml	Anti-human IgG conjugate (ten times the concentration) From rabbit, contains NaN ₃

5.2 recomLine Campylobacter IgA

Additionally to the components listed under point 5 each reagent set contains:

TESTSTR	1 piece	Test tubes with 25 consecutively numbered test strips, coated with recombinant Campylobacter -antigens
CONJ IgA	6 ml	Anti-human IgA conjugate (ten times the concentration) From rabbit, contains NaN ₃

6. Additional reagents and accessory equipment required

Deionised water, vacuum extraction system with disinfection trap, micro pipettes, plastic forceps, shaker, metric cylinder with graduations, scales.

Incubation trays (can be obtained from MIKROGEN GmbH)

7. Information on test and reagents

7.1 Precautions

- ☞ Suitable single-use gloves must be worn during the entire test procedure.
- ☞ The conjugates contain sodium azide. Avoid contact with skin or mucosa.
- ☞ All reagents and materials contaminated with potentially infectious samples must be treated with suitable disinfectants or autoclaved at 121°C for at least 1 hour.
- ☞ Use incubation trays only once.
- ☞ Handle strips carefully with a plastic forceps.

7.2 Handling information

Store reagents before and after use at 2°C - 8°C, **do not freeze**. Temper all components before starting the test for at least 30 minutes to 18°C - 25°C (room temperature). Both the test and incubation procedures are carried out at room temperature.

Mix the concentrated conjugates well before use. Mix the patient sera well before use as well.

Do not open the tube containing the test strips until just before use to avoid condensation of water. The strips that are not used must be left in the tube and stored further at 2°C - 8°C (reclose tube tightly, test strips must not become moist before testing!).

A quality guarantee can only be given up to the expiry date on the packages.

Only reagents with lot numbers corresponding to the respective lot number on the label of the kit package may be used.

The test must be performed by well-trained and authorised qualified personnel.

In case of substantial modifications of the product or the instructions for use, the application of the test might differ from the purpose intended by MIKROGEN.

Automation is possible, please refer to MIKROGEN for details.

7.3 Preparation of solutions

7.3.1 Preparation of ready-to-use wash buffer

This buffer is used for both serum and conjugate dilution as well as for the washing steps.

Prior to dilution, the volume of wash buffer required for the corresponding number of tests must be determined.

The skim milk powder is first dissolved in wash buffer concentrate. This mixture is then filled up with deionised water to the final volume (dilution: 1 + 9). See Table 1 for the amounts required. Volumes for numbers of test strips not listed in the table have to be determined by calculation.

Ready-to-use wash buffer can be stored at 2°C - 8°C for four weeks.

Table 1: Wash buffer per test strip used

Test strips used	Skim milk powder	Wash buffer concentrate	Deionised water	Ready-to-use wash buffer
1	0,1 g	2 ml	18 ml	20 ml
2	0,2 g	4 ml	36 ml	40 ml
3	0,3 g	6 ml	54 ml	60 ml
5	0,5 g	10 ml	90 ml	100 ml
10	1 g	20 ml	180 ml	200 ml
20	2 g	40 ml	360 ml	400 ml
25	2,5 g	50 ml	450 ml	500 ml

7.3.2 Preparation of conjugate solutions

The conjugate solution for the IgG resp. IgA test is to be prepared immediately before use. It is not possible to store the ready-to-use conjugate solution.

One part of IgG resp. IgA conjugate concentrate is diluted with 9 parts of ready-to-use wash buffer (1 + 9).

See Table 2 for the amounts required. Volumes for numbers of test strips not listed in the table have to be determined by calculation.

Table 2: Volumes for anti-human IgG resp. IgA conjugate dilution

Test strips used *	IgG / IgA conjugate concentrate	Ready-to-use wash buffer	Ready-to-use conjugate dilution
1	0,2 ml	1,8 ml	2 ml
2	0,4 ml	3,6 ml	4 ml
3	0,6 ml	5,4 ml	6 ml
5	1 ml	9 ml	10 ml
10	2 ml	18 ml	20 ml
20	4 ml	36 ml	40 ml
25	5 ml	45 ml	50 ml

* The volumes have been calculated without dead volume. Depending on handling (manually or by machine processing) please prepare conjugate solution for additional 1 to 3 strips.

7.3.3 Substrate solution

The substrate is ready for use! Bring to room temperature (18°C - 25°C) before starting the colour reaction.

Avoid contamination of the unused substrate solution by nonsterile pipette tips, etc. at all costs since this may affect test sensitivity.

7.4 Storage and stability

Store reagents at 2°C - 8°C before and after use.

Ready-to-use wash buffer can be stored at 2°C - 8°C for four weeks.

The conjugate solution is to be prepared always freshly.

8. Sample material

The sample material can be serum or plasma that is separated from the coagulum as soon as possible after sampling. A microbial contamination of the sample has to be avoided at all costs. Insoluble substances are to be removed from the sample prior to incubation by means of centrifugation.

Use of lipaemic, haemolytic or turbid samples may result in a dark background in the *recomLine* Campylobacter IgG/IgA.

Important!

If the tests are not carried out immediately, the samples can be stored for up to 2 weeks at 2°C-8°C. Longer storage of the samples is possible at -20°C or colder. Repeated freezing and thawing of the samples is not recommended because this may affect the quality of the results.

9. Test procedure

9.1 General

The reproducibility of the results depends to a great extent on consistent washing of the strips. The washing frequencies described under 9.3. and 9.5 should therefore always be maintained.

9.2 Incubation of samples

1. One well in the incubation tray (see 6) is required per sample to be tested. **2 ml** of ready-to-use wash buffer is pipetted into each well. Then one test strip is carefully dipped into each of the wells filled with ready-to-use wash buffer using a plastic forceps. The strip number must face upwards.

Important!

The strip must be completely wet and immersed in the ready-to-use wash buffer.

Record the tube and strip numbers used on the evaluation sheet.

2. Adding the samples

IgG/IgA test procedure: **20 µl** of an undiluted sample (human serum or plasma) is pipetted into the proper wells (**dilution 1 + 100**).

Please be sure to add the sample at one end of the immersed strips into the wash buffer and mix as soon as possible by shaking the tray carefully.

Record the sample numbers and the immunoglobulin class (IgG, IgA) on the evaluation sheet.

Cover the incubation tray with the plastic lid and incubate for **1 hour** at room temperature while shaking gently.

Important!

Make sure the incubation solutions are not carried over into other wells; be particularly careful to avoid splashing when opening and closing the lid.

9.3 Washing

1. Following incubation, the plastic lids are carefully removed from the incubation trays.
2. The serum dilution is carefully aspirated from the individual wells.

Important!

When the solutions have been aspirated from a well, the pipette tip has to be changed or rinsed well with deionised water after each aspiration procedure to prevent cross-contamination.

In the case of machine processing the references of the equipment manufacturer have to be considered.

3. Then place **2 ml** of the ready-to-use wash buffer in each well and wash on the shaker while shaking gently for **5 minutes**. The wash buffer is aspirated after the washing procedure.
4. Carry out the washing step under 3. a total of **three times**.

9.4 Incubation with peroxidase conjugate

After washing the strips, place **2 ml** of the appropriately prepared **conjugate solution** (see Table 2) into each well and incubate for **45 minutes** while shaking gently at room temperature, whereby the incubation tray is covered with the plastic lid.

9.5 Washing

The conjugate solutions are aspirated from the wells and the strips are washed again (see 9.3).

9.6 Substrate reaction

1. Add **1.5 ml substrate solution** into each well and incubate for **5 - 10 minutes**, shaking gently and under observation, at room temperature.
2. You should terminate the reaction as soon as the cutoff control band can be seen.

9.7 Stopping the reaction

1. After the solution is aspirated, wash the strips briefly **three times with deionised water**.
2. Use a plastic forceps to remove the strips carefully from the water and place them between 2 layers of absorbent paper to dry for about 2 hours. The strips may be adhesively attached to the enclosed evaluation sheet and the results may be recorded.
3. The strips should be stored protected from exposure to light.

9.8 Test procedure of avidity determination

The avidity test procedure is described in detail in the instructions for use for the avidity reagent (Article No.: 11010).

10. Summary of the test procedure

1.	bring all reagents to room temperature
2.	deposit the strips in 2 ml ready-to-use wash buffer and take care that they are completely immersed in the ready-to-use wash buffer
3.	pipette 20 µl of the sample
4.	incubate for 1 hour at room temperature while shaking gently
5.	wash on the shaker three times with 2 ml of ready-to-use wash buffer for 5 minutes each time
6.	add 2 ml appropriately prepared conjugate solution
7.	incubate at room temperature for 45 minutes while shaking gently
8.	wash on the shaker three times with 2 ml of ready-to-use wash buffer for 5 minutes each time
9.	add 1.5 ml of the substrate solution, incubate while shaking at room temperature for 5 - 10 minutes
10.	wash the strips at least three times with deionised water
11.	dry the strips between 2 layers of absorbent paper for 2 hours and read off the result

11. Evaluation

11.1 Evaluation of band intensity

1. On the enclosed evaluation sheet, record the date, lot and tube number along with the antibody class detected.
2. Enter the sample identification number on the protocol sheet.
3. Now attach the corresponding test strips with a glue stick into the corresponding fields of the evaluation sheet. To do this, place the test strips with the reaction control band on the marking lines. Then use clear adhesive tape to attach the test strips left from the marking line. Complete sticking of the test strip with glue or adhesive tape may lead to aberrations of the colouring.
4. Now identify the bands of the developed test strips based on the printed control strip on the evaluation sheet and record them in the protocol sheet (Table 3). Identify the bands on the test strips separately for the respective immunoglobulin classes.

11.2 Control results

The test can be evaluated if the following criteria are met:

- Reaction control band (upper line) pronounced colour, dark band
- Antibody class (second and third band): The corresponding conjugate control band must show a clear colour reaction. The other conjugate control band may develop a weak, non-specific colouring.
- Cutoff control (fourth band): weak, but clearly evident colouring.

Table 3: Intensity of bands in relation to the cutoff band

Bands	Intensity
No reaction	-
Intensity weaker than the cutoff band	±
Same intensity as the cutoff band	+
Strong intensity (stronger than the cutoff band)	++
Very strong intensity	+++

Important!

The band patterns in *recomLine* Campylobacter IgG and IgA detection tests may show varying intensities. It is possible that the IgG *recomLine* will show more intense and darker bands than the IgA *recomLine*. The intensity of the protein bands depends on the concentration of the Campylobacter-specific antibodies.

11.3 Test results and test interpretation

On the basis of clinical testing and a mathematical analysis a point evaluation of Campylobacter antigens in **recomLine Campylobacter** was carried out for safe and easy test evaluation. According to this, the test result can be obtained by adding up the points and subsequent evaluation.

The test result is obtained by the addition of the point values of the separate bands, that are determined as +, ++ or +++ (Table 4). The resulting total of points is entered into the column, marked with the sum symbol. Finally the positive, borderline or negative result are determined on the basis of the point value and recorded into the column "interpretation".

Table 4: Point evaluation of Campylobacter antigens in *recomLine* Campylobacter

Antigen	Points in IgG	Points in IgA
MOMP	3	3
PEB4	5	5
PEB2	1	1
PEB1	2	2
OMP18	5	5

Table 5: Testevaluation in *recomLine* Campylobacter

Sum of points	Evaluation
≤ 3	negative
4	borderline
≥ 5	positive

11.4 Notes on interpretation of the results

When evaluating the test results, IgG and IgA findings should always be considered together.

Samples with unclear or questionable results should be rechecked after 2 - 3 weeks in keeping with the clinical situation.

In all test interpretations, especially of weak positive results, it is important to include any other clinical record data available. Close cooperation between laboratory and the physician in charge of treatment is recommended here as well.

A negative result does not exclude the possibility of a campylobacter infection. False negative results may occur with serum samples obtained at a very early stage after infection in which no antibodies to *C.jejuni/C.coli* have been produced as yet.

Dark test strips: Some patient sera may cause a dark, continuous or patterned coloration on the entire nitrocellulose strip (e.g. sera from patients with milk protein allergies). Various different factors in each patient serum are responsible for this effect. Evaluation of these strips is usually feasible in a restricted sense only. For instance, "inverse" bands (white bands on a dark background) must be evaluated as negative. The corresponding serum should in any case be tested using other serological methods.

12. Clinical results

Patient samples from different sources were tested to assess the performance capability of recomLine Campylobacter IgG/IgA. The results are listed in the corresponding tables, correlated with additional data on the frequency of positive reactivity for the individual recombinant antigens.

12.1 Sera from patients with positive *C.jejuni* / *C.coli* stool culture

This study considered a total of 100 serum samples from patients with positive Campylobacter stool cultures. Blood samples were taken 0 to 40 days after the stool diagnostics, the beginning of the disease was unknown. The IgG results are shown in Table 6 and the IgA results in Table 8.

Following confirmed Campylobacter infections, 76% of patient samples were found to contain IgG antibodies and 22% IgA antibodies.

12.2 Blood donors

130 non-selected blood donor plasma samples were tested.

The results were a seroprevalence in the IgG test of 14% and in the IgA test of 1.5%, see Table 6 (IgG) and Table 8 (IgA). The only questionable result in this study was among the positive IgG results in the blood donor samples.

12.3 Samples from patients with suspected Guillain-Barré-Syndrome (GBS) of unclarified aetiology

Within the framework of this study, 30 samples from patients with suspected Guillain-Barré syndrome of unclarified aetiology were tested. These samples were characterised by the clinical picture of GBS, whereby no data on previous illnesses involving diarrhoea were known. IgA antibodies were detected 6.7 times more frequently in this sample collective than among the blood donors. No significant frequency percentiles of IgG antibodies were observed. The results are presented in Table 7 (IgG) and Table 9 (IgA).

12.4 Sera from patients with suspected Reactive Arthritis (ReA) of unclarified aetiology

Samples from 60 patients with suspected Reactive Arthritis of unclarified aetiology were tested for the presence of IgG and IgA antibodies to *C. jejuni* / *C. coli*. This collective was characterised by the clinical picture of Reactive Arthritis, whereby no data on previous illnesses involving diarrhoea were known. IgA antibodies were detected in this collective 5.3 times more frequently than in the blood donor collective. No significant frequency percentiles of IgG antibodies were observed. The results are presented in Table 7 (IgG) and Table 9 (IgA).

Table 6: recomLine Campylobacter IgG – Samples from patients with positive *C.jejuni* / *C. coli* stool culture and blood donors

	n	frequency of the rec. antigens (%)					Results
		MOMP	PEB4	PEB2	PEB1	OMP18	IgG-positive
<i>C.jejuni</i>, <i>C.coli</i> pos. stool culture	100	9	30	4	9	70	76%
Blood donors	130	5	4	2	2	13	14%

Table 7: recomLine Campylobacter IgG – Samples from patients with suspected Guillain-Barré-Syndrom resp. Reactive Arthritis unclarified aetiology

	n	frequency of the rec. antigens (%)					Results
		MOMP	PEB4	PEB2	PEB1	OMP18	IgG-positive
GBS-samples with unclarified aetiology	30	13	13	10	7	10	20%
ReA-samples with unclarified aetiology	60	0	17	8	2	12	23%

Table 8: recomLine Campylobacter IgA – Samples from patients with positive *C.jejuni* / *C. coli* stool culture and blood donors

	n	frequency of the rec. antigens (%)					Results
		MOMP	PEB4	PEB2	PEB1	OMP18	IgG-positive
<i>C.jejuni</i>, <i>C.coli</i> pos. stool culture	100	3	4	0	2	21	22%
Blood donors	130	1	1	0	0	1	1,5%

Table 9: recomLine Campylobacter IgA – Samples from patients with suspected Guillain-Barré-Syndrom resp. Reactive Arthritis unclarified aetiology

	n	frequency of the rec. antigens (%)					Results
		MOMP	PEB4	PEB2	PEB1	OMP18	IgG-positive
GBS-samples with unclarified aetiology	30	7	7	0	0	3	10%
ReA-samples with unclarified aetiology	60	2	7	0	2	2	8%

12.5 Potentially falsifying sera

Investigations of potential interference factors revealed no evidence that test results were affected in cases of haemolytic patient samples or in samples from patients that were positive for rheumatism factor or showed high total immunoglobulin levels.

False positive results may occur with lipaemic or icteric samples or samples from patients with infectious mononucleosis (EBV infection).

A total of 69 patient samples with potentially interfering factors were tested. The results are presented in Table 10.

Table 10: recomLine Campylobacter IgA /IgG, results of potentially falsifying sera

	IgG-positive	IgA-positive
lipaemic samples	2/10	0/10
icteric samples	1/10	1/10
haemolytic samples	1/10	0/10
Fresh EBV infections	2/10	0/10
high total immunoglobulin levels samples	0/9	0/10
rheuma factor positive samples	2/20	1/20

12.6 Serum-plasma comparison

Analysis of samples from patients revealed no differences between serum and plasma analyses .

13. Literature

- (1) Skirrow MB, Blaser MJ (2000), Clinical Aspects of Campylobacter Infection, in: Nachamkin I, Blaser MJ (eds), Campylobacter, 2nd Ed. ASM, Washington.
- (2) Loch H, Krogfeldt KA (2002), Comparison of rheumatological and gastrointestinal symptoms after infection with *Campylobacter jejuni/coli* and enterotoxigenic *Escherichia coli*, Ann Rheum Dis 61: 448-52.
- (3) Cox CJ, Kempell KE, Gaston JS (2003), Investigation of infectious agents associated with arthritis by reverse transcription PCR of bacterial rRNA, Arthritis Res Ther 5: R1-8.
- (4) Hannu T, Kauppi M, Tuomala M, Laaksonen I, Klements P, Kuusi M (2004), Reactive arthritis following an outbreak of Campylobacter jejuni infection, J Rheumatol 31: 528-30.
- (5) Nachamkin I et al. (2000), *Campylobacter jejuni* Infection and the Association with Guillain-Barré Syndrome, in: Nachamkin I, Blaser MJ (eds), Campylobacter, 2nd Ed. ASM, Washington.
- (6) Prendergast MM, Tribble DR, Baqar S, Scott DA, Ferris JA, Walker RI, Moran AP (2004), In Vivo Phase Variation and Serologic Response to Lipooligosaccharide of *Campylobacter jejuni* in Experimental Human Infection, Infect Immun 72: 916-22.
- (7) Gilbert M, Godschalk PC, Karwaski MF, Ang CW, van Belkum A, Li J, Wakarchuk WW, Endtz HP, (2004), Evidence for Acquisition of the Lipooligosaccharide Biosynthesis Locus in *Campylobacter jejuni* GB11, a Strain Isolated from a Patient with Guillain-Barré Syndrome, by Horizontal Exchange, Infect Immun 72: 1162-65.
- (8) Leonard EE 2nd, Tompkins LS, Falkow S, Nachamkin I, (2004), Comparison of *Campylobacter jejuni* Isolated Implicated in Guillain-Barré Syndrome and Strains That Cause Enteritis by a DNA Microarray, Infect Immun 72: 1199-203.
- (9) Moran AP, Penner JL, Aspinall GO (2000), Campylobacter Lipopolysaccharides, in: Nachamkin I, Blaser MJ (eds), Campylobacter, 2nd Ed. ASM, Washington.
- (10) Yuki N, Susu+ki K, Koga M, Nishimoto Y, Odaka M, Hirata K, Tagushi K, Miyatake T, Furukawa K, Kobata T, Yamada M, (2004), Carbohydrate mimicry between human ganglioside GM1 and Campylobacter jejuni lipooligosaccharide causes Guillain-Barré syndrome, Proc Natl Acad Sci USA 101: 11404-9.
- (11) Kist M, Bockemühl J (Hersg.), (2000), Gastrointestinale Infektionen, MIQ (Qualitätsstandards in der mikrobiologischen Diagnostik, Urban & Fischer, München.
- (12) Heintschel v. Heinegg E. (1997), Diagnostische Bibliothek Campylobacter, Blackwell Wissenschafts-Verlag Laboratoriumsmedizin 52.
- (13) Kist M. (2000), Lebensmittelbedingte Infektionen durch Campylobacter, Springer-Verlag, Bundesgesundheitsbl-Gesundheitsforsch-Gesundheitsschutz 45, 497-506.
- (14) Altekruse S. F., Stern N. J., Fields P. I., Swerdlow D. L. (1999), Campylobacter jejuni-An Emerging Foodborne Pathogen, Emerging Infectious Diseases 5/1.
- (15) Zhang Q., Meitzler J. C., Huang S., Morishita T. (2000), Sequence Polymorphism, Predicted Secondary Structures, and Surface-Exposed Conformational Epitopes of Campylobacter Major Outer Membrane Protein, Infection and Immunity 68/10, 5679-5689.
- (16) Burnens A., Stucki U., Nicolet J., Frey J. (1995), Identification and Characterization of an Immunogenic Outer Membrane Protein of Campylobacter jejuni, J Clin Microbiol 33/11 2826-2832.
- (17) Pei Z., Ellison R. T., Blaser M. J. (1991), Identification, Purification, and Characterization of Major Antigenic Proteins of *Campylobacter jejuni*, J Biological Chemistry 266/25, 16363-16369.
- (18) Ketley J. M. (1997), Pathogenesis of enteric infection by Campylobacter, Microbiology 143, 5-21.
- (19) Nachamkin I., Allos B. M., Ho T. (1998), Campylobacter Species and Guillain-Barré-Syndrome, Clinical Microbiology Reviews 11/3, 555-567.
- (20) Nachamkin I. (2002), Chronic effects of Campylobacter infection, Microbes and Infection 4, 399-403.
- (21) Colmegna I., Cuchacovich R., Espinoza L.R., (2004), HLA-B27-Associated Reactive Arthritis: Pathogenetic and Clinical Considerations, Clinical Microbiology Reviews 17/2, 348-369
- (22) Lecuit et al. (2004), Immunoproliferative Small Intestinal Disease Associated with Campylobacter jejuni, N Engl J Med 350: 239-48
- (23) Yli-Kerttula T., Luukkainen R., Yli-Kerttula U., Möttönen T., Hakola M., Korpela M. Sanila M., Uksila J., Toivanen A., (2003), Effect of a three month course of ciprofloxacin on the late prognosis of reactive arthritis, Ann Rheum Dis, 62:880-884,

We will be pleased to send you further literature on Campylobacter diagnostics at your request.

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