recomWell Campylobacter IgG

recomWell Campylobacter IgA

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recomWell Campylobacter IgG

recomWell Campylobacter IgA

Enzyme immunoassay with recombinant antigens for the detection of IgG or IgA antibodies against *Campylobacter jejuni* and *Campylobacter coli* in human serum or plasma

1. General aspects, intended use

recomWell Campylobacter is a quantitative in vitro test for the detection and safe identification of IgG or IgA antibodies against *Campylobacter jejuni* and *Campylobacter coli* in human serum or plasma. **recomWell Campylobacter** is a screening test based on the principle of an indirect sandwich ELISA.

2. Campylobacter jejuni and Campylobacter coli

The genus *Campylobacter* comprises gram-negative, spiral-shaped, microaerophilic, mesophilic to thermophilic bacteria with bipolar flagella. In 1963, Sebald and Vernon named a bacterium Escherich had described as early as 1889 *Campylobacter jejuni*. Isolation from stool samples was achieved in 1972 by Dekeyser et al. Jones et al. made out a food association in 1981/1982. In taxonomic terms, *Campylobacter* is classified with the epsilon subdivision of the Proteobacteria.

The pathogen reservoir is mainly the intestinal tract of warm-blooded wild, domestic and pet animals (birds and mammals).

Intestinal *Campylobacter* infections are the second most frequent enteric bacterial infections reported in Germany after enteric salmonellosis (Robert Koch Institute 2004: *RKI Infektionsepidemiologisches Jahrbuch meldepflichtiger Krankheiten 2004: Salmonella 56.947, Campylobacter 55.745*), whereby *Campylobacter jejuni* is much more frequent, accounting for over 90% of cases as compared to *Campylobacter coli* at approx. 9%. According to the Robert Koch Institute, *Campylobacter* infections have now become more frequent than *Salmonella* infections in 2005, an indication that incidence of such infections is increasing.

Contaminated foods (mainly poultry) and drinking water (tropical countries) constitute the main sources of infection. Unreported cases not reflected in the statistics probably outnumber reported cases many times over. Annual incidence rates vary seasonally (higher rates in the summer months) as well as regionally, with diagnoses following travel to warm countries throughout the year due to considerably higher incidence in Third World countries.

Acute disease

The course of infection with *Campylobacter* is almost exclusively oral, the pathogen is highly adapted to the intestinal tract. Systemic infections up to and including meningitis are observed less frequently (especially in immunocompromised patients, infants and the elderly) (Skirrow & Blaser, 2000). The incubation period is brief (1-7 days) and presumably depends on the infective dose. Besides nearly asymptomatic (clinically inapparent) courses, infected persons suffer from painful gastrointestinal symptoms with sometimes bloody diarrhoea, fever, meningism and myalgias. The acute clinical picture persists for only a few days in most patients. The course of the infection is self-limiting in most cases, although severe losses of fluids and electrolytes must be replaced as required.

Sequelae:

Reactive arthritis (ReA) is among the known sequelae to enteral, and frequently urogenital infections as well caused by certain bacterial and viral pathogens, with onset a few weeks after the primary infection. In addition to urogenital mycoplasmas and chlamydiae, enteral salmonellae, shigellae and yersiniae, *Campylobacter spp.* are prominent among secondary pathogen diagnoses (Locht & Krogfeldt, 2002; Cox et al., 2003; Hannu et al., 2004). The background pathophysiological mechanism is presumably molecular mimicry of cross-reacting antibodies to antigens of the synovial membrane. Some cases of arthritis of unclear genesis can be correctly diagnosed as post-infection Reactive Arthritis, and not rheumatoid arthritis, on the basis of the appropriate laboratory parameters.

Similarly to Reactive Arthritis, **Guillain-Barré-Syndrom (GBS)** may also develop a few weeks after infections, usually enteral or respiratory, caused by certain bacterial and viral pathogens. The presumed level of incidence is between 1:1,000 and 1:10,000. GBS is an acute, immunomediated polyradicular neuropathy caused by an abnormal humoral immune response to the peripheral myelin sheath and/or the neural axon. A number of different viruses (e.g. CMV, EBV, VZV, measles virus, mumps virus), *Mycoplasma pneumoniae* and bacteria (*Borrelia burgdorferi*, *Haemophilus influenzae*) can cause a GBS,

but various case-control studies have implicated that mainly *Campylobacter spp.*, and in particular *Campylobacter jejuni* (Nachamkin et al., 2000; Prendergast et al., 2004; Gilbert et al., 2004; Leonard et al., 2004) as prior agent with a frequency level of 30%, based mainly on serological finding. Direct pathogen detection methods produced negative results in most cases. In pathophysiological terms, pathogen-specific antibodies cross-react with neuronal antigens (e.g. Moran et al., 2000; Yuki et al., 2004), thereby inducing inflammation mediators to invade macrophages, resulting in subsequent local micro-oedemas with transitory or permanent failure of the affected neurones. Neurological failures in GBS therefore include purely motor ("AMAN" = acute motor axonal neuropathy, "AIDP" = acute inflammatory demyelinating polyneuropathy) as well as sensory and mixed types ("ANSAM" = acute motor-sensory axonal neuropathy). Mild and reversible courses are observed as well as severe courses with permanent paralysis. About 5% of GBS cases are terminal. The course of disease appears to be directly proportional to the specific antibody titre. Clinical observations and initial findings from surveillance studies on "flaccid paralysis" (differential diagnosis: poliomyelitis) suggest that, in pathophysiological terms, discrete, transitory paralysis symptoms that take the form of signs of exhaustion in the extremities correspond much more frequently to an inapparent GBS than has been assumed to date.

3. Diagnosis

In suspected *Campylobacter* infections, the bacteriological stool culture, preferably determined on selective mediums, is the standard diagnostic method (MIQ, Kist et al., 2000). However, the result or species-specific diagnosis may require two to three days due to the special demands of the pathogens. As is the case with other enteral bacterial infections, the results of serological antibody detection methods are not very definitive in acutely ill patients and are at most useful in aetiological clarification after the pathogens are no longer detectable. As a rule, however, *Campylobacter spp.* are still excreted sporadically during convalescence, in some cases over a period of several weeks, whereby negative and positive stool samples may alternate. Chronic shedders over periods of many months are observed, who may form a partial collective among patients with infection sequelae.

Campylobacter infections are notifiable according to the (German) Infection Protection Law (IfSG).

Due to the relatively high frequency of occurrence of these pathogens among the causes of post-infection complications, serological identification of antibodies to *Campylobacter spp.* is a necessary standard procedure in patients suffering from infection sequelae. The necessary tests can be based on complement-binding reactions (CBR), ELISA methods or immunoblots.

Campylobacter infections are still underestimated. Testing rates are generally low and usually involve CBR tests, whereby the use of the total protein does not ensure a very high level of specificity in view of the cross-reactivity involved. CBR can detect IgG and IgM antibodies. IgA antibodies cannot be detected, nor can IgG and IgM antibodies be differentiated.

4. Test principle

The antigens used in **recomWell Campylobacter** are made by genetic engineering techniques. In this way, an optimal presentation without other interfering and cross reacting proteins can be achieved. Only specific antigens with importance for the sensitive serological diagnosis are used in comparison to other ELISA-systems using lysate antigens.

The following recombinant antigens are presented in this **recomWell Campylobacter**: PEB4 (cytoplasmic protein, membrane protein) and OMP18 (outer membrane protein, peptidoglycan-associated lipoprotein).

The determination of IgG and IgA antibodies will be practised separately.

5. Package contents

The reagents in one pack are sufficient for 96 determinations.

Each reagent set contains:

| | | |
|---------------------|----------|---|
| WASHBUF 10 X | 100 ml | Wash buffer (ten times the concentration) Contains phosphate buffer, NaCl and detergent, Preservatives: methylisothiazolone and oxyprirone |
| DILUBUF | 125 ml | Dilution buffer (ready-for-use) Contains protein, detergent and blue dye Preservatives: methylisothiazolone and oxyprirone |
| SUBS TMB | 12 ml | Chromogenic substrate tetramethylbenzidine (TMB, ready-for-use) |
| SOLN STOP | 12 ml | Stop solution 25 % Phosphoric acid (H ₃ PO ₄) |
| INSTRU | 1 | Instructions for use |
| EVALFORM | 1 | Evaluation form |
| TAPE | 2 pieces | Sealing tape |

5.1. recomWell Campylobacter IgG

Additionally to the components listed under point 5 each reagent set contains:

| | | |
|----------------------|--------------|---|
| MTP | 12 x 8 wells | Microplate (section marked in red) coated with recombinant Campylobacter antigens in vacuum-pressure sealed bag |
| CONTROL + IgG | 150 µl | Positive control (violet screw cap) Preservatives: methylisothiazolone and oxyprirone |
| CONTROL ± IgG | 150 µl | Cutoff control (yellow screw cap) Preservatives: methylisothiazolone and oxyprirone |
| CONTROL - IgG | 150 µl | Negative control (white screw cap) Preservatives: methylisothiazolone and oxyprirone |
| CONJ IgG | 250 µl | Anti-human IgG conjugate (101 times the concentration, red screw cap) Preservatives: methylisothiazolone, sodium azide and chloroacetamide |

5.2. recomWell Campylobacter IgA

Additionally to the components listed under point 5 each reagent set contains:

| | | |
|----------------------|--------------|--|
| MTP | 12 x 8 wells | Microplate (section marked in blue) coated with recombinant Campylobacter antigens in vacuum-pressure sealed bag |
| CONTROL + IgA | 150 µl | Positive control (brown screw cap) Preservatives: methylisothiazolone and oxyprirone |
| CONTROL ± IgA | 150 µl | Cutoff control (orange screw cap) Preservatives: methylisothiazolone and oxyprirone |
| CONTROL - IgA | 150 µl | Negative control (white screw cap) Preservatives: methylisothiazolone and oxyprirone |
| CONJ IgA | 250 µl | Anti-human IgA conjugate (101 times the concentration, blue screw cap) Preservatives: methylisothiazolone, sodium azide and chloroacetamide |

6. Additional reagents and accessory equipment required

Deionised water, test tubes, micro pipettes, incubator 37°C, microplate photometer

7. Information on test and reagents

7.1. Precautions

- ☞ Control sera are from blood donors verified for the absence of antibodies to HIV 1/2, HCV and no Hepatitis Bs-antigen. Since an infection cannot be excluded with absolute certainty despite this precaution, the product must be treated with the same care as the patient sample.
- ☞ Suitable single-use gloves must be worn during the entire test procedure.
- ☞ The conjugates, controls, dilution buffer and wash buffer contain preservative agents. Avoid contact with skin or mucosa.
- ☞ Phosphoric acid is an irritant. Avoid all contact with skin or mucosa.
- ☞ All reagents and materials contaminated with potentially infectious samples must be treated with suitable disinfectants or autoclaved at 121°C for at least 1 hour.

7.2. Handling information

A quality guarantee can only be given up to the expiry date on the packages.

Only reagents with lot numbers corresponding to the respective lot number on the label of the kit package may be used.

The test must be performed by well-trained and authorised qualified personnel.

In case of substantial modifications of the product or the instructions for use, the application of the test might differ from the purpose intended by MIKROGEN.

Automation is possible, please refer to MIKROGEN for details.

7.3. Preparation of the solutions

The test reagents are sufficient for 96 IgG or IgA tests. The amounts indicated below refer to processing of a microplate strip with 8 wells. If several microplate strips are used at the same time, the amounts indicated must be multiplied by the number of microplate strips used.

Substrate and stop solution are ready to use.

7.3.1. Preparation of ready-to-use wash buffer

The wash buffer concentrate is diluted **1 + 9** with deionised H₂O. 5 ml concentrate are mixed with 45 ml deionised H₂O per microplate strip with 8 wells.

7.3.2. Preparation of conjugate solution

Per each microplate strip with 8 wells, 10 µl anti-human IgG peroxidase conjugate (red cap) or IgA peroxidase conjugate (blue cap) are added to 1 ml dilution buffer in a clean vessel and mixed well (dilution **1 + 100**).

7.4. Storage and stability

Store the reagents at **2°C - 8°C** before and after use.

The ready-to-use wash buffer can be prepared in larger amounts. Ready-to-use wash buffer may be stored at room temperature for one week for use in further tests.

The sample dilutions, controls and conjugate solution must always be prepared freshly.

8. Sample material

The sample material can be serum or plasma that is separated from the coagulum as soon as possible after sampling. Heat-inactivated samples will result in raised background reaction levels and are therefore unsuitable for use. A microbial contamination of the sample has to be avoided at all costs. Insoluble substances must be removed from the sample prior to incubation by centrifugation.

Important!

If the tests are not carried out immediately, the samples can be stored for up to 2 weeks at 2°C - 8°C. Longer storage of the samples is possible at -20°C or below. Repeated freezing and thawing of the samples is not recommended because this may affect the quality of the results.

9. Test procedure

9.1. Test preparations

Temper all reagents to room temperature (18°C - 25°C) before use for about 30 minutes. To avoid condensation of water in the microplate, it must be brought up to room temperature **in the closed bag**. After the required section is removed, the bag containing the plate must be reclosed and placed in the refrigerator.

Before use, the control and patient sera and concentrated conjugates must be mixed well and briefly centrifuged if practicable to collect the liquid at the bottom of the tubes.

9.2. Preparation of samples and controls

Pipette 10 µl of sample or control into 1 ml dilution buffer each and mix well (dilution **1 + 100**). A negative control, cutoff control and positive control must be run parallel to each test run and diluted just like the patient samples.

9.3. Incubation of the samples

The microplate is removed from the pressure-sealed bag.

Pipette **100 µl** per well of the diluted samples and diluted controls. One value is tested for each negative control, positive control and patient sample, whereas the cutoff control must be double-tested. It is preferable to pipette one cutoff control at the beginning and again at the end of the series.

In manual processing, the microplate is **carefully taped** over with sealing tape and incubated for **1 hour** at **37°C**.

9.4. Washing procedure

The wells are emptied completely, then washed **four times**, each time with 300 µl ready-to-use wash buffer per well. We recommend performing this step with the appropriate ELISA washing equipment. Make absolutely sure that the wash buffer is removed completely between the washing steps.

After the last washing step is completed (even if washing equipment is being used) tap the plate over a paper towel to remove any residual liquid from the wells.

9.5. Incubation with peroxidase conjugate

Pipette **100 µl** per well.

In manual processing, the microplate is **carefully taped** over with sealing tape and incubated for **30 minutes** at **37°C**.

9.6. Washing procedure

The wells are emptied and washed as described under 9.4.

9.7. Substrate reaction

The substrate solution is ready to use.

Pipette **100 µl** per well.

The microplate is incubated for **30 minutes** at **room temperature** while protecting it from direct sunlight. The time is counted from pipetting of the first well.

9.8. Stopping the reaction

To stop the reaction, pipette **100 µl** of stop solution per well, using the same pipetting scheme as for the substrate solution.

11.2.2. Quantitative evaluation

The antibody activity levels in units per ml are assigned to the extinction values using a formula.

$U/ml \text{ sample} = (\text{extinction sample} / \text{extinction cutoff}) \times 20$

Grey range: Lower limit = 20 U/ml

Upper limit = 24 U/ml

- U/ml sample > 24 **positive** test result
- U/ml sample < 20 **negative** test result
- $20 \leq U/ml \text{ sample} \leq 24$ **borderline** test result

11.3. Directions for the interpretation of test results

In all test interpretations it is important to include any other clinical record data available. Close cooperation between laboratory and the physician in charge of treatment is recommended here as well.

Samples with unclear or questionable results should be rechecked after 2 - 3 weeks in keeping with the clinical situation.

A negative *recomWell* Campylobacter result does not exclude the possibility of a campylobacter infection. False negative results may occur with serum samples obtained at a very early stage after infection in which no antibodies to *C.jejuni/C.coli* have been produced as yet.

When evaluating the test results, IgG and IgA findings should always be considered together.

We generally recommend to check positive and borderline ELISA results by a confirmation test (e.g. *recomLine* Campylobacter).

12. Clinical results

12.1. Sera from patients with positive *C.jejuni* / *C.coli* stool culture

This study considered a total of 292 serum samples from patients with positive Campylobacter stool cultures. Blood samples were taken 0 to 81 days after the stool diagnostics, the beginning of the disease was unknown.

Following confirmed Campylobacter infections, 73% of patient samples were found to contain IgG antibodies and 24% IgA antibodies (see Table 1).

This 292 samples were tested with a reference ELISA (extract-antigen) also. Only 43% of the samples were IgG- positive or borderline in the reference ELISA. 109 samples (37%) showed a positive or borderline result in both testing systems.

89 samples of 292 (30%) samples were positive in the reference IgA-ELISA. Only 28 samples (13%) showed a positive or borderline result in both testing systems.

12.2. Samples from patients with suspected Guillain-Barré-Syndrome (GBS) of unclarified aetiology

Within the framework of this study, 30 samples from patients with suspected Guillain-Barré syndrome of unclarified aetiology were tested. These samples were characterised by the clinical picture of GBS, whereby no data on previous illnesses involving diarrhoea were known. IgA antibodies were detected 5 times more frequently in this sample collective than among the blood donors. No significant frequency percentiles of IgG antibodies were observed. The results are presented in Table 2.

12.3. Sera from patients with suspected Reactive Arthritis (ReA) of unclarified aetiology

Samples from 38 patients with suspected Reactive Arthritis of unclarified aetiology were tested for the presence of IgG and IgA antibodies to *C. jejuni* / *C. coli*. This collective was characterised by the clinical picture of Reactive Arthritis, whereby no data on previous illnesses involving diarrhoea were known. IgA antibodies were detected in this collective 4 times more frequently than in the blood donor collective. No significant frequency percentiles of IgG antibodies were observed. The results are presented in Table 2.

12.4. Blood donors

220 non-selected blood donor plasma samples were tested.

The results were a seroprevalence in the IgG test of 19% and in the IgA test of 2,7% (see Table 1).

Table 1: recomWell Campylobacter - Samples from patients with positive *C.jejuni* / *C. coli* stool culture and blood donors

| | n | positive or borderline interpretation | |
|--|------------|---------------------------------------|-----------------|
| | | IgG | IgA |
| <i>C.jejuni</i>, <i>C.coli</i> positive stool culture | 292 | 213 (73%) | 69 (24%) |
| Blood donors | 220 | 41 (19%) | 6 (2,7%) |

Table 2: recomWell Campylobacter - Samples from patients with suspected Guillain-Barré-Syndrom resp. Reactive Arthritis unclarified aetiology

| | n | positive or borderline interpretation | |
|---|-----------|---------------------------------------|----------------|
| | | IgG | IgA |
| GBS-samples with unclarified aetiology | 30 | 7 (23%) | 4 (13%) |
| ReA-samples with unclarified aetiology | 38 | 7 (18%) | 4 (11%) |

12.5. Potentially falsifying sera

Sera of patients (n=100) infected by *Yersinia enterocolitica*/*Y. pseudotuberculosis*, *Helicobacter pylori*, *Borrelia burgdorferi*, *Parvovirus B19* and *Toxoplasma gondii* did not show a higher number of increased Campylobacter antibody titers.

12.6. Serum-plasma comparison

Analysis of samples from patients revealed no differences between serum and plasma analyses.

12.7. Precision

12.7.1. Intra-Assay-Variance:

A serum was examined on 32 cavities of a microtiterplate. The variation coefficient (VC) was calculated:
VC (IgG) = 6%, VC (IgA) = 7%.

12.7.2. Inter-Assay-Variance:

7 sera of different ODs were examined in 8 different determinations. The variation coefficient was calculated for each single serum (U/ml):

VC (IgG) < 12%, VC (IgA) < 12 %.

13. Literature

- (1) Skirrow MB, Blaser MJ (2000), Clinical Aspects of Campylobacter Infection, in: Nachamkin I, Blaser MJ (eds), Campylobacter, 2nd Ed. ASM, Washington.
- (2) Loch H, Krogfeldt KA (2002), Comparison of rheumatological and gastrointestinal symptoms after infection with *Campylobacter jejuni/coli* and enterotoxigenic *Escherichia coli*, Ann Rheum Dis 61: 448-52.
- (3) Cox CJ, Kempell KE, Gaston JS (2003), Investigation of infectious agents associated with arthritis by reverse transcription PCR of bacterial rRNA, Arthritis Res Ther 5: R1-8.
- (4) Hannu T, Kauppi M, Tuomala M, Laaksonen I, Klements P, Kuusi M (2004), Reactive arthritis following an outbreak of Campylobacter jejuni infection, J Rheumatol 31: 528-30.
- (5) Nachamkin I et al. (2000), *Campylobacter jejuni* Infection and the Association with Guillain-Barré Syndrome, in: Nachamkin I, Blaser MJ (eds), Campylobacter, 2nd Ed. ASM, Washington.
- (6) Prendergast MM, Tribble DR, Baqar S, Scott DA, Ferris JA, Walker RI, Moran AP (2004), In Vivo Phase Variation and Serologic Response to Lipooligosaccharide of *Campylobacter jejuni* in Experimental Human Infection, Infect Immun 72: 916-22.
- (7) Gilbert M, Godschalk PC, Karwaski MF, Ang CW, van Belkum A, Li J, Wakarchuk WW, Endtz HP, (2004), Evidence for Acquisition of the Lipooligosaccharide Biosynthesis Locus in *Campylobacter jejuni* GB11, a Strain Isolated from a Patient with Guillain-Barré Syndrome, by Horizontal Exchange, Infect Immun 72: 1162-65.
- (8) Leonard EE 2nd, Tompkins LS, Falkow S, Nachamkin I, (2004), Comparison of *Campylobacter jejuni* Isolated Implicated in Guillain-Barré Syndrome and Strains That Cause Enteritis by a DNA Microarray, Infect Immun 72: 1199-203.
- (9) Moran AP, Penner JL, Aspinall GO (2000), Campylobacter Lipopolysaccharides, in: Nachamkin I, Blaser MJ (eds), Campylobacter, 2nd Ed. ASM, Washington.
- (10) Yuki N, Susu+ki K, Koga M, Nishimoto Y, Odaka M, Hirata K, Tagushi K, Miyatake T, Furukawa K, Kobata T, Yamada M, (2004), Carbohydrate mimicry between human ganglioside GM1 and Campylobacter jejuni lipooligosaccharide causes Guillain-Barré syndrome, Proc Natl Acad Sci USA 101: 11404-9.
- (11) Kist M, Bockemühl J (Hersg.), (2000), Gastrointestinale Infektionen, MIQ (Qualitätsstandards in der mikrobiologischen Diagnostik, Urban & Fischer, München.
- (12) Heintschel v. Heinegg E. (1997), Diagnostische Bibliothek Campylobacter, Blackwell Wissenschafts-Verlag Laboratoriumsmedizin 52.
- (13) Kist M. (2000), Lebensmittelbedingte Infektionen durch Campylobacter, Springer-Verlag, Bundesgesundheitsbl-Gesundheitsforsch-Gesundheitsschutz 45, 497-506.
- (14) Altekruse S. F., Stern N. J., Fields P. I., Swerdlow D. L. (1999), Campylobacter jejuni-An Emerging Foodborne Pathogen, Emerging Infectious Diseases 5/1.
- (15) Zhang Q., Meitzler J. C., Huang S., Morishita T. (2000), Sequence Polymorphism, Predicted Secondary Structures, and Surface-Exposed Conformational Epitopes of Campylobacter Major Outer Membrane Protein, Infection and Immunity 68/10, 5679-5689.
- (16) Burnens A., Stucki U., Nicolet J., Frey J. (1995), Identification and Characterization of an Immunogenic Outer Membrane Protein of Campylobacter jejuni, J Clin Microbiol 33/11 2826-2832.
- (17) Pei Z., Ellison R. T., Blaser M. J. (1991), Identification, Purification, and Characterization of Major Antigenic Proteins of *Campylobacter jejuni*, J Biological Chemistry 266/25, 16363-16369.
- (18) Ketley J. M. (1997), Pathogenesis of enteric infection by Campylobacter, Microbiology 143, 5-21.
- (19) Nachamkin I., Allos B. M., Ho T. (1998), Campylobacter Species and Guillain-Barré-Syndrome, Clinical Microbiology Reviews 11/3, 555-567.
- (20) Nachamkin I. (2002), Chronic effects of Campylobacter infection, Microbes and Infection 4, 399-403.
- (21) Colmegna I., Cuchacovich R., Espinoza L.R., (2004), HLA-B27-Associated Reactive Arthritis: Pathogenetic and Clinical Considerations, Clinical Microbiology Reviews 17/2, 348-369
- (22) Lecuit et al. (2004), Immunoproliferative Small Intestinal Disease Associated with Campylobacter jejuni, N Engl J Med 350: 239-48
- (23) Yli-Kerttula T., Luukkainen R., Yli-Kerttula U., Möttönen T., Hakola M., Korpela M. Sanila M., Uksila J., Toivanen A., (2003), Effect of a three month course of ciprofloxacin on the late prognosis of reactive arthritis, Ann Rheum Dis, 62:880-884,

We will be pleased to send you further literature on the diagnosis of Campylobacter at your request.

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MIKROGEN GmbH Tel: +49 (0)89 54 80 1-0